

Equipment Covers

ECO² Setup

CODE	DESCRIPTION (LxW)	CASE QTY	CASE DIMENSIONS
E41-4470	ECO ² Setup Cover, 44x70" (112x178cm)	50	25x14x9" (63x35x23cm)
E41-4470-S	ECO ² Setup Cover, 44x70" (112x178cm)	40	25x14x9" (63x35x23cm)
E41-4477	ECO ² Setup Cover, 44x77" (112x196cm)	50	25x14x9" (63x35x23cm)
E41-4477-S	ECO ² Setup Cover, 44x77" (112x196cm)	40	25x14x9" (63x35x23cm)
E41-4490	ECO ² Setup Cover, 44x90" (112x229cm)	100	25x14x9" (63x35x23cm)
E41-4490-S	ECO ² Setup Cover, 44x90" (112x229cm)	25	25x14x9" (63x35x23cm)
E41-5470	ECO ² Setup Cover, 54x70" (137x178cm)	100	20x15x20" (52x39x50cm)
E41-5470-S	ECO ² Setup Cover, 54x70" (137x178cm)	25	20x15x20" (52x39x50cm)
E41-5480	ECO ² Setup Cover, 54x80" (137x203cm)	50	20x15x20" (52x39x50cm)
E41-5480-S	ECO ² Setup Cover, 54x80" (137x203cm)	40	20x15x20" (52x39x50cm)
E41-5490	ECO ² Setup Cover, 54x90" (137x229cm)	40	18x12x6" (47x30x15cm)
E41-5490-S	ECO ² Setup Cover, 54x90" (137x229cm)	15	16x12x16" (41x31x41cm)

Materials

- ECO² Clear Polyethylene (PE) Film

Product Attributes

- Clear film allows for visibility
- Certain components of this product incorporate ECO² sustainable materials which enhance product biodegradability

Sterilization

Sterile: Sterilized by EO (ethylene oxide). Product is sterile unless the package has been damaged or opened.
Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packer. This includes further packaging and sterilization according to the validated processes of the re-packer. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

Shelf Life

Sterile: 5 years
Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

Product Certifications/Standards


- Sterile: CE certified as Class I, according to MDD93/42/EEC Annex IX, Rule 1
- Non-Sterile: CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1
- ASTM D5511-12 Standard Test Method for Determining Anaerobic Biodegradation of Plastics Materials Under High-Solids Anaerobic Digestion Condition

Manufacturing Certifications/Registrations

- EN ISO 13485:2016
- Medical Device Single Audit Process Certified

Storage

- General warehouse conditions



Alleset
 4142 Industry Way
 Flowery Branch, GA 30542,
 USA
www.alleset.com
 SM-TDS-003 2023
 Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.

GRI Medical & Electronic Technology Co., LTD
 1805 Honggao Road,
 314031 Jiaxing, Zhejiang
 PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
 US FDA Establishment Registration #3004911384

Emergo Europe B.V.
 Westervoortsedijk 60,
 6827 AT Arnhem,
 The Netherlands
www.emergobyul.com

