

# Equipment Covers

## Setup

CODE	DESCRIPTION (LxW)	CASE QTY	CASE DIMENSIONS
41-4470	Setup Cover, 44x70" (112x178cm)	50	25x14x9" (63x35x23cm)
41-4470-S	Setup Cover, 44x70" (112x178cm)	40	25x14x9" (63x35x23cm)
41-4477	Setup Cover, 44x77" (112x196cm)	50	25x14x9" (63x35x23cm)
41-4477-S	Setup Cover, 44x77" (112x196cm)	40	25x14x9" (63x35x23cm)
41-4490	Setup Cover, 44x90" (112x229cm)	100	25x14x9" (63x35x23cm)
41-4490-S	Setup Cover, 44x90" (112x229cm)	25	25x14x9" (63x35x23cm)
41-5470	Setup Cover, 54x70" (137x178cm)	100	20x15x20" (52x39x50cm)
41-5470-S	Setup Cover, 54x70" (137x178cm)	25	20x15x20" (52x39x50cm)
41-5480	Setup Cover, 54x80" (137x203cm)	50	20x15x20" (52x39x50cm)
41-5480-S	Setup Cover, 54x80" (137x203cm)	40	20x15x20" (52x39x50cm)
41-5490	Setup Cover, 54x90" (137x229cm)	40	18x12x6" (47x30x15cm)
41-5490-S	Setup Cover, 54x90" (137x229cm)	15	16x12x16" (41x31x41cm)

### Materials

- Clear Polyethylene (PE) Film

### Product Attributes

- Clear film allows for visibility

### Sterilization

**Sterile:** Sterilized by EO (ethylene oxide). Product is sterile unless the package has been damaged or opened.  
**Non-Sterile:** If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packager. This includes further packaging and sterilization according to the validated processes of the re-packager. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

### Shelf Life

**Sterile:** 5 years  
**Non-Sterile:** Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

### Product Certifications/Standards


- Sterile: CE certified as Class I, according to MDD93/42/EEC Annex IX, Rule 1
- Non-Sterile: CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1

### Manufacturing Certifications/Registrations

- EN ISO 13485:2016
- Medical Device Single Audit Process Certified

### Storage

- General warehouse conditions



Alleset  
 4142 Industry Way  
 Flowery Branch, GA 30542,  
 USA  
[www.alleset.com](http://www.alleset.com)  
 SM-TDS-003 2023  
 Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.

GRI Medical & Electronic Technology Co., LTD  
 1805 Honggao Road,  
 314031 Jiaxing, Zhejiang  
 PEOPLE'S REPUBLIC OF CHINA  
[www.alleset.com](http://www.alleset.com)  
 US FDA Establishment Registration #3004911384

Emergo Europe B.V.  
 Westervoortsedijk 60,  
 6827 AT Arnhem,  
 The Netherlands  
[www.emergobyul.com](http://www.emergobyul.com)

