

# Surgical Accessories

## Surgical Skin Markers

CODE	DESCRIPTION (LxD)	CASE QTY	CASE DIMENSIONS
<span>N</span> 25-S001	Skin Marker, Standard Tip, 5.3in/135mm x 0.45in/11.2mm	2000	19x11x9" (49x27x22cm)
<span>S</span> 25-S001-S	Skin Marker, Standard Tip, 5.3in/135mm x 0.45in/11.2mm	100	13x11x9" (32x27x22cm)
<span>S</span> 25-S002-S	Skin Marker, Standard Tip w/ Ruler, 5.3in/135mm x 0.45in/11.2mm	100	12x11x8" (31x27x20cm)
<span>S</span> 25-S003-S	Skin Marker, Standard Tip w/ Ruler and Labels, 5.3in/135mm x 0.45in/11.2mm	100	13x11x9" (32x27x22cm)
<span>N</span> 25-F001	Skin Marker, Fine Tip, 5.3in/135mm x 0.45in/11.2mm	2000	19x11x9" (49x27x22cm)
<span>S</span> 25-F001-S	Skin Marker, Fine Tip, 5.3in/135mm x 0.45in/11.2mm	100	13x11x9" (32x27x22cm)
<span>S</span> 25-F002-S	Skin Marker, Fine Tip w/ Ruler, 5.25in/133mm x 0.45in/12.7mm	100	12x11x8" (31x27x20cm)
<span>S</span> 25-F003-S	Skin Marker, Fine Tip w/ Ruler and Labels, 5.3in/135mm x 0.45in/11.2mm	100	13x11x9" (32x27x22cm)
<span>N</span> 25-D001	Skin Marker, Dual Tip, 6.75in/171mm x 0.45in/11.2mm	1800	19x11x9" (49x27x22cm)
<span>S</span> 25-D001-S	Skin Marker, Dual Tip, 6.75in/171mm x 0.45in/11.2mm	100	13x11x9" (32x27x22cm)
<span>N</span> 25-D002	Skin Marker, Dual Tip w/ Ruler, 6.75in/171mm x 0.45in/11.2mm	500	14X12X12" (35x31x31cm)
<span>S</span> 25-D002-S	Skin Marker, Dual Tip w/ Ruler, 6.75in/171mm x 0.45in/11.2mm	100	12x11x8" (31x27x20cm)
<span>S</span> 25-D003-S	Skin Marker, Dual Tip w/ Ruler and Labels, 6.75in/171mm x 0.45in/11.2mm	100	13x11x9" (32x27x22cm)
<span>N</span> 25-C001	Skin Marker, Combo Tip (skin/utility), 6.75in/171mm x 0.45in/11.2mm	1800	19x11x9" (49x27x22cm)
<span>S</span> 25-C001-S	Skin Marker, Combo Tip (skin/utility), 6.75in/171mm x 0.45in/11.2mm	100	13x11x9" (32x27x22cm)
<span>S</span> 25-C002-S	Skin Marker, Combo Tip (skin/utility) w/ Ruler, 6.75in/171mm x 0.45in/11.2mm	100	13x11x9" (32x27x22cm)
<span>S</span> 25-C003-S	Skin Marker, Combo Tip (skin/utility) w/ Ruler Labels, 6.75in/171mm x 0.45in/11.2mm	100	13x11x9" (32x27x22cm)
<span>N</span> 25-FX001	Skin Marker, Fine Tip, XL, Prep-Resistant	2000	13x11x9" (32x27x22cm)
<span>S</span> 25-FX001-S	Skin Marker, Fine Tip, XL, Prep-Resistant	100	13x11x9" (32x27x22cm)
<span>S</span> 25-FX002-S	Skin Marker, Fine Tip, w/ Ruler, XL, Prep-Resistant	100	13x11x9" (32x27x22cm)
<span>S</span> 25-FX003-S	Skin Marker, Fine Tip w/ Ruler and Labels, XL, Prep-Resistant	100	13x11x9" (32x27x22cm)
<span>N</span> 25-SX001	Skin Marker, Standard Tip, XL, Prep-Resistant	2000	13x11x9" (32x27x22cm)
<span>S</span> 25-SX001-S	Skin Marker, Standard Tip, XL, Prep-Resistant	100	13x11x9" (32x27x22cm)
<span>S</span> 25-SX002-S	Skin Marker, Standard Tip w/ Ruler, XL, Prep-Resistant	100	13x11x9" (32x27x22cm)
<span>S</span> 25-SX003-S	Skin Marker, Standard Tip w/ Ruler and Labels, XL, Prep-Resistant	100	13x11x9" (32x27x22cm)
<span>N</span> 80226	Ruler, White, 6x0.5" (157x15mm)	1000	7x5x2" (19x13x7cm)

### Materials

- Barrels and Caps – Polypropylene (PP)
- Filler – Polyester
- Nib (Standard Tip) – Polyester Oriented Fiber
- Nib (Fine Tip) – Polypropylene (PP)
- Ink – Gentian Violet
- Ruler – Polyvinyl Chloride (PVC)

### Product Attributes

- Skin safe ink
- Available in standard tip or fine tip
- Color coded barrels

### Sterilization

Sterile: Sterilized by EO (ethylene oxide). Product is sterile unless the package has been damaged or opened.  
 Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packer. This includes further packaging and sterilization according to the validated processes of the re-packer. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

### Shelf Life

Sterile: 5 years  
 Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

### Product Certifications/Standards


- CE certified as Class I, according to MDR 2017/745, Annex VIII, Rule 1

### Storage


- General warehouse conditions

# Surgical Accessories


## Surgical Skin Markers



Fine/Standard Tip



Dual Tip Skin/Utility



Combo Tip Skin/Utility



Alleset  
4142 Industry Way  
Flowery Branch, GA 30542,  
USA

[www.alleset.com](http://www.alleset.com)

SM-TDS-0032023

Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD  
1805 Honggao Road,  
314031 Jiaxing, Zhejiang  
PEOPLE'S REPUBLIC OF CHINA  
[www.alleset.com](http://www.alleset.com)

US FDA Establishment Registration #3004911384



Emergo Europe B.V.  
Westervoortsedijk 60,  
6827 AT Arnhem,  
The Netherlands  
[www.emergobyul.com](http://www.emergobyul.com)

