

# Surgical Accessories

## Surgical Skin Markers

| CODE                      | DESCRIPTION (LxD)   | CASE QTY | CASE DIMENSIONS        |
|---------------------------|---|----------|------------------------|
| <span>N</span> 25-S001    | Skin Marker, Standard Tip, 5.3in/135mm x 0.45in/11.2mm                              | 2000     | 19x11x9" (49x27x22cm)  |
| <span>S</span> 25-S001-S  | Skin Marker, Standard Tip, 5.3in/135mm x 0.45in/11.2mm                              | 100      | 13x11x9" (32x27x22cm)  |
| <span>S</span> 25-S002-S  | Skin Marker, Standard Tip w/ Ruler, 5.3in/135mm x 0.45in/11.2mm                     | 100      | 12x11x8" (31x27x20cm)  |
| <span>S</span> 25-S003-S  | Skin Marker, Standard Tip w/ Ruler and Labels, 5.3in/135mm x 0.45in/11.2mm          | 100      | 13x11x9" (32x27x22cm)  |
| <span>N</span> 25-F001    | Skin Marker, Fine Tip, 5.3in/135mm x 0.45in/11.2mm                                  | 2000     | 19x11x9" (49x27x22cm)  |
| <span>S</span> 25-F001-S  | Skin Marker, Fine Tip, 5.3in/135mm x 0.45in/11.2mm                                  | 100      | 13x11x9" (32x27x22cm)  |
| <span>S</span> 25-F002-S  | Skin Marker, Fine Tip w/ Ruler, 5.25in/133mm x 0.45in/12.7mm                        | 100      | 12x11x8" (31x27x20cm)  |
| <span>S</span> 25-F003-S  | Skin Marker, Fine Tip w/ Ruler and Labels, 5.3in/135mm x 0.45in/11.2mm              | 100      | 13x11x9" (32x27x22cm)  |
| <span>N</span> 25-D001    | Skin Marker, Dual Tip, 6.75in/171mm x 0.45in/11.2mm                                 | 1800     | 19x11x9" (49x27x22cm)  |
| <span>S</span> 25-D001-S  | Skin Marker, Dual Tip, 6.75in/171mm x 0.45in/11.2mm                                 | 100      | 13x11x9" (32x27x22cm)  |
| <span>N</span> 25-D002    | Skin Marker, Dual Tip w/ Ruler, 6.75in/171mm x 0.45in/11.2mm                        | 500      | 14X12X12" (35x31x31cm) |
| <span>S</span> 25-D002-S  | Skin Marker, Dual Tip w/ Ruler, 6.75in/171mm x 0.45in/11.2mm                        | 100      | 12x11x8" (31x27x20cm)  |
| <span>S</span> 25-D003-S  | Skin Marker, Dual Tip w/ Ruler and Labels, 6.75in/171mm x 0.45in/11.2mm             | 100      | 13x11x9" (32x27x22cm)  |
| <span>N</span> 25-C001    | Skin Marker, Combo Tip (skin/utility), 6.75in/171mm x 0.45in/11.2mm                 | 1800     | 19x11x9" (49x27x22cm)  |
| <span>S</span> 25-C001-S  | Skin Marker, Combo Tip (skin/utility), 6.75in/171mm x 0.45in/11.2mm                 | 100      | 13x11x9" (32x27x22cm)  |
| <span>S</span> 25-C002-S  | Skin Marker, Combo Tip (skin/utility) w/ Ruler, 6.75in/171mm x 0.45in/11.2mm        | 100      | 13x11x9" (32x27x22cm)  |
| <span>S</span> 25-C003-S  | Skin Marker, Combo Tip (skin/utility) w/ Ruler Labels, 6.75in/171mm x 0.45in/11.2mm | 100      | 13x11x9" (32x27x22cm)  |
| <span>N</span> 25-FX001   | Skin Marker, Fine Tip, XL, Prep-Resistant   | 2000     | 13x11x9" (32x27x22cm)  |
| <span>S</span> 25-FX001-S | Skin Marker, Fine Tip, XL, Prep-Resistant   | 100      | 13x11x9" (32x27x22cm)  |
| <span>S</span> 25-FX002-S | Skin Marker, Fine Tip, w/ Ruler, XL, Prep-Resistant                                 | 100      | 13x11x9" (32x27x22cm)  |
| <span>S</span> 25-FX003-S | Skin Marker, Fine Tip w/ Ruler and Labels, XL, Prep-Resistant                       | 100      | 13x11x9" (32x27x22cm)  |
| <span>N</span> 25-SX001   | Skin Marker, Standard Tip, XL, Prep-Resistant                                       | 2000     | 13x11x9" (32x27x22cm)  |
| <span>S</span> 25-SX001-S | Skin Marker, Standard Tip, XL, Prep-Resistant                                       | 100      | 13x11x9" (32x27x22cm)  |
| <span>S</span> 25-SX002-S | Skin Marker, Standard Tip w/ Ruler, XL, Prep-Resistant                              | 100      | 13x11x9" (32x27x22cm)  |
| <span>S</span> 25-SX003-S | Skin Marker, Standard Tip w/ Ruler and Labels, XL, Prep-Resistant                   | 100      | 13x11x9" (32x27x22cm)  |
| <span>N</span> 80226      | Ruler, White, 6x0.5" (157x15mm)   | 1000     | 7x5x2" (19x13x7cm)     |

### Materials

- Barrels and Caps – Polypropylene (PP)
- Filler – Polyester
- Nib (Standard Tip) – Polyester Oriented Fiber
- Nib (Fine Tip) – Polypropylene (PP)
- Ink – Gentian Violet
- Ruler – Polyvinyl Chloride (PVC)

### Product Attributes

- Skin safe ink
- Available in standard tip or fine tip
- Color coded barrels

### Sterilization

**Sterile:** Sterilized by EO (ethylene oxide). Product is sterile unless the package has been damaged or opened.

**Non-Sterile:** If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packager. This includes further packaging and sterilization according to the validated processes of the re-packager. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

### Shelf Life

**Sterile:** 5 years

**Non-Sterile:** Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

### Product Certifications/Standards




- Sterile: CE certified as Class I, according to MDD93/42/EEC Annex IX, Rule 1
- Non-Sterile: CE certified as Class I, according to MDR 2017/745, Annex VIII, Rule 1

### Storage

- General warehouse conditions

# Surgical Accessories

## Surgical Skin Markers



Alleset  
4142 Industry Way  
Flowery Branch, GA 30542,  
USA  
www.alleset.com  
SM-TDS-003 2023  
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD  
1805 Honggao Road,  
314031 Jiaxing, Zhejiang  
PEOPLE'S REPUBLIC OF CHINA  
www.alleset.com  
US FDA Establishment Registration #3004911384

|    |     |
|----|-----|
| EC | REP |
|----|-----|

Emergo Europe B.V.  
Westervoortsedijk 60,  
6827 AT Arnhem,  
The Netherlands  
www.emergobyul.com

