

Surgical Accessories

Solidifier

CODE	DESCRIPTION	CASE QTY	CASE DIMENSIONS
70-S1000	Dissolvable Pouch, Gels 1000cc	100	28x18x7" (25x30x23cm)
70-S1500	Bottle, Gels 1500cc	96	11x9x5" (43x28x15cm)
70-S3000	Bottle, Gels 3000cc	96	(36x53x18cm)

Materials

- Super Absorbent Polymer (SAP)

Product Attributes

- Sized according to suction canister capacity
- No shaking required, reduces chance of spills or splashes
- Funnel style cap for delivery through suction ports when necessary

Sterilization

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packer. This includes further packaging and sterilization according to the validated processes of the re-packer. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

Shelf Life

Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

Product Certifications/Standards


- CE certified as Class I, according to MDR 2017/745 Annex IX, Rule 1

Manufacturing Certifications/Registrations

- EN ISO 13485:2016

Storage

- General warehouse conditions



Alleset
4142 Industry Way
Flowery Branch, GA 30542,
USA
www.alleset.com
SM-TDS-002 2022
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.

GRI Medical & Electronic Technology Co., LTD
1805 Honggao Road,
314031 Jiaxing, Zhejiang
PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
US FDA Establishment Registration #3004911384

Emergo Europe B.V.
Westervoortsedijk 60,
6827 AT Arnhem,
The Netherlands
www.emergobyul.com

