

Prepping Sponge Forceps

| CODE | DESCRIPTION | CASE QTY | CASE DIMENSIONS |
|---------|---|----------|-------------------------|
| 25-0001 | Forceps, White (Tweezer-Style), 5" (12.5cm) | 2000 | 24x15x16" (61x38x40cm) |
| 25-0106 | Forceps, Towel, Blue, 5" (13cm) | 500 | 10x16x16" (25x40x40cm) |
| 25-0540 | Forceps, Sponge, Green, 5.5" (14cm) | 500 | 19x13x7" (48x33x18.5cm) |
| 25-0550 | Forceps, Sponge, Green, 7.7" (19.5cm) | 500 | 12x13x20" (30x33x49cm) |
| 25-0562 | Forceps, Sponge, Green, 9" (23cm) | 400 | 14x12x18" (36x30x45cm) |

Materials

- Polypropylene (PP)
- Acrylonitrile Butadiene tyrene (ABS)
- Nylon

Product Attributes

- Features locking mechanism for clamping to desired materials

Sterilization

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packer. This includes further packaging and sterilization according to the validated processes of the re-packer. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

Shelf Life

Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

Product Certifications/Standards


- CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1

Manufacturing Certifications/Registrations


- EN ISO 13485:2016

Storage


- General warehouse conditions



25-0001



25-0106




25-0540, 25-0550, 25-0562



Alleset
4142 Industry Way
Flowery Branch, GA 30542,
USA
www.alleset.com
SM-TDS-003 2023
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD
1805 Honggao Road,
314031 Jiaxing, Zhejiang
PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
US FDA Establishment Registration #3004911384



Emergo Europe B.V.
Westervoortsedijk 60,
6827 AT Arnhem,
The Netherlands
www.emergobyul.com

