

# Prepping

## Sponge Forceps

CODE	DESCRIPTION	CASE QTY	CASE DIMENSIONS
25-0001	Forceps, White (Tweezer-Style), 5" (12.5cm)	2000	24x15x16" (61x38x40cm)
25-0106	Forceps, Towel, Blue, 5" (13cm)	500	10x16x16" (25x40x40cm)
25-0540	Forceps, Sponge, Green, 5.5" (14cm)	500	19x13x7" (48x33x18.5cm)
25-0550	Forceps, Sponge, Green, 7.7" (19.5cm)	500	12x13x20" (30x33x49cm)
25-0562	Forceps, Sponge, Green, 9" (23cm)	400	14x12x18" (36x30x45cm)

### Materials

- Polypropylene (PP)
- Acrylonitrile Butadiene tyrene (ABS)
- Nylon

### Product Attributes

- Features locking mechanism for clamping to desired materials

### Sterilization

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packager. This includes further packaging and sterilization according to the validated processes of the re-packager. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

### Shelf Life


Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

### Product Certifications/Standards


- Non-Sterile: CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1

### Manufacturing Certifications/Registrations


- EN ISO 13485:2016



25-001



25-0106



25-0540,  
25-0550,  
25-0562



Alleset  
4142 Industry Way  
Flowery Branch, GA 30542,  
USA  
www.alleset.com  
SM-TDS-003 2023  
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD  
1805 Honggao Road,  
314031 Jiaxing, Zhejiang  
PEOPLE'S REPUBLIC OF CHINA  
www.alleset.com  
US FDA Establishment Registration #3004911384



Emergo Europe B.V.  
Westervoortsedijk 60,  
6827 AT Arnhem,  
The Netherlands  
www.emergobyul.com

