

Prepping Sponge Sticks

CODE	DESCRIPTION	CASE QTY	CASE DIMENSIONS
25-2008	Sponge Stick, Blue, 6" (15cm)	1000	26x16x14" (67x40x35cm)
25-2009	Sponge Stick, Orange, 8" (20cm)	1000	31x20x14" (78x52x35cm)
25-2010	Sponge Stick, Rounded Top, Blue, 6" (15cm)	1000	26x16x14" (67x40x35cm)
25-2011	Sponge Stick, Angled, Blue, 6.3" (16cm)	1000	26x16x14" (67x40x35cm)
25-2012	Sponge Stick, Rounded Top, Blue, 6.3" (16cm)	1000	26x16x14" (67x40x35cm)
25-2020	Sponge Stick, Vaginal, Orange, 8.6" (22cm)	800	31x20x14" (78x52x35cm)

Materials

- Sponge: Polyurethane (PU) Foam
- Sponge: Polyurethane (PU) Foam
- Stick: Polystyrene

Product Attributes

- Ultrasonic welding is used to adhere foam to the stick

Sterilization

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packager. This includes further packaging and sterilization according to the validated processes of the re-packager. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

Shelf Life

Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

Product Certifications/Standards


- Non-Sterile: CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1

Manufacturing Certifications/Registrations


- EN ISO 13485:2016

Storage


- General warehouse conditions




25-2008




25-2009




25-2010



25-2011



25-2012



25-2020



Alleset
4142 Industry Way
Flowery Branch, GA 30542,
USA
www.alleset.com
SM-TDS-003 2023
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.

GRI Medical & Electronic Technology Co., LTD
1805 Honggao Road,
314031 Jiaxing, Zhejiang
PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
US FDA Establishment Registration #3004911384

Emurgo Europe B.V.
Westervoortsedijk 60,
6827 AT Arnhem,
The Netherlands
www.emurgobyul.com