

Orthopedic Dressings

Impervious Stockinette

CODE	DESCRIPTION	Size	CASE QTY	CASE DIMENSIONS
N 85-07340	Impervious Stockinette, w/o Pull Tab, 7x34" (18x86cm)	S	105	20x17x9" (53x43x22cm)
S 85-07340-S	Impervious Stockinette, w/o Pull Tab, 7x34" (18x86cm)	S	30	18x13x9" (45x34x24cm)
N 85-07341	Impervious Stockinette, w/ Pull Tab, 7x34" (18x86cm)	S	84	20x17x9" (53x43x22cm)
S 85-07341-S	Impervious Stockinette, w/ Pull Tab, 7x34" (18x86cm)	S	30	18x13x9" (45x34x24cm)
N 85-07342	Impervious Stockinette, w/o Pull Tab, 7x34" (18x86cm)	S	84	20x17x9" (53x43x22cm)
N 85-09480	Impervious Stockinette, w/o Pull Tab, 9x48" (23x122cm)	M	68	20x17x9" (53x43x22cm)
S 85-09480-S	Impervious Stockinette, w/o Pull Tab, 9x48" (23x122cm)	M	20	23x9x12" (60x22x30cm)
N 85-09481	Impervious Stockinette, w/ Pull Tab, 9x48" (23x122cm)	M	40	23x9x12" (60x22x30cm)
S 85-09481-S	Impervious Stockinette, w/ Pull Tab, 9x48" (23x122cm)	M	36	20x17x13" (50x42x32cm)
N 85-09482	Impervious Stockinette, w/o Pull Tab, 9x48" (23x122cm)	M	70	20x17x13" (50x42x32cm)
S 85-09482-S	Impervious Stockinette, w/o Pull Tab, 9x48" (23x122cm)	M	40	23x9x12" (60x22x30cm)
N 85-09300	Impervious Stockinette, w/o Pull Tab, 9x39" (23x75cm)	M	70	20x17x13" (50x42x32cm)
S 85-09300-S	Impervious Stockinette, w/o Pull Tab, 9x39" (23x75cm)	M	20	14x11x8" (36x29x21cm)
N 85-12480	Impervious Stockinette, w/o Pull Tab, 12x48" (30x122cm)	L	44	17x11x18" (42x28x45cm)
S 85-12480-S	Impervious Stockinette, w/o Pull Tab, 12x48" (30x122cm)	L	44	23x15x15" (59x37x38cm)
N 85-12481	Impervious Stockinette, w/ Pull Tab, 12x48" (30x122cm)	L	44	17x11x18" (42x28x45cm)
S 85-12481-S	Impervious Stockinette, w/ Pull Tab, 12x48" (30x122cm)	L	40	23x14x17" (58x36x43cm)
N 85-12482	Impervious Stockinette, w/o Pull Tab, 12x48" (30x122cm)	L	44	17x11x18" (42x28x45cm)
N 85-16480	Impervious Stockinette, w/o Pull Tab, 16x48" (41x122cm)	XL	27	18x13x9" (45x34x24cm)
N 85-16481	Impervious Stockinette, w/ Pull Tab, 16x48" (41x122cm)	XL	27	14x14x11" (36x36x27cm)
S 85-16481-S	Impervious Stockinette, w/ Pull Tab, 16x48" (41x122cm)	XL	40	23x14x17" (58x36x43cm)
N 85-16482	Impervious Stockinette, w/o Pull Tab, 16x48" (41x122cm)	XL	27	14x14x11" (36x36x27cm)

Materials

- Inner Layer: Polyester, White
- Impervious Outer Layer: Blue Elastic Polymer
- Pull Tab: Tyvek® High-Density Polyethylene Fiber Strips, White

Product Attributes

- Protective, impervious outer layer
- Woven polyester inner layer for comfort

Sterilization

Sterile: Sterilized by EO (ethylene oxide). Product is sterile unless the package has been damaged or opened.

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packager. This includes further packaging and sterilization according to the validated processes of the re-packager. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

Shelf Life

Sterile: 5 years

Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

Storage

- General warehouse conditions




Product Certifications/Standards

- Sterile: CE certified as Class I, according to MDD93/42/EEC Annex IX, Rule 4
- Non-Sterile: CE certified as Class I, according to MDR 2017/745, Annex VIII, Rule 1

Manufacturing Certifications/ Registrations

- EN ISO 13485:2016



Alleset
4142 Industry Way
Flowery Branch, GA 30542,
USA
www.alleset.com
SM-TDS-003 2023
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD
1805 Honggao Road,
314031 Jiaxing, Zhejiang
PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
US FDA Establishment Registration #3004911384



Emergo Europe B.V.
Westervoortsedijk 60,
6827 AT Arnhem,
The Netherlands
www.emergobyul.com