

























Surgical Accessories

Suction Tubings

CODE	DESCRIPTION	CASE QTY	CASE DIMENSIONS
 70-0506C	Suction Tubing, 3/16inx6ft (5mmx1.8m)	25	14x14x15" (36x36x39cm)
 70-0506C-S	Suction Tubing, 3/16inx6ft (5mmx1.8m)	50	11x11x16" (29x28x41cm)
 70-0510C	Suction Tubing, 3/16inx10ft (5mmx3.0m)	25	14x14x15" (36x36x39cm)
 70-0510C-S	Suction Tubing, 3/16inx10ft (5mmx3.0m)	50	11x11x16" (29x28x41cm)
 70-0512C	Suction Tubing, 3/16inx12ft (5mmx3.7m)	25	14x14x15" (36x36x39cm)
 70-0512C-S	Suction Tubing, 3/16inx12ft (5mmx3.7m)	50	11x11x16" (29x28x41cm)
 70-0520C	Suction Tubing, 3/16inx20ft (5mmx6.1m)	25	14x16x15" (35x40x37cm)
 70-0520C-S	Suction Tubing, 3/16inx20ft (5mmx6.1m)	50	13x11x18" (34x28x46cm)
 70-0606C	Suction Tubing, 1/4inx6ft (6mmx1.8m)	25	14x16x15" (35x40x37cm)
 70-0606C-S	Suction Tubing, 1/4inx6ft (6mmx1.8m)	50	11x11x16" (29x28x41cm)
 70-0610C	Suction Tubing, 1/4inx10ft (6mmx3.0m)	25	14x16x15" (35x40x37cm)
 70-0610C-S	Suction Tubing, 1/4inx10ft (6mmx3.0m)	50	11x11x16" (29x28x41cm)
 70-0612C	Suction Tubing, 1/4inx12ft (6mmx3.7m)	25	14x16x15" (35x40x37cm)
 70-0612C-S	Suction Tubing, 1/4inx12ft (6mmx3.7m)	50	11x11x16" (29x28x41cm)
 70-0620C	Suction Tubing, 1/4inx20ft (6mmx6.1m)	25	14x16x15" (35x40x37cm)
 70-0620C-S	Suction Tubing, 1/4inx20ft (6mmx6.1m)	50	13x11x18" (34x28x46cm)
 70-0706C	Suction Tubing, 9/32inx6ft (7mmx1.8m)	25	16x19x17" (40x48x42cm)
 70-0706C-S	Suction Tubing, 9/32inx6ft (7mmx1.8m)	50	11x11x16" (29x28x41cm)
 70-0710C	Suction Tubing, 9/32inx10ft (7mmx3.0m)	25	16x19x17" (40x48x42cm)
 70-0710C-S	Suction Tubing, 9/32inx10ft (7mmx3.0m)	50	11x11x16" (29x28x41cm)
 70-0712C	Suction Tubing, 9/32inx12ft (7mmx3.7m)	25	16x19x17" (40x48x42cm)
 70-0712C-S	Suction Tubing, 9/32inx12ft (7mmx3.7m)	50	13x21x17" 34x54x44cm
 70-0720C	Suction Tubing, 9/32inx20ft (7mmx6.1m)	25	13x20x21" 34x52x54cm
 70-0720C-S	Suction Tubing, 9/32inx20ft (7mmx6.1m)	50	13x21x17" 34x54x44cm

Materials

- Polyvinyl Chloride (PVC)

Product Attributes

- Low coil memory and greater flexibility
- Increase wall thickness helps prevent tubing collapse
- Color-coded female fittings for easy size identification

Sterilization

Sterile: Sterilized by EO (ethylene oxide). Product is sterile unless the package has been damaged or opened.

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packager. This includes further packaging and sterilization according to the validated processes of the re-packager. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

Shelf Life

Sterile: 5 years

Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions.


After 5 years, the product should be performance tested to ensure compliance to specifications.

Product Certifications/Standards

- Sterile: CE certified as Class I, according to MDD93/42/EEC Annex IX, Rule 1
- Non-Sterile: CE certified as Class I, according to MDR 2017/745, Annex VIII, Rule 1

Surgical Accessories

Suction Tubings



Alletes
4142 Industry Way
Flowery Branch, GA 30542,
USA
www.alletes.com
SM-TDS-0032023

Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD
1805 Honggao Road,
314031 Jiaxing, Zhejiang
PEOPLE'S REPUBLIC OF CHINA
www.alletes.com
US FDA Establishment Registration #3004911384



Emergo Europe B.V.
Westervoortsedijk 60,
6827 AT Arnhem,
The Netherlands.
www.emergobyul.com

