

Surgical Gowns

Film-Reinforced

CODE	DESCRIPTION	SIZE	CASE QTY	CASE DIMENSIONS
V90-2210	Surgical Gown, Film-Reinforced	M	44	18x15x12 (46x38x29.5cm)
V90-2210-S2	Surgical Gown, Film-Reinforced, 2 Hand Towels	M	26	18x12x19" (46x30x49cm)
V90-2220	Surgical Gown, Film-Reinforced	L	44	19x16x12.5" (50x42x32cm)
V90-2220-S2	Surgical Gown, Film-Reinforced, 2 Hand Towels	L	26	18x12x19" (46x30x49cm)
V90-2221	Surgical Gown, Film-Reinforced	L, X-Long	44	19x16x12.5" (50x42x32cm)
V90-2221-S2	Surgical Gown, Film-Reinforced, 2 Hand Towels	L, X-Long	26	18x12x19" (46x30x49cm)
V90-2230	Surgical Gown, Film-Reinforced	XL	44	19x16x12.5" (50x42x32cm)
V90-2230-S2	Surgical Gown, Film-Reinforced, 2 Hand Towels	XL	24	18x12x19" (46x30x49cm)
V90-2231	Surgical Gown, Film-Reinforced	XL, X-Long	44	19x16x12.5" (50x42x32cm)
V90-2231-S2	Surgical Gown, Film-Reinforced, 2 Hand Towels	XL, X-Long	24	18x12x19" (46x30x49cm)
V90-2240	Surgical Gown, Film-Reinforced	2XL	36	19x16x12.5" (50x42x32cm)
V90-2240-S2	Surgical Gown, Film-Reinforced, 2 Hand Towels	2XL	24	18x12x19" (46x30x49cm)

Materials

- Multilayered SMS Material
- Non-Woven Polypropylene (PP) and Polyethylene (PE) Film

Product Attributes

- Adjustable hook and loop fenestration
- Knit cuffs
- Set-in sleeves
- Sonically sealed on all construction points

Sterilization

Sterile: Sterilized by EO (ethylene oxide). Product is sterile unless the package has been damaged or opened.
Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packager. This includes further packaging and sterilization according to the validated processes of the re-packager. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

Shelf Life

Sterile: 5 years
Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

Product Certifications/Standards


- ANSI/AAMI PB70:2012 Level 3 Liquid Barrier Protection
- EN 13795:2011 High Performance
- Sterile: CE certified as Class I, according to MDD93/42/EEC Annex IX, Rule 4
- Non-Sterile: CE certified as Class I, according to MDR 2017/745, Annex VIII, Rule 1

Manufacturing Certifications/ Registrations

- EN ISO 13485:2016

Storage

- General warehouse conditions



CODE	A: Back Length	B: Sleeve Length	C: Body Reinf Top Width	D: Body Reinf Bottom Width	E: Body Reinf Height
V90-2210 (S2/NS)	41" (104cm)	23" (58cm)	12" (30cm)	17" (43cm)	22" (56cm)
V90-2220 (S2/NS)	45" (115cm)	24" (62cm)	13" (33cm)	17" (43cm)	25" (64cm)
V90-2221 (S2/NS)	54" (138cm)	24" (62cm)	13" (33cm)	17" (43cm)	25" (64cm)
V90-2230 (S2/NS)	49" (124cm)	25" (64cm)	13" (33cm)	17" (43cm)	28" (71cm)
V90-2231 (S2/NS)	58" (147cm)	25" (64cm)	13" (33cm)	17" (43cm)	28" (71cm)
V90-2240 (S2/NS)	51" (130cm)	30" (75cm)	16" (41cm)	21" (53cm)	30" (76cm)
Tolerance	±0.8" (2cm)	±0.4" (1cm)	±0.4" (1cm)	±0.4" (1cm)	±0.8" (2cm)



Alleset
 4142 Industry Way
 Flowery Branch, GA 30542,
 USA
 www.alleset.com
 SM-TDS-002 2022
 Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.

GRI Medical & Electronic Technology Co., LTD
 1805 Honggao Road,
 314031 Jiaxing, Zhejiang
 PEOPLE'S REPUBLIC OF CHINA
 www.alleset.com
 US FDA Establishment Registration #3004911384

Emergo Europe B.V.
 Prinsessegracht 20
 2514 AP The Hague
 The Netherlands
 www.emergobyul.com

