


Surgical Drapes

Ophthalmic

CODE	DESCRIPTION	CASE QTY	CASE DIMENSIONS
 80-04406-S	Ophthalmic Drape, w/ Full Incise and Fluid Collection Pouch, 56x52" (142x132cm)	50	23x14x17" (58x36x43cm)

Materials

- Multilayered SMS Material with Antistatic Properties and Alcohol Repellency
- Fluid Collection Pouch, Clear Polyethylene (PE) Film

Product Attributes

- Non-linting
- Eye wick
- Fluid collection pouch
- Adhesive incise area, non fenestrated

Sterilization

Sterile: Sterilized by EO (ethylene oxide). Product is sterile unless the package has been damaged or opened.

Shelf Life

Sterile: 5 years

Product Certifications/Standards


- CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1

Manufacturing Certifications/Registrations

- EN ISO 13485:2016
- Medical Device Single Audit Process Certified

Storage

- General warehouse conditions



Alleset
4142 Industry Way
Flowery Branch, GA 30542,
USA

www.alleset.com

SM-TDS-003 2023

Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.

THINK.MAKE.PROTECT.



GRI Medical & Electronic Technology Co., LTD
1805 Honggao Road,
314031 Jiaxing, Zhejiang
PEOPLE'S REPUBLIC OF CHINA
www.alleset.com

US FDA Establishment Registration #3004911384



Emergo Europe B.V.
Westervoortsedijk 60,
6827 AT Arnhem,
The Netherlands
www.emergobyul.com

A GRI COMPANY

