

# Surgical Accessories

## Towels

CODE	DESCRIPTION	CASE QTY	CASE DIMENSIONS
25-0005	OR Towel, Cotton, Blue, 17x24" (43x61cm)	400	23x17x11" (58x44x29cm)
25-0006	OR Towel, Cotton, White, 17x24" (43x61cm)	400	23x17x11" (58x44x29cm)
25-0007	OR Towel, Cotton, Green, 17x24" (43x61cm)	400	23x17x11" (58x44x29cm)
25-0008	OR Towel, Cotton, Fenestrated, 17x24" (43x61cm)	400	23x17x11" (58x44x29cm)

### Materials

Cotton

### Product Attributes

- Highly absorbent and durable
- Color fast

### Sterilization

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packer. This includes further packaging and sterilization according to the validated processes of the re-packer. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

### Shelf Life

Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

### Product Certifications/Standards


- Non-Sterile: CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1

### Manufacturing Certifications/Registrations


- EN ISO 13485:2016

### Storage

- General warehouse conditions



25-005, 25-006, 25-007



25-0008



Alleset  
4142 Industry Way  
Flowery Branch, GA 30542,  
USA  
www.alleset.com  
SM-TDS-002 2023  
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD  
1805 Honggao Road,  
314031 Jiaxing, Zhejiang  
PEOPLE'S REPUBLIC OF CHINA  
www.alleset.com  
US FDA Establishment Registration #3004911384



Emergo Europe B.V.  
Westervoortsedijk 60,  
6827 AT Arnhem,  
The Netherlands  
www.emergobyul.com

