

WasteStop™

Waste Containment Bowl


	CODE	DESCRIPTION	CASE QTY	CASE DIMENSIONS
	19-0500-BL	WasteStop Waste Containment Basin, Blue, Assembled	20	TBD
	19-0500-BL-S	WasteStop Waste Containment Basin, Blue, Assembled	20	TBD
	19-0500B-BL	WasteStop Waste Containment Bowl, Blue, Bulk	200	TBD

Materials

- Basin: Polypropylene (PP)
- Lid: Polyvinyl Chloride (PVC)
 - Adhesive

Product Attributes

- Conical insertion point directs fluid into the bowl
- Biohazard adhesive tape closure
- SAP pad absorbs liquid contents



Sterilization

Sterile: Sterilized by EO (ethylene oxide). Product is sterile unless the package has been damaged or opened.

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packager. This includes further packaging and sterilization according to the validated processes of the re-packager. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

Shelf Life

Sterile: 5 years

Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

Product Certifications/Standards

- Sterile: CE certified as Class I, according to MDD93/42/EEC Annex IX, Rule 1
- Non-Sterile: CE certified as Class I, according to MDR 2017/745, Annex VIII, Rule 1

Manufacturing Certifications/Registrations

- EN ISO 13485:2016
- Medical Device Single Audit Process Certified

Storage


- General warehouse conditions



Alleset
 4142 Industry Way
 Flowery Branch, GA 30542,
 USA
www.alleset.com
 SM-TDS-003 2023
 Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD
 1805 Honggao Road,
 314031 Jiaxing, Zhejiang
 PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
 US FDA Establishment Registration #3004911384



Emergo Europe B.V.
 Westervoortsedijk 60,
 6827 AT Arnhem,
 The Netherlands
www.emergobyul.com

