

# Surgical Accessories

## Yankauer Suction Tips

CODE	DESCRIPTION	CASE QTY	CASE DIMENSIONS
70-50130	Bulb Tip, Non-Vented, 0.42x10.78" (1.07x27.4cm)	130	12x6x10" (30x15x26cm)
70-50130-S	Bulb Tip, Non-Vented, 0.42x10.78" (1.07x27.4cm)	50	12x6x10" (30x15x26cm)
70-50132	Bulb Tip, Non-Vented, 0.42x10.78" (1.07x27.4cm)	130	12x6x10" (30x15x26cm)
70-50132-S	Bulb Tip, Non-Vented, 0.42x10.78" (1.07x27.4cm)	50	12x6x10" (30x15x26cm)
70-50140	Straight Tip, Vented, 0.25x9.78" (0.65x24.85cm)	130	12x6x10" (30x15x26cm)
70-50140-S	Straight Tip, Vented, 0.25x9.78" (0.65x24.85cm)	50	12x6x10" (30x15x26cm)
70-50142	Straight Tip, Vented, 0.25x9.78" (0.65x24.85cm)	130	12x6x10" (30x15x26cm)
70-50142-S	Straight Tip, Vented, 0.25x9.78" (0.65x24.85cm)	50	12x6x10" (30x15x26cm)

### Materials

- K-Resin

### Product Attributes

- Clear material for ideal visualization during suction
- Optional vents allow for controlled suction

### Sterilization

**Sterile:** Sterilized by EO (ethylene oxide). Product is sterile unless the package has been damaged or opened.

**Non-Sterile:** If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packager. This includes further packaging and sterilization according to the validated processes of the re-packager. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

### Shelf Life

**Sterile:** 5 years

**Non-Sterile:** Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

### Product Certifications/Standards


- Sterile: CE certificated as Class I, according to MDD93/42/EEC Annex IX, Rule 1
- Non-Sterile: CE certified as Class I, according to MDR 2017/745, Annex VIII, Rule 1

### Manufacturing Certifications/Registrations

- EN ISO 13485:2016

### Storage

- General warehouse conditions



Yankauer Suction, Straight Tip

Yankauer Suction, Bulb Tip



Alleset  
4142 Industry Way  
Flowery Branch, GA 30542,  
USA  
www.alleset.com  
SM-TDS-003 2023  
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.

GRI Medical & Electronic Technology Co., LTD  
1805 Honggao Road,  
314031 Jiaxing, Zhejiang  
PEOPLE'S REPUBLIC OF CHINA  
www.alleset.com  
US FDA Establishment Registration #3004911384

Emergo Europe B.V.  
Westervoortsedijk 60,  
6827AT Arnhem,  
The Netherlands.  
www.emergobyul.com

