


EYE SHIELDS

High-Performance



Alleset Performance Eye Shields are a stylish way to provide complete eye protection during medical and surgical procedures. These protective eye shields offer sleek, updated styling that are comfortable and fit over prescription glasses, thus promoting compliance.



The Performance Eye Shields have a special coating to prevent fogging during use, providing for a clear view throughout the procedure.

The lens attachment holes are completely occluded by the frame, offering complete protection from splashes and sprays.

- | | | |
|----------|------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------|
| 1 | Fits Over Prescription Glasses | Unique design provides a comfortable and secure fit over prescription glasses so special eye wear is not needed. |
| 2 | Colored Strips on Lens Top Edge Available in Four Stylish Colors | Allows for finger-proof handling. |
| 3 | Non-Perforated Lens | Helps prevent fluids from passing through lens attachment holes. |
| 4 | Protective Film Covering | Protective film release liner is easily removed with pull tab and prevents scratching. |
| 5 | Innovative Extra Wide Inner Frame Ridge | Fills the gap at the temples, providing better top-down fluid and splash protection. |
| 6 | Anti-Fog, Anti-Glare & Anti-Static Clear Lens | Inhibits fog from forming on the lens, provide a glare-free view, and anti-static properties for safety in the operating theater. |
| 7 | Wrap Around Frame | Provides for a more secure fit, preventing slipping, and comfortable placement behind the ears. |
| 8 | Reusable Frame | Easily replaceable lenses promote reuse of frame helping to be more cost effective. |
| | Does not contain natural rubber latex or phthalates | Safe for patients and staff with possible allergies and prevention of leaching of potentially harmful compounds. |

CODE	COLOR	DESCRIPTION	PKG
33-0500	Multicolor	Eye Shields, Pre-Assembled, Clear	200
33-0501	Multicolor	Eye Shields, Lens, Clear	2500
33-0502	Black	Eye Shields, Frames, Polypropylene	200