



SKIN MARKER

Leave Your Mark, Safely



Alleset offers a growing range of skin markers for use before and during surgery. Utilizing a special water-based formulation, the Gentian-violet ink is safe for use on the skin. It is non-toxic and will not tattoo like some other markers may do.

The dual-tip marker provides both a fine tip and a standard tip in one convenient skin marker. This dual-tip marker is filled with the same skin-safe Gentian-violet ink. GRI-Alleset also offers a combo-tip marker, with one end safe for use on skin, while the other "utility" end is formulated with a solvent-based permanent ink for marking plastics, metals, or other items.



- | | | |
|----------|---|---|
| 1 | Skin-Safe Ink | Specially formulated ink is non-toxic and will not tattoo |
| 2 | Available in Standard Tip or Fine Tip | Preferences are met, largely depending on the type of procedure (e.g., plastics versus general use) |
| 3 | Color-Coded Barrels | Blue for standard tip and white for fine tip makes differentiation simple and clear |
| 4 | Available Utility Ink | Dual tip marker includes permanent ink for writing on plastic, metal, or other objects |
| 5 | Available Prep-Resistant Ink | Prep-Resistant ink markers are available that better withstand common alcohol-based solutions such as Chloraprep or Duraprep™ |
| | Does not contain natural rubber latex or phthalates | Safe for patients and staff with possible allergies and prevention of leaching of potentially harmful compounds |

STERILE	PKG	NON-STERILE	PKG	DESCRIPTION
25-C001-S	100/CS	25-C001	1800/CS	Combo Tip*
25-C002-S	100/CS			Combo Tip* with Ruler
25-C003-S	100/CS			Combo Tip* with Ruler and Labels
25-S001-S	100/CS	25-S001	2000/CS	Standard Tip
25-S002-S	100/CS			Standard Tip with Ruler
25-S003-S	100/CS			Standard Tip with Ruler and Labels
25-F001-S	100/CS	25-F001	2000/CS	Fine Tip
25-F002-S	100/CS			Fine Tip with Ruler
25-F003-S	100/CS			Fine Tip with Ruler and Labels
25-D001-S	100/CS	25-D001	1800/CS	Dual Tip
25-D002-S	100/CS			Dual Tip with Ruler
25-D003-S	100/CS			Dual Tip with Ruler and Labels
25-SX001-S	100/CS	25-SX001	2000/CS	Standard Tip, "XL" Prep-Resistant
25-SX002-S	100/CS			Standard Tip with Ruler, "XL" Prep-Resistant
25-SX003-S	100/CS			Standard Tip with Ruler and Labels, "XL" Prep-Resistant
25-FX001-S	100/CS	25-FX001	2000/CS	Fine Tip, "XL" Prep-Resistant
25-FX002-S	100/CS			Fine Tip with Ruler, "XL" Prep-Resistant
25-FX003-S	100/CS			Fine Tip with Ruler and Labels, "XL" Prep-Resistant

*Combo tip has skin ink in one end and utility ink in the other end.