

ROOM TURNAROUND KITS

Cleanup, Preparation, Protection

Introducing **CORE (Clean Operating Room Environment) Room Turnaround Kits®**, a new and more efficient way to change your Operating Theatre environment between procedures, customised to the unique needs of your hospital.



BETTER SYSTEMS. BETTER EFFICIENCY. BETTER PROTECTION.

Reducing the risk of hospital-acquired, or nosocomial infections, is a major focus for today's healthcare professionals. With the increase of multi drug resistant organisms, minimising the risk of transmission is critical, particularly in an environment where infections can be spread by various means; via the healthcare staff, contaminated equipment, bed linens, air droplets, or fluids.

Choosing the right products can help to lessen the chances of hospital-acquired infections, improving patient outcomes and keeping patients, staff, and equipment protected. Our kits and components can play an important role in the reduction of such risks.

Better Systems

At Alleset, we are focused on helping hospitals achieve their goals through our comprehensive line of disposable operating room turnaround kits. Each component is specifically designed to minimise cross-contamination and work with other kit components to improve room cleaning protection and efficiency.

Hospitals can purchase these non-sterile disposable items as individual components or as a fully assembled kit, each containing the products necessary to prepare the room for the next surgical procedure. Our comprehensive portfolio includes a variety of impervious absorbent operating room table covers, arm board covers, headrest covers, pillow covers, patient safety straps, waste/liner bags, and mops. Kits can easily be tailored to every surgical department's needs.

Waste/Liner Bags



Mop Heads and Kick Bucket Liner



Disposable, Protective OR Linens



Standard or Custom Kits



Better Efficiency

Alleset's room turnaround kits enable staff to turn and prepare the operating theatre for every procedure with maximum efficiency, helping reduce set up time, improve patient and staff experience, and maintain a consistent approach to cleaning practices. Each customisable kit can be opened as a unit, rather than as individual pieces, making preparation quick and easy, and allowing staff to spend more time with patients while potentially increasing the volume of surgical procedures conducted on a daily or weekly basis.


Better Protection

The need for protection encompasses the patient, staff, and operating theatre equipment, with each component working together to help reduce cross-contamination, improve patient outcomes and help maintain a safe work environment.

Fluid Management

Providing a dry environment for patients is important for maintaining skin integrity and reducing the potential for breakdown and ulcers. Our super absorbent table covers are designed to quickly wick fluids away from the patient, trapping them in the absorbent core and away from the dry, soft, and smooth outer layer.


- Protection from blood and other body fluids is increased with our absorbent table sheet, which moves excess fluid from the surface, trapping it in the middle core away from the staff and equipment.
- Reducing exposure to fluids helps to extend the life of the operating table and makes clean-up easier.
- Our impermeable bottom layer prevents strike through, lowering the risk of cross-contamination.
- The impermeable bottom layer of the Premium Absorbent Table Cover was designed with improved grip to help reduce shifting and bunching.



Cross-Contamination

When it comes to cleaning the operating theatre mattress, even the smallest tear, cut, or hole can make disinfection difficult, potentially leaving harmful microorganisms, which can be transferred through contact. An impervious, disposable table cover provides a protective barrier between equipment, patient, and staff.


Cross-contamination can also occur through stitch holes, which are typically found in arm board covers, pillow cases, and headrest covers. Alleset is the first to incorporate sonic bonding into our linens, eliminating holes which can provide a pathway for microbial contamination.



Lint Reduction

Another important consideration in minimising the potential for post-operative complications is to decrease the possibility of a foreign body reaction. Lint can be shed from many items including woven reusable table covers, sheets, and pillow cases, which continue to be used during surgical procedures. This lint can act as a vehicle to carry unwanted particulates.

The more lint generated, the greater the possibility of post-operative complications caused by micro-organism transfer or by a foreign body. Our single patient use products, made with non-woven fabrics and film, are designed for minimal lint generation, helping to reduce this risk.



Better Results, Done More Effectively

Reducing hospital-acquired infections leads to better outcomes for patients, for staff, and for hospitals and surgical facilities. Alleset's CORE Room Turnaround Kits® provide a simple yet effective solution for ensuring those outcomes are achieved in advance of every procedure.

AVAILABLE KITS & COMPONENTS

Following is a list of components available to build your room turnaround kit. Each kit can be customised to suit your needs.

DISPOSABLE LINENS		SIZE
50-HRC001	Head Rest Cover, Blue	14 x 15 in / 36 x 38 cm
50-IPC001	Impervious Pillow Cover, Blue	19.5 x 27 in / 50 x 69 cm
50-LFT001	Lift/Draw Sheet, White	34 x 66 in / 86 x 168 cm
50-LFT002	Lift/Draw Sheet, XL, White	40 x 92 in / 102 x 234 cm
50-LFT003	Lift/Draw Sheet, L, White	40 x 80 in / 102 x 203 cm
99-UPIMP4090	Absorbent Table Cover, Blue	40 x 90 in / 102 x 229 cm
99-UPIMP4060	Absorbent Table Cover, Blue	40 x 60 in / 102 x 152 cm
99-UPIMP4090G	Premium Absorbent Table Cover	40 x 90 in / 102 x 229 cm
99-UPIMP4060G	Premium Absorbent Table Cover	40 x 60 in / 102 x 152 cm
50-STRPW1	Patient Safety Strap, White	43 x 3 in / 109 x 8 cm
50-AMBC01	Armboard Cover, Blue	13 x 30 in / 33 x 76 cm
50-AMBC02	Armboard Cover w/2 Straps, Blue	13 x 30 in / 33 x 76 cm
50-TOWL01	Disposable Towel, White	12 x 15 in / 30 x 38 cm
50-TOWL02	Disposable Towel, Microfiber, White	12 x 12 in / 30 x 30 cm

KICK BUCKET BAGS		SIZE
50-KBBR01	Kick Bucket Bag, Red	23 x 24 in / 58 x 61 cm
50-KBBR02	Kick Bucket Bag, Red	24 x 31 in / 61 x 79 cm
50-KBBC01	Kick Bucket Bag, Clear	16 x 12 x 24 in / 41 x 30 x 61 cm
50-KBBB01	Kick Bucket Bag, Black	40 x 46 in / 102 x 117 cm

MOP HEADS		SIZE
50-MOPR12	Mop Head, White	12 oz / 340 g
50-MOPF12	Mop Head, Fantail, White	12 oz / 340 g
50-MOPR16	Mop Head, White	16 oz / 454 g
50-MOPF16	Mop Head, Fantail, White	16 oz / 454 g
50-MOPC01	Mop Head, Non-Cellulose, White	16 oz / 454 g
50-MOPM01	Mop, Microfiber, White	5.75 x 16 in / 15 x 40 cm

DRAW TAPE BAGS		SIZE
50-DTCL01	Draw Tape Bag, Clear	40 x 40 in / 102 x 102 cm
50-DTRD01	Draw Tape Bag, Red	40 x 40 in / 102 x 102 cm
50-DTBL01	Draw Tape Bag, Blue	40 x 40 in / 102 x 102 cm

LINERS/MISCELLANEOUS BAGS		
		SIZE
50-CTB01	Bag, Clear	14 x 18 in / 36 x 46 cm
50-CTB02	Bag, Clear	14 x 24 in / 36 cm x 61 cm
50-CTB03	Bag, Clear	16 x 9 x 30 in / 41 x 23 x 76 cm
50-CTB05	Bag, Clear	38 x 60 in / 97 x 152 cm
50-CTB06	Trash Liner Bag, Clear	33 x 38 in / 84 x 97 cm
50-CTB07	Bag, Red	21 x 19 x 46 in / 53 x 48 x 117 cm
50-CTB08	Bag, Clear	40 x 46 in / 102 x 117 cm
50-CTB09	Bag, Brown	38 x 38 in / 97 x 97 cm
50-CTB10	Bag, Green	38 x 38 in / 97 x 97 cm
50-CTB11	Bag, Brown	38 x 58 in / 97 x 147 cm

DRAWCORD BAGS		SIZE
50-DCCL01	Drawcord Bag, Clear	33 x 36 in / 84 x 91 cm
50-DCCL02	Drawcord Bag, Clear	38 x 38 in / 97 x 97 cm
50-DCCL03	Drawcord Bag, Clear	48 x 42 in / 122 x 107 cm
50-DCCL04	Drawcord Bag, Clear	40 x 40 in / 102 x 102 cm
50-DCCL05	Drawcord Bag, Clear	54 x 42 in / 137 x 107 cm
50-DCRD01	Drawcord Bag, Red	29 x 38 in / 74 x 97 cm
50-DCRD02	Drawcord Bag, Red	33 x 38 in / 84 x 97 cm
50-DCRD03	Drawcord Bag, Red	38 x 38 in / 97 x 97 cm
50-DCRD04	Drawcord Bag, Red	40 x 40 in / 102 x 102 cm
50-DCRD05	Drawcord Bag, Red	40 x 44 in / 102 x 112 cm
50-DCBL01	Drawcord Bag, Blue	29 x 38 in / 74 x 97 cm
50-DCBL02	Drawcord Bag, Blue	33 x 38 in / 84 x 97 cm
50-DCBL03	Drawcord Bag, Blue	38 x 38 in / 97 x 97 cm
50-DCBRN01	Drawcord Bag, Brown	38 x 38 in / 97 x 97 cm
50-DCBRN02	Drawcord Bag, Brown	54 x 42 in / 137 x 107 cm
50-DCGRN02	Drawcord Bag, Green	38 x 38 in / 97 x 97 cm
50-DCBYL01	Drawcord Bag, Yellow	29 x 38 in / 74 x 97 cm
50-DCBW01	Drawcord Bag, White	38 x 38 in / 97 x 97 cm

STANDARD CORE KITS

Kit Code: 49-RTK00

Absorbent Table Cover, Lift Sheet, Armboard Covers

Kit Code: 49-RTK01

Absorbent Table Cover, Lift Sheet

Kit Code: 49-RTK02

Absorbent Table Cover, Lift Sheet, Headrest Cover, Armboard Covers

Kit Code: 49-RTK04

Absorbent Table Cover, Lift Sheet, Headrest Cover, Armboard Covers, Drawcord Bag (Clear), Drawcord Bag (Red), Drawcord Bag (Blue), Kickbucket Bag (Red)