

PrimaPad™

TRENDELENBURG

Positioning System

With the PrimaPad system, patients can be safely placed in the Trendelenburg position during surgery. PrimaPad helps to limit patient movement during the procedure while also helping to reduce nerve damage, skin abrasions, and pressure ulcers.



Because PrimaPad is simple and fast to set up, it improves room turnover times and helps to reduce infections. It also eliminates the need for other position systems such as shoulder restraints and beanbags.

Different Procedures, Same System

PrimaPad is ideal for any surgery where the patient is placed in the Trendelenburg position, including:

- Robotic surgery
- Laparoscopy
- Gynecology
- Urology
- Colorectal procedures

Kit Components

- Memory foam pad
- Lift sheet
- Set of Body Straps

Head cradles and arm protectors are also available as kit components



Viscoelastic Foam Construction

Constructed with a special viscoelastic foam (sometimes called memory foam), PrimaPad is vital to many procedures:

- Limits Trendelenburg patient movement
- Conforms to the patient's anatomy
- Helps reduce pressure points
- Assists in maintaining normothermia through body heat retention
- Moisture wicking
- Hypoallergenic
- Not made with natural rubber latex

Features & Benefits

1 High Coefficient of Friction

Our special memory foam limits movement of the patient when placed in the Trendelenburg position.

1 High Density Viscoelastic Foam


The high density nature of the viscoelastic foam helps reduce tissue damage, skin abrasions, and pressure ulcers. The natural properties of the foam retain body heat, thus helping to maintain normothermia.

2 Included Body Straps

A set of hook-and-loop body straps is included that easily attach to the accessory rail, helping to further secure the patient, especially when tilted.

3 Strong Non-Woven Lift Sheet

The included non-woven lift sheet can be used for traditional patient transfer, then utilized to wrap and secure the patient arms.



Does not contain natural rubber latex, nano-silver additives or phthalates


Safe for patients and staff with possible allergies and prevention of leaching of potentially harmful compounds.

PrimaPad™ Trendelenburg Positioning Kits

CODE	DESCRIPTION	PACKAGING
50-582	PrimaPad™ Trendelenburg Positioning Kit, Universal, Non-Sterile	6 kits/cs
50-580	PrimaPad™ Trendelenburg Positioning Kit, Large, Non-Sterile	6 kits/cs
50-583	PrimaPad™ Trendelenburg Positioning Kit, X-Large, Non-Sterile	6 kits/cs
50-433	PrimaPad™ Trendelenburg Arm Protectors, Standard, Pair, Non-Sterile	4 kits/cs



PrimaPad™ Components A La Carte



Trendelenburg Positioning Pad

Packaging: 6 ea/cs

CODE	SIZE
50-602	20 x 35 x 1 in / 51 x 89 x 2.5 cm



Body Straps

Packaging: 60 pr/cs

CODE	SIZE
50-KB1023	4 x 61 in / 10 x 155 cm



Lift Sheet

Packaging: 50 ea/cs

CODE	SIZE
50-LFT008	12 x 80 in / 30 x 205 cm