

# Equipment Covers

## Endo Camera

| CODE                                | DESCRIPTION (WxL)  | CASE QTY | CASE DIMENSIONS        |
|-------------------------------------|--|----------|------------------------|
| <input type="radio"/> N 44-0596PR   | Camera Cover, Perforated Tip, w/ Plastic Ring Insert, 5x96" (13x244cm) | 100      | 19x15x12" (47x38x30cm) |
| <input type="radio"/> S 44-0596PR-S | Camera Cover, Perforated Tip, w/ Plastic Ring Insert, 5x96" (13x244cm) | 25       | 15x10x14" (37x25x36cm) |
| <input type="radio"/> N 44-0796PR   | Camera Cover, Perforated Tip, w/ Plastic Ring Insert, 7x96" (18x244cm) | 100      | 20x13x13" (52x32x32cm) |
| <input type="radio"/> S 44-0796PR-S | Camera Cover, Perforated Tip, w/ Plastic Ring Insert, 7x96" (18x244cm) | 25       | 19x15x12" (47x38x30cm) |
| <input type="radio"/> N 44-0996PR   | Camera Cover, Perforated Tip, w/ Plastic Ring Insert, 9x96" (23x244cm) | 100      | 25x13x17" (64x33x43cm) |
| <input type="radio"/> S 44-0996PR-S | Camera Cover, Perforated Tip, w/ Plastic Ring Insert, 9x96" (23x244cm) | 25       | 19x17x11" (48x43x27cm) |

### Materials

- Clear Polyethylene (PE) Film
- Plastic Ring Insert

### Product Attributes

- Adhesive securement strap
- Plastic insertion ring

### Sterilization

**Sterile:** Sterilized by EO (ethylene oxide). Product is sterile unless the package has been damaged or opened.  
**Non-Sterile:** If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packager. This includes further packaging and sterilization according to the validated processes of the re-packager. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

### Shelf Life

**Sterile:** 5 years  
**Non-Sterile:** Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

### Product Certifications/Standards


- Sterile: CE certified as Class I, according to MDD93/42/EEC Annex IX, Rule 1
- Non-Sterile: CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1

### Manufacturing Certifications/Registrations

- EN ISO 13485:2016
- Medical Device Single Audit Process Certified

### Storage

- General warehouse conditions



Alleset  
 4142 Industry Way  
 Flowery Branch, GA 30542,  
 USA  
[www.alleset.com](http://www.alleset.com)  
 SM-TDS-002 2023  
 Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.

GRI Medical & Electronic Technology Co., LTD  
 1805 Honggao Road,  
 314031 Jiaxing, Zhejiang  
 PEOPLE'S REPUBLIC OF CHINA  
[www.alleset.com](http://www.alleset.com)  
 US FDA Establishment Registration #3004911384

Emergo Europe B.V.  
 Westervoortsedijk 60,  
 6827 AT Arnhem,  
 The Netherlands  
[www.alleset.com](http://www.alleset.com)