

CORE Room Turnaround Kits®

Drawcord Bag

CODE	DESCRIPTION (LxW)	CASE QTY	CASE DIMENSIONS
50-DCCL01	Drawcord Bag, Clear, 33x36" (84x91cm)	100	16x12.5x10" (40x32x25cm)
50-DCCL02	Drawcord Bag, Clear, 38x38" (97x97cm)	100	16x12.5x10" (40x32x25cm)
50-DCCL03	Drawcord Bag, Clear, 48x42" (122x107cm)	100	22x12x12" (56x30.5x30.5cm)
50-DCCL04	Drawcord Bag, Clear, 40x40" (102x102cm)	100	16x12.5x10" (40x32x25cm)
50-DCCL05	Drawcord Bag, Clear, 54x42" (137x107cm)	100	16x15x5" (41x39x12cm)

Materials


- Clear Polyethylene (PE) Film
- Black Polypropylene (PP) Cord

Product Attributes

- Drawcord closure with cord lock

Sterilization

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packer. This includes further packaging and sterilization according to the validated processes of the re-packer. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.



Shelf Life

Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

Product Certifications/Standards

- Non-Sterile: CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1

Manufacturing Certifications/Registrations

- EN ISO 13485:2016
- Medical Device Single Audit Process Certified

Storage


- General warehouse conditions



Alleset
4142 Industry Way
Flowery Branch, GA 30542,
USA
www.alleset.com
SM-TDS-002 2022
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD
1805 Honggao Road,
314031 Jiaxing, Zhejiang
PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
US FDA Establishment Registration #3004911384



GRI-Alleset Limited B.V.
Kingsfordweg 151,
1043 GR Amsterdam,
The Netherlands
www.alleset.com

