








Allecare

Foam Wheelchair Cushions

CODE	DESCRIPTION (L x W x H)	CASE QTY	CASE DIMENSIONS
 FWC1616-2	Foam Wheelchair Cushion, Black 16x16x2" (41x41x5cm)	4	18x13x20.5" (46x33x52cm)
 FWC1816-2	Foam Wheelchair Cushion, Black 18x16x2" (46x41x5cm)	4	18x13x20.5" (46x33x52cm)
 FWC2016-2	Foam Wheelchair Cushion, Black 20x16x2" (51x41x5cm)	4	18x13x20.5" (46x33x52cm)
 FWC1616-3	Foam Wheelchair Cushion, Black 16x16x3" (41x41x8cm)	4	18x13x20.5" (46x33x52cm)
 FWC1816-3	Foam Wheelchair Cushion, Black 18x16x3" (46x41x8cm)	4	18x13x20.5" (46x33x52cm)
 FWC2016-3	Foam Wheelchair Cushion, Black 20x16x3" (51x41x8cm)	4	18x13x20.5" (46x33x52cm)
 FWC2218-3	Foam Wheelchair Cushion, Black 22x18x3" (56x46x8cm)	4	24x18x20" (61x46x51cm)

Materials

- Polyurethane (PU) Foam
- Polyester Cover
- Hook and Loop

Product Attributes

- High density foam
- Hook and loop strap
- Passes California Technical Bulletin 117-2013

Shelf Life


- 5 years

Manufacturing Certifications/Registrations

- EN ISO 13485:2016
- Medical Device Single Audit Process Certified

Storage

- General warehouse conditions



Alleset
4142 Industry Way
Flowery Branch, GA 30542,
USA

www.alleaset.com

SM-TDS-003 2023

Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



TENNESSEE FOAM, LLC
910 Callahan Dr. Ste 108,
Knoxville, TN 37912
USA

www.alleaset.com

US FDA Establishment Registration #3017728970

