

Protective Apparel

Extra Protective Lab Coat

CODE	DESCRIPTION	SIZE	CASE QTY	CASE DIMENSIONS
60-4031	Heavy Weight Tri-Layer Fabric, Knit Collar and Cuffs, White	M	30	16x10x12" (41x25x30cm)
60-4032	Heavy Weight Tri-Layer Fabric, Knit Collar and Cuffs, White	L	30	16x10x12" (41x25x30cm)
60-4033	Heavy Weight Tri-Layer Fabric, Knit Collar and Cuffs, White	XL	30	16x10x12" (41x25x30cm)
60-4034	Heavy Weight Tri-Layer Fabric, Knit Collar and Cuffs, White	2XL	30	16x10x12" (41x25x30cm)

Materials

- Multilayered SMS Material with Antistatic Properties and Alcohol Repellency

Product Attributes

- Heavy weight material
- Knit collar and cuffs
- Three pockets
- Hanging loop

Sterilization

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packer. This includes further packaging and sterilization according to the validated processes of the re-packer. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

Shelf Life

Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

Product Certifications/Standards


- CE certified as Category 1 PPE per Regulation (EU) 2016/42

Manufacturing Certifications/ Registrations

- EN ISO 13485:2016

Storage

- General warehouse conditions



Extra Protection Lab Coat with Knit Collar, White

Code	A: Length	B: Width
60-4031	40" (101.6cm)	24.5" (62cm)
60-4032	42" (106.7cm)	26.5" (67cm)
60-4033	44" (112cm)	28.5" (72cm)
60-4034	46" (116.8cm)	30.5" (77cm)
Tolerance	±1" (2.54cm)	±1" (2.54cm)



Alleset
4142 Industry Way
Flowery Branch, GA 30542,
USA
www.alleset.com
SM-TDS-002 2022
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.

GRI Medical & Electronic Technology Co., LTD
1805 Honggao Road,
314031 Jiaxing, Zhejiang
PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
US FDA Establishment Registration #3004911384

Emergo Europe B.V.
Westervoortsedijk 60,
6827 AT Arnhem,
The Netherlands
www.emergobyul.com

