

CORE Room Turnaround Kits®

Disposable Linens

CODE	DESCRIPTION (WxL)	CASE QTY	CASE DIMENSIONS
50-LFT001	Lift Sheet, 34x66" (86x168cm)	100	23x17x13" (58x43x34cm)
50-LFT002	Lift Sheet, 40x92" (102x234cm)	50	24x14x15" (62x35x38cm)
50-LFT003	Lift Sheet, 40x80" (102x203cm)	50	22x13x13" (55x33x32cm)
50-LFT004	Lift Sheet, 34x66" (86x168cm)	100	23x17x13" (58x43x34cm)
50-LFT005	Lift Sheet, 34x66" (86x168cm)	50	19x17x13" (48x43x34cm)
50-LFT006	Lift Sheet, 79x98" (200x250cm)	30	21x14x18" (54x35x46cm)

Materials

- Polyethylene Terephthalate (PET)

Product Attributes

- Non-woven PET spunlace construction

Sterilization

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packer. This includes further packaging and sterilization according to the validated processes of the re-packer. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

Shelf Life

Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

Product Certifications/Standards


CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1

Manufacturing Certifications/Registrations

- EN ISO 13485:2016
- Medical Device Single Audit Process Certified

Storage

- General warehouse conditions



Alleset
4142 Industry Way
Flowery Branch, GA 30542,
USA
www.alleset.com
SM-TDS-002 2022
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD
1805 Honggao Road,
314031 Jiaxing, Zhejiang
PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
US FDA Establishment Registration #3004911384



Alleset Limited B.V.
Kingsfordweg 151,
1043 GR Amsterdam,
The Netherlands
www.alleset.com

