

# Underpads

## Allesorb® STANDARD

CODE	DESCRIPTION	CASE QTY	CASE DIMENSIONS
99-UPST2436	Underpad Allesorb® STANDARD, 24x36" (61x91cm)	10/Pack, 6/CS	19x13x11" (48x33x28cm)
99-UPST3136	Underpad Allesorb® STANDARD, 31X36" (79x91cm)	10/Pack, 4/CS	16x14x11" (40x35x29cm)

### Materials

- Breathable Polyethylene (PE) Film and Hydrophilic Spunbond Polypropylene (PP)
- Super Absorbent Polymer Core

### Product Attributes

- Unique fabric construction wicks fluids away quickly
- Middle layer traps fluids
- Bottom layer allows air permeability for heat transfer while maintaining a moisture barrier
- Easy to use

### Sterilization

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packager. This includes further packaging and sterilization according to the validated processes of the re-packager. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

### Shelf Life

Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

### Product Certifications/Standards


- CE certified as class I, according to MDR 2017/745, Annex VIII, Rule 1

### Manufacturing Certifications/Registrations

- EN ISO 13485:2016

### Storage

- General warehouse conditions



Alleset  
4142 Industry Way  
Flowery Branch, GA 30542,  
USA  
www.alleset.com  
SM-TDS-002 2022  
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD  
1805 Honggao Road,  
314031 Jiaxing, Zhejiang  
PEOPLE'S REPUBLIC OF CHINA  
www.alleset.com  
US FDA Establishment Registration #3004911384



Emergo Europe B.V.  
Prinsessegracht 20  
2514 AP The Hague  
The Netherlands  
www.emergobyul.com

