

Surgical Accessories

Suction Tubing

CODE	DESCRIPTION	CASE QTY	CASE DIMENSIONS
70-0506EC-S	Suction Tubing w/ EU Connector, 3/16inx6ft (5mmx1.8m)	25	16x11x11" (41x29x28cm)
70-0510EC-S	Suction Tubing w/ EU Connector, 3/16inx10ft (5mmx3.0m)	25	16x11x11" (41x29x28cm)
70-0512EC-S	Suction Tubing w/ EU Connector, 3/16inx12ft (5mmx3.7m)	25	16x11x11" (41x29x28cm)
70-0520EC-S	Suction Tubing w/ EU Connector, 3/16inx20ft (5mmx6.1m)	25	18x13x11" (46x34x28cm)
70-0606EC-S	Suction Tubing w/ EU Connector, 1/4inx6ft (6mmx1.8m)	25	16x11x11" (41x29x28cm)
70-0610EC-S	Suction Tubing w/ EU Connector, 1/4inx10ft (6mmx3.0m)	25	15x14x14" (39x36x36cm)
70-0612EC-S	Suction Tubing w/ EU Connector, 1/4inx12ft (6mmx3.7m)	25	15x14x14" (39x36x36cm)
70-0620EC-S	Suction Tubing w/ EU Connector, 1/4inx20ft (6mmx6.1m)	25	16x15x18" (41x38x45cm)
70-0706EC-S	Suction Tubing w/ EU Connector, 9/32inx6ft (7mmx1.8m)	25	16x11x11" (41x29x28cm)
70-0710EC-S	Suction Tubing w/ EU Connector, 9/32inx10ft (7mmx3.0m)	25	15x14x14" (39x36x36cm)
70-0712EC-S	Suction Tubing w/ EU connector, 9/32inx12ft (7mmx3.7m)	25	15x14x14" (39x36x36cm)
70-0720EC-S	Suction Tubing w/ EU connector, 9/32inx20ft (7mmx6.1m)	25	18x16x18" (45x41x45cm)

Materials

- Polyvinyl Chloride (PVC)

Product Attributes

- Low coil memory and greater flexibility
- Increased wall thickness helps prevent tubing collapse

Sterilization

Sterile: Sterilized by EO (ethylene oxide). Product is sterile unless the package has been damaged or opened.

Shelf Life

Sterile: 5 years

Product Certifications/Standards



- Sterile: CE certified as Class IIa, according to MDD93/42/EEC Annex IX, Rule 2

Manufacturing Certifications/Registrations


- EN ISO 13485:2016
- Medical Device Single Audit Process Certified

Storage


- General warehouse conditions



3/16" (5mm)



1/4" (6mm)



9/32" (7mm)



Alleset
4142 Industry Way
Flowery Branch, GA 30542,
USA
www.alleset.com
SM-TDS-003 2023
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD
1805 Honggao Road,
314031 Jiaxing, Zhejiang
PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
US FDA Establishment Registration #3004911384



Emergo Europe B.V.
Westervoortsedijk 60,
6827 AT Arnhem,
The Netherlands.
www.emergobyul.com

