

# Allecare

## Back Support Cushions

CODE	DESCRIPTION (L x H)	CASE QTY	CASE DIMENSIONS
WBS1617	Contoured Wheelchair Back Support, 16x17" (41x43cm)	4	18x13x20.5" (46x33x52cm)
WBS1817	Contoured Wheelchair Back Support, 18x17" (46x43cm)	4	18x13x20.5" (46x33x52cm)
WBS2017	Contoured Wheelchair Back Support, 20x17" (51x43cm)	4	18x13x20.5" (46x33x52cm)
WBS2217	Contoured Wheelchair Back Support, 22x17" (56x43cm)	4	24x18x20" (61x46x51cm)

### Materials

- Polyurethane (PU) Foam
- Polyester Cover
- Hook and Loop

### Product Attributes

- High density foam
- Contoured cut for back support
- Hook and loop securement strap
- Passes California Technical Bulletin 117-2013

### Shelf Life


- 5 years

### Manufacturing Certifications/Registrations

- EN ISO 13485:2016
- Medical Device Single Audit Process Certified

### Storage

- General warehouse conditions



Alleset  
4142 Industry Way  
Flowery Branch, GA 30542,  
USA

[www.alleaset.com](http://www.alleaset.com)

SM-TDS-003 2023

Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.

THINK.MAKE.PROTECT.



TENNESSEE FOAM, LLC  
910 Callahan Dr. Ste 108,  
Knoxville, TN 37912  
USA

[www.alleaset.com](http://www.alleaset.com)

US FDA Establishment Registration #3017728970

A GRI COMPANY

