

NEEDLE COUNTERS


Our needle counter product portfolio offers a wide variety of configurations and colors. The line is composed of our SECURE and SELECT product lines which feature many design options to facilitate sharps counting and proper disposal of surgical sharps during the surgical procedure.



Available in yellow or red, configurations include: slide locks for additional safety, counters with foam strips, foam blocks, magnets, with or without blade removers, adhesive securement strips. Counts ranging from 10 to 100.

SELECT Needle Counters

- Specially selected foam does not shrink during sterilization
- Adhesive tabs secure the counter to desired workspace
- Easy-to-read numbers allow easy count and placement of blades and needles
- Clear lid available, providing visualization when closed
- Does not contain natural rubber latex, nano-silver additives or phthalates



SECURE Needle Counters

Our SECURE Needle Counters maintain all of the features of our SELECT Needle Counters plus the following:

- Built-in blade remover ensures safety
- Three (3) hinge design prevents inadvertent opening and separation
- Passive locking latch engages automatically and securely
- Interlocking shells prevent sharps from escaping when closed