

Allecare

Gel Wheelchair Cushions

CODE	DESCRIPTION (L x W x H)	CASE QTY	CASE DIMENSIONS
GWC1616-2	Gel Wheelchair Cushion, 16x16x2" (41x41x5cm)	4	18x13x20.5" (46x33x52cm)
GWC1816-2	Gel Wheelchair Cushion, 18x16x2" (46x41x5cm)	4	18x13x20.5" (46x33x52cm)
GWC2016-2	Gel Wheelchair Cushion, 20x16x2" (51x41x5cm)	4	18x13x20.5" (46x33x52cm)
GWC2218-3	Gel Wheelchair Cushion, 22x18x3" (41x41x8cm)	4	24x18x20" (61x46x51cm)
GWC2418-3	Gel Wheelchair Cushion, 24x18x3" (46x41x8cm)	4	24x18x20" (61x46x51cm)
GWC2618-3	Gel Wheelchair Cushion, 26x18x3" (51x41x8cm)	4	24x20x36" (61x51x91cm)
GWC2620-3	Gel Wheelchair Cushion, 26x20x3" (56x46x8cm)	4	24x20x36" (61x51x91cm)
GWC2820-3	Gel Wheelchair Cushion, 28x20x3" (56x46x8cm)	4	24x20x36" (61x51x91cm)
GWC3020-3	Gel Wheelchair Cushion, 30x20x3" (56x46x8cm)	4	24x20x36" (61x51x91cm)
GWC1416-3	Gel Wheelchair Cushion, 14x16x3" (56x46x8cm)	4	18x13x20.5" (46x33x52cm)
GWC1616-3	Gel Wheelchair Cushion, 16x16x3" (56x46x8cm)	4	18x13x20.5" (46x33x52cm)
GWC1816-3	Gel Wheelchair Cushion, 18x16x3" (56x46x8cm)	4	18x13x20.5" (46x33x52cm)
GWC2016-3	Gel Wheelchair Cushion, 20x16x3" (56x46x8cm)	4	18x13x20.5" (46x33x52cm)
GWC2018-3	Gel Wheelchair Cushion, 20x18x3" (56x46x8cm)	4	18x13x20.5" (46x33x52cm)

Materials

- Polyurethane (PU) Foam
- Water Based Gel
- Polyethylene (PE) Bladder Bag
- Polyester Cover

Product Attributes

- High density foam
- Hook and loop Strap
- Multi-channel gel bladder
- Passes California Technical Bulletin 117-2013

Shelf Life


- 5 years

Manufacturing Certifications/Registrations

- EN ISO 13485:2016
- Medical Device Single Audit Process Certified

Storage

- General warehouse conditions



Alleset
4142 Industry Way
Flowery Branch, GA 30542,
USA

www.alleaset.com

SM-TDS-003 2023

Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



TENNESSEE FOAM, LLC
910 Callahan Dr. Ste 108,
Knoxville, TN 37912
USA

www.alleaset.com

US FDA Establishment Registration #3017728970

