

Protective Apparel

ComfortGuard® Tri-Layer SMS Isolation Gown

CODE	DESCRIPTION	SIZE	CASE QTY	CASE DIMENSIONS
60-3070	Isolation Gown w/ Elastic Wrist, Tab Neck Closure, Closed Back, Yellow	L	100	21x14x22" (53x36x56cm)
60-3071	Isolation Gown w/ Elastic Wrist, Tab Neck Closure, Closed Back, Yellow	XL	100	21x14x22" (53x36x56cm)
60-3072	Isolation Gown w/ Elastic Wrist, Tab Neck Closure, Closed Back, Blue	L	100	21x14x22" (53x36x56cm)
60-3073	Isolation Gown w/ Elastic Wrist, Tab Neck Closure, Closed Back, Blue	XL	100	21x14x22" (53x36x56cm)
60-3074	Isolation Gown w/ Knit Cuff, Tab Neck Closure, Closed Back, Yellow	L	100	21x14x28" (53x35x70cm)
60-3075	Isolation Gown w/ Knit Cuff, Tab Neck Closure, Closed Back, Yellow	XL	100	21x14x28" (53x35x70cm)
60-3076	Isolation Gown, Thumb Loop w/ Elastic, Over the Head Style, Closed Back, Yellow	L	100	21x14x22" (53x36x56cm)
60-3077	Isolation Gown, Thumb Loop w/ Elastic, Over the Head Style, Closed Back, Yellow	XL	100	21x14x22" (53x36x56cm)
60-3078	Isolation Gown w/ Elastic Wrist, Hook and Loop Closure, Closed Back, Yellow	L	100	21x14x22" (53x36x56cm)
60-3079	Isolation Gown w/ Elastic Wrist, Hook and Loop Closure, Closed Back, Yellow	XL	100	21x14x22" (53x36x56cm)

Materials

- Multilayered SMS Material

Product Attributes

- Full coverage
- Regular and extra-long lengths
- Tab neck or over the head design
- Tie waist
- Elastic cuff, knit cuff or thumb hook options
- Sonic bonded seam construction

Sterilization

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packager. This includes further packaging and sterilization according to the validated processes of the re-packager. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

Shelf Life

Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

Product Certifications / Standards


- EN13795: 2011 Standard Performance
- ANSI/AAMI PB70:2022 Level 3 Liquid Barrier Protection
- CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1

Manufacturing Certifications / Registrations


- EN ISO 13485:2016
- Medical Device Single Audit Process Certified

Storage

- General warehouse conditions



60-3070, 60-3071




60-3072, 60-3073




Protective Apparel


ComfortGuard® Tri-Layer SMS Isolation Gown



60-3074, 60-3075



60-3076, 60-3077



60-3078, 60-3079

Code	A: Chest Width	B: Total Length	B: Sleeve Down Length
60-3070	28.5" (72cm)	46" (117cm)	22" (57cm)
60-3071	30.5" (77cm)	48" (122cm)	23" (60cm)
60-3072	28.5" (72cm)	46" (117cm)	22" (57cm)
60-3073	30.5" (77cm)	48" (122cm)	22" (57cm)
60-3074	28.5" (72cm)	46" (117cm)	22" (57cm)
60-3075	30.5" (77cm)	48" (122cm)	23" (60cm)
60-3078	28.5" (72cm)	46" (117cm)	22" (57cm)
60-3079	30.5" (77cm)	48" (122cm)	22" (57cm)
Tolerance	±1" (2.54cm)	±1" (2.54cm)	±1" (2.54cm)

Code	A: Chest Width	B: Total Length	B: Sleeve Up Length
60-3076	28.5" (72cm)	46" (117cm)	27.5" (70cm)
60-3077	30.5" (77cm)	48" (122cm)	29.5" (75cm)
Tolerance	±0.5" (1.27cm)	±0.5" (1.27cm)	±0.5" (1.27cm)



Alleset
4142 Industry Way
Flowery Branch, GA 30542,
USA

www.alleset.com

SM-TDS-003 2023

Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.

THINK.MAKE.PROTECT.



GRI Medical & Electronic Technology Co., LTD
1805 Honggao Road,
314031 Jiaxing, Zhejiang
PEOPLE'S REPUBLIC OF CHINA
www.alleset.com

US FDA Establishment Registration #3004911384



Emergo Europe B.V.
Westervoortsedijk 60,
6827 AT Arnhem,
The Netherlands
www.emergobyul.com

