

# Protective Apparel

## Isolation Gown

CODE	DESCRIPTION	SIZE	CASE QTY	CASE DIMENSIONS
60-3030	Single Layer Spunbond Isolation Gown, Non-Reinforced, Yellow	L	100	21x14x14" (53x35x36cm)
60-3032	Single Layer Spunbond Isolation Gown, Non-Reinforced, Yellow	XL	100	21x14x14" (53x35x36cm)
60-3021	Single Layer Spunbond Isolation Gown, Non-Reinforced, Yellow	XL, X-Long	100	21x14x14" (53x35x36cm)
60-3033	Single Layer Spunbond Isolation Gown, Non-Reinforced, Blue	L	100	21x14x14" (53x35x36cm)
60-3034	Single Layer Spunbond Isolation Gown, Non-Reinforced, Blue	XL	100	21x14x14" (53x35x36cm)
60-3035	Single Layer Spunbond Isolation Gown, Non-Reinforced, Blue	XL, X-Long	100	21x14x14" (53x35x36cm)

### Materials

Hydrophilic Spunbond Polypropylene with Antistatic Properties

### Product Attributes

- Full coverage
- Tie neck closure
- Tie waist
- Elastic cuffs

### Sterilization

**Non-Sterile:** If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packer. This includes further packaging and sterilization according to the validated processes of the re-packer. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

### Shelf Life

**Non-Sterile:** Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

### Product Certifications/Standards



- CE certified as class I, according to MDR 2017/745, Annex VIII, Rule 1

### Manufacturing Certifications/ Registrations

- EN ISO 13485:2016
- Medical Device Single Audit Process Certified

### Storage

- General warehouse conditions



Code	A: Total Length	B: Chest Length	C: Sleeve Length
60-3030	46" (117cm)	28" (72cm)	22" (57cm)
60-3032	48" (122cm)	30" (77cm)	24" (60cm)
60-3021	56" (141cm)	30" (77cm)	24" (60cm)
60-3033	46" (117cm)	28" (72cm)	22" (57cm)
60-3034	48" (122cm)	30" (77cm)	24" (60cm)
60-3035	50" (127cm)	28" (72cm)	22" (57cm)
Tolerance	±1" (2.54cm)	±1" (2.54cm)	±1" (2.54cm)



Alleset  
4142 Industry Way  
Flowery Branch, GA 30542,  
USA  
www.alleset.com  
SM-TDS-003 2023  
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.

GRI Medical & Electronic Technology Co., LTD  
1805 Honggao Road,  
314031 Jiaxing, Zhejiang  
PEOPLE'S REPUBLIC OF CHINA  
www.alleset.com  
US FDA Establishment Registration #3004911384

Emergo Europe B.V.  
Westervoortsewijk 60,  
6827 AT Arnhem,  
The Netherlands  
www.emergobyul.com

