

Protective Apparel

Procedure Gown

CODE	DESCRIPTION	SIZE	CASE QTY	CASE DIMENSIONS
60-3212	Tri-Layer SMS Procedure Gown	L	60	21x14x18" (54x35x46cm)
60-3213	Tri-Layer SMS Procedure Gown	XL	60	21x14x18" (54x35x46cm)

Materials

- Multilayered SMS Material with Antistatic Properties and Alcohol Repellency

Product Attributes

- Full-coverage
- Alcohol-repellent
- Tie waist
- Hook and loop neck closure
- Raglan sleeve
- Knit cuff

Sterilization

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packer. This includes further packaging and sterilization according to the validated processes of the re-packer. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

Shelf Life

Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

Product Certifications/Standards


- ANSI/AAMI PB70:2012 Level 2 Liquid Barrier Protection
- CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1

Manufacturing Certifications/ Registrations

- EN ISO 13485:2016
- Medical Device Single Audit Process Certified

Storage

- General warehouse conditions



Code	A: Total Length	B: Chest Length	C: Sleeve Length
60-3212	46" (117cm)	28.5" (72cm)	22.5" (57cm)
60-3213	48" (122cm)	30.5" (77cm)	23.5" (60cm)
Tolerance	±1" (2.54cm)	±1" (2.54cm)	±1" (2.54cm)



Alleset
4142 Industry Way
Flowery Branch, GA 30542,
USA
www.alleset.com
SM-TDS-002 2022
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.

THINK.MAKE.PROTECT.



GRI Medical & Electronic Technology Co., LTD
1805 Honggao Road,
314031 Jiaxing, Zhejiang
PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
US FDA Establishment Registration #3004911384



Emergo Europe B.V.
Westervoortsedijk 60,
6827 AT Arnhem,
The Netherlands
www.emergobyul.com

A GRI COMPANY

