

Protective Apparel

Coveralls

CODE	DESCRIPTION	SIZE	CASE QTY	CASE DIMENSIONS
60-6001	Regular Protection Coverall w/ Open Cuff and Ankles, Blue	L	30	16x10x12" (41x25x30cm)
60-6002	Regular Protection Coverall w/ Open Cuff and Ankles, Blue	XL	30	16x10x12" (41x25x30cm)

Materials

- Multilayered SMS Material with Antistatic Properties

Product Attributes

- Open sleeve and cuff
- Full zipper
- Traditional collar

Sterilization

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packer. This includes further packaging and sterilization according to the validated processes of the re-packer. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

Shelf Life

Non-Sterile: Establishes at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

Product Certifications/Standards


- CE certified as Category 1 PPE per Regulation (EU) 2016/425

Manufacturing Certifications/Registrations

- EN ISO 13485:2016

Storage

- General warehouse conditions



Item	Description	L	XL	Tolerance
A	Total Length	65.5" (166cm)	68.5" (174cm)	±1"
B	Chest	25.5" (65cm)	27.5" (70cm)	±1"
C	Sleeve	22.5" (57cm)	23" (58cm)	±1"
D	Inseam Length	29.5" (75cm)	30.5" (77cm)	±1"



Alleset
 4142 Industry Way
 Flowery Branch, GA 30542,
 USA
 www.alleset.com
 SM-TDS-002 2022
 Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.

GRI Medical & Electronic Technology Co., LTD
 1805 Honggao Road,
 314031 Jiaxing, Zhejiang
 PEOPLE'S REPUBLIC OF CHINA
 www.alleset.com
 US FDA Establishment Registration #3004911384

Emergo Europe B.V.
 Westervoortsedijk 60,
 6827 AT Arnhem,
 The Netherlands
 www.emergobyul.com

