

Protective Apparel

Isolation Gown

CODE	DESCRIPTION	SIZE	CASE QTY	CASE DIMENSIONS
60-3036	Tri-Layer SMS Isolation Gown, Non-Reinforced	L	100	21x14x14" (53x35x36cm)
60-3037	Tri-Layer SMS Isolation Gown, Non-Reinforced	XL	100	21x14x14" (53x35x36cm)
60-3025	Tri-Layer SMS Isolation Gown, Non-Reinforced	XL, X-Long	100	21x14x14" (53x35x36cm)

Materials

- Multilayered SMS Material with Antistatic Properties

Product Attributes

- Full coverage
- Regular and extra-long lengths
- Tie waist
- Tie neck closure
- Elastic cuffs

Sterilization

Non-Sterile: If a product is delivered non-sterile and is ultimately for use in sterile environments, it is intended to be sterilized by the re-packer. This includes further packaging and sterilization according to the validated processes of the re-packer. These products are to be sterilized by ethylene oxide (EO). Other methods of sterilization have not been validated by Alleset.

Shelf Life

Non-Sterile: Established at 5 years based upon sterile packaged product aging studies and normal warehouse storage conditions. After 5 years, the product should be performance tested to ensure compliance to specifications.

Product Certifications/Standards


- CE certified as class I, according to MDR 2017/745, Annex VIII, Rule 1

Manufacturing Certifications/ Registrations

- EN ISO 13485:2016
- Medical Device Single Audit Process Certified

Storage

- General warehouse conditions



Code	A: Total Length	B: Chest Length	C: Sleeve Length
60-3036	46" (117cm)	28" (72cm)	22" (57cm)
60-3037	48" (122cm)	30" (77cm)	24" (60cm)
60-3025	56" (141cm)	30" (77cm)	24" (60cm)
Tolerance	±1" (2.54cm)	±1" (2.54cm)	±1" (2.54cm)



Alleset
4142 Industry Way
Flowery Branch, GA 30542,
USA
www.alleset.com
SM-TDS-003 2023
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.



GRI Medical & Electronic Technology Co., LTD
1805 Honggao Road,
314031 Jiaxing, Zhejiang
PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
US FDA Establishment Registration #3004911384



Emergo Europe B.V.
Westervoortsedijk 60,
6827 AT Arnhem,
The Netherlands
www.emergobyul.com

