



Patient Positioning

PrimaPad - Trendelenburg Position System

CODE	DESCRIPTION	CASE QTY	CASE DIMENSIONS
 50-582	PrimaPad - Trendelenburg Positioning System, 35x20" (89x51cm)	6 kits/case	37x22x9" (95x55x24cm)

Materials

PrimaPad

-  Polyurethane (PU) Memory Foam, Green

Lift Sheet

-  Spunlace

Body Straps

-  Composite Foam

Product Attributes

- 1" (2.5 cm) thick viscoelastic "memory" foam pad
- Integrated pad securement straps
- Spunlace lift sheet for patient transfer
- Foam patient positioning straps

Sterilization

Non-Sterile: This product is intended to be used non-sterile. If the repackager intends to sterilize this product, a validation must be conducted to demonstrate their sterilization process will meet all the ISO 11135 and EN 10993 requirements to demonstrate sterility and EO/ECH residual compliance. It must be done according to their validated sterilization process.

Shelf Life

Not Available

Product Certifications/Standards

- CE certified as Class I, according to MDR 2017/745 Annex VIII, Rule 1

Manufacturing Certifications/Registrations


- EN ISO 13485:2016

Storage


- General warehouse conditions

Components per Kit


Description	Quantity
PrimaPad	1 each
Lift Sheet	1 each
Body Straps	1 pair



PrimaPad



Lift Sheet



Body Straps



Alleset
4142 Industry Way
Flowery Branch, GA 30542,
USA
www.alleset.com
SM-TDS-002 2022
Country of Origin Abbreviations per ISO 3166-1, Alpha-2 Codes.

THINK.MAKE.PROTECT.



GRI Medical & Electronic Technology Co., LTD
1805 Honggao Road,
314031 Jiaxing, Zhejiang
PEOPLE'S REPUBLIC OF CHINA
www.alleset.com
US FDA Establishment Registration #3004911384



Emergo Europe B.V.
Westervoortsedijk 60,
6827 AT Arnhem,
The Netherlands
www.emergobyul.com

A GRI COMPANY

