

# SURGICAL GOWNS & DRAPES

 COMFORT GUARD<sup>®</sup>

 STERIGUARD<sup>®</sup>



# ABOUT GRI-ALLESET

## DELIVERING PRODUCT SOLUTIONS™

GRI-Alleset is your “one-stop shop” for the widest variety of self-manufactured surgical disposable products. At the same time other companies are rationalizing their product offering, GRI-Alleset approaches each market and its customers with one question in mind ... “What is essential to the healthcare provider?”

As a wholly-owned subsidiary of Global Resources International, Inc. (GRI), a market leading medical device manufacturer, the mission of GRI-Alleset is to develop and deliver high-quality, cost-effective product solutions to healthcare providers and supply chain partners as a means to improving the value proposition across the global healthcare marketplace. Initiating at product ideation and development, through a wide variety of manufacturing capabilities, and offering several cost-effective distribution options, we maintain control of cost and quality throughout the entire product lifecycle.

As an organization we pride ourselves on delivering “true value” - The perfect combination of quality products, competitive pricing, and unparalleled support services.

### Reducing Risks

- ▶ Head-to-toe protection for clinicians and patients
- ▶ Barrier fabrics that leverage the latest, most advanced protective technologies
- ▶ Innovative design features that maximize risk prevention

### Striving For Better

- ▶ Innovative designs that help provide superior protection and performance while controlling cost
- ▶ Solutions that match the needs of the latest surgical methods, technologies and equipment
- ▶ Flexibility to customize solutions based on the evolving needs of surgeons and clinicians worldwide

### Affording The Best

- ▶ Self-manufactured solutions versus leveraging middle-men or “brokers”
- ▶ Vertically integrated solutions – from design to manufacturing to sterilization
- ▶ Optimally designed manufacturing and supply solutions minimize costs within the production process



DELIVERING PRODUCT SOLUTIONS™

# TABLE OF CONTENTS

|                                                                                     |                               |    |
|-------------------------------------------------------------------------------------|-------------------------------|----|
|    | Surgical Gowns . . . . .      | 4  |
|    | Angiography . . . . .         | 16 |
|    | Special Procedures . . . . .  | 24 |
|    | Basic Drapes . . . . .        | 28 |
|    | Cardiovascular . . . . .      | 32 |
|   | EENT / Neurology . . . . .    | 36 |
|  | General Drapes . . . . .      | 40 |
|  | Laparotomy . . . . .          | 46 |
|  | Lithotomy / Urology . . . . . | 52 |
|  | OB / GYN . . . . .            | 56 |
|  | Ophthalmic . . . . .          | 62 |
|  | Orthopedic . . . . .          | 68 |
|  | Equipment Covers . . . . .    | 78 |



# SURGICAL GOWNS

GRI-Alleset's ComfortGuard® Surgical Gowns have considered every element of the surgical gown as an opportunity to improve performance. From the unique formation of our fibers to the innovative design of our sleeves, we have built our gowns from the ground up with the challenges of different surgical procedures in mind.



## PROTECTION

Excellent barrier and strength properties help protect clinicians from bodily fluids. Low flammability and low lint help eliminate one source of post-operative patient complications.



## COMFORT

Excellent breathability created by advanced fabric technologies and innovative design, utilizing a lightweight, highly breathable cool-back panel to maximize comfort. Soft fabrics feel comfortable and stay quiet.



## FIT

Generous cuts with raglan and set-in sleeves and pleats on the back provide increased range of motion. Added gown length improves lower-body protection. Extra long hook-and-loop creates a more customized fit around the neckline. Extra soft, extended cuffs enhance comfort.



## SUSTAINABILITY


Innovative cool-back feature provides the comfort you demand, while helping reduce fossil fuels and carbon footprint by up to 20%. While your doctors take care of their patients, we are doing our part to ensure efficient and responsible manufacturing processes.



DELIVERING PRODUCT SOLUTIONS™

# GRI-Alleaset understands that every cut and stitch in surgery is critical. This is why we apply the same thinking to the design and construction of our surgical gowns.

GRI-Alleaset's ComfortGuard® Surgical Gowns provide all the protection needed, paired with comfort and fit that hours of surgery demands. For every clinician and for each procedure performed. ComfortGuard® Surgical Gowns deliver superior performance at maximum value throughout our product line.



## Breathable Protection

Triple-density materials transfer heat out of the garment while keeping fluids from getting in.



## Fiber-Perfect Softness

Superior flex, strength, and fine fibers create a unique soft construction for unrivaled comfort.



## Protective Length

Full coverage with low movement drag and maximum comfort.



## Cool Back Panel

Unique design combination of protection and ventilation.



## Generous Hook-and-Loop

Secure hook-and-loop aids donning, enhances fit, and maximizes comfort.



## Surgical-Fit Cuffs


Extended length and unique construction ensure a better fit and softer feel.



## Full-Range Sleeves

Raglan or Set-In design gives full range of motion in any position.





## 8000 SERIES

### Surgical Gown, Non-Reinforced, Set-In Sleeves

- ▶ Premium SMS Barrier Fabric
- ▶ Set-In Sleeves
- ▶ Adjustable Comfort-Soft Hook-and-Loop Closures
- ▶ Non-Reinforced

#### BASIC


| Size | <b>S</b> Sterile | Case Qty | <b>N</b> Non-Sterile | Case Qty | <b>N</b> Non-Sterile, Handi-Bin | Case Qty |
|------|------------------|----------|----------------------|----------|---------------------------------|----------|
| S    | 90-8000-S1       | 22       | 90-8000              | 32       |                                 |          |
| M    | 90-8010-S1       | 32       | 90-8010              | 40       |                                 |          |
| L    | 90-8020-S1       | 30       | 90-8020              | 40       | 90-8020-HB                      | 432      |
| XL   | 90-8030-S1       | 28       | 90-8030              | 36       | 90-8030-HB                      | 378      |
| XXL  | 90-8040-S1       | 26       | 90-8040              | 32       | 90-8040-HB                      | 324      |

#### SPECIALTY

| Size       | <b>S</b> Sterile | Case Qty | <b>N</b> Non-Sterile | Case Qty | <b>N</b> Non-Sterile, Handi-Bin | Case Qty |
|------------|------------------|----------|----------------------|----------|---------------------------------|----------|
| L, XLong   | 90-8022-S1       | 28       | 90-8022              | 36       |                                 |          |
| XL, XLong  | 90-8032-S1       | 26       | 90-8032              | 32       |                                 |          |
| XXL, XLong | 90-4042-S1       | 24       | 90-8042              | 32       |                                 |          |

#### SCRUB NURSE

| Size | <b>N</b> Non-Sterile, Handi-Bin       | Case Qty |
|------|---------------------------------------|----------|
| L    | 90-8020-W1-HB (Hand Towel & Overwrap) | 432      |
| XL   | 90-8030-W1-HB (Hand Towel & Overwrap) | 378      |

**1200 SERIES****Surgical Gown, Film-Reinforced, Set-In Sleeves**

- ▶ Premium SMS Barrier Fabric
- ▶ Set-In Sleeves
- ▶ Adjustable Comfort-Soft Hook-and-Loop Closures
- ▶ Film Reinforcement
- ▶ Sleeve-N-Sleeve Construction
- ▶ Extra Long, Extra Wide Film Reinforcement

## BASIC


| Size | <b>S</b> Sterile | Case Qty | <b>N</b> Non-Sterile | Case Qty | <b>N</b> Non-Sterile, Handi-Bin | Case Qty |
|------|------------------|----------|----------------------|----------|---------------------------------|----------|
| S    | 90-1200-S1       | 20       | 90-1200              | 32       |                                 |          |
| M    | 90-1210-S1       | 26       | 90-1210              | 28       |                                 |          |
| L    | 90-1220-S1       | 24       | 90-1220              | 36       |                                 |          |
| XL   | 90-1230-S1       | 24       | 90-1230              | 32       |                                 |          |
| XXL  | 90-1240-S1       | 22       | 90-1240              | 32       |                                 |          |

## SPECIALTY

| Size       | <b>S</b> Sterile | Case Qty | <b>N</b> Non-Sterile | Case Qty | <b>N</b> Non-Sterile, Handi-Bin | Case Qty |
|------------|------------------|----------|----------------------|----------|---------------------------------|----------|
| L, XLong   | 90-1222-S1       | 22       | 90-1222              | 32       |                                 |          |
| XL, XLong  | 90-1232-S1       | 20       | 90-1232              | 28       |                                 |          |
| XXL, XLong | 90-1242-S1       | 18       | 90-1242              | 24       |                                 |          |

## SCRUB NURSE

| Size | <b>N</b> Non-Sterile, Handi-Bin       | Case Qty |
|------|---------------------------------------|----------|
| L    | 90-1220-W1-HB (Hand Towel & Overwrap) | 378      |
| XL   | 90-1230-W1-HB (Hand Towel & Overwrap) | 324      |
| XXL  | 90-1240-W1-HB (Hand Towel & Overwrap) | 324      |



## 4000 SERIES

### Surgical Gown, Non-Reinforced, Set-In Sleeves

- ▶ SMS Barrier Fabric
- ▶ Set-In Sleeves
- ▶ Adjustable Comfort-Soft Hook-and-Loop Closures
- ▶ Non-Reinforced

#### BASIC


| Size | <b>S</b> Sterile | Case Qty | <b>N</b> Non-Sterile | Case Qty | <b>N</b> Non-Sterile, Handi-Bin | Case Qty |
|------|------------------|----------|----------------------|----------|---------------------------------|----------|
| S    | 90-4000-S1       | 22       | 90-4000              | 40       |                                 |          |
| M    | 90-4010-S1       | 30       | 90-4010              | 40       |                                 |          |
| L    | 90-4020-S1       | 28       | 90-4020              | 40       | 90-4020-HB                      | 432      |
| XL   | 90-4030-S1       | 26       | 90-4030              | 40       | 90-4030-HB                      | 378      |
| XXL  | 90-4040-S1       | 24       | 90-4040              | 36       | 90-4040-HB                      | 378      |

#### SPECIALTY

| Size       | <b>S</b> Sterile | Case Qty | <b>N</b> Non-Sterile | Case Qty | <b>N</b> Non-Sterile, Handi-Bin | Case Qty |
|------------|------------------|----------|----------------------|----------|---------------------------------|----------|
| L, XLong   | 90-4022-S1       | 24       | 90-4022              | 40       |                                 |          |
| XL, XLong  | 90-4032-S1       | 24       | 90-4032              | 36       |                                 |          |
| XXL, XLong | 90-4042-S1       | 22       | 90-4042              | 32       |                                 |          |

#### SCRUB NURSE

| Size | <b>N</b> Non-Sterile, Handi-Bin       | Case Qty |
|------|---------------------------------------|----------|
| L    | 90-4020-W1-HB (Hand Towel & Overwrap) | 432      |
| XL   | 90-4030-W1-HB (Hand Towel & Overwrap) | 378      |
| XXL  | 90-4040-W1-HB (Hand Towel & Overwrap) | 378      |

**4200 SERIES****Surgical Gown, Film-Reinforced, Set-In Sleeves**

- ▶ SMS Barrier Fabric
- ▶ Set-In Sleeves
- ▶ Adjustable Comfort-Soft Hook-and-Loop Closures
- ▶ Film Reinforcement
- ▶ Sleeve-N-Sleeve Construction
- ▶ Extra Long, Extra Wide Film Reinforcement

## BASIC


| Size | <b>S</b> Sterile | Case Qty | <b>N</b> Non-Sterile | Case Qty | <b>N</b> Non-Sterile, Handi-Bin | Case Qty |
|------|------------------|----------|----------------------|----------|---------------------------------|----------|
| S    | 90-4200-S1       | 22       | 90-4200              | 40       |                                 |          |
| M    | 90-4210-S1       | 26       | 90-4210              | 40       |                                 |          |
| L    | 90-4220-S1       | 24       | 90-4220              | 36       | 90-4220-HB                      | 378      |
| XL   | 90-4230-S1       | 24       | 90-4230              | 32       | 90-4230-HB                      | 324      |
| XXL  | 90-4240-S1       | 22       | 90-4240              | 32       | 90-4240-HB                      | 324      |

## SPECIALTY

| Size       | <b>S</b> Sterile | Case Qty | <b>N</b> Non-Sterile | Case Qty | <b>N</b> Non-Sterile, Handi-Bin | Case Qty |
|------------|------------------|----------|----------------------|----------|---------------------------------|----------|
| L, XLong   | 90-4222-S1       | 22       | 90-4222              | 32       |                                 |          |
| XL, XLong  | 90-4232-S1       | 20       | 90-4232              | 32       |                                 |          |
| XXL, XLong | 90-4242-S1       | 20       | 90-4242              | 28       |                                 |          |

## SCRUB NURSE

| Size | <b>N</b> Non-Sterile, Handi-Bin       | Case Qty |
|------|---------------------------------------|----------|
| L    | 90-4220-W1-HB (Hand Towel & Overwrap) | 378      |
| XL   | 90-4230-W1-HB (Hand Towel & Overwrap) | 324      |
| XXL  | 90-4240-W1-HB (Hand Towel & Overwrap) | 324      |



**4300 SERIES**  
**Surgical Gown, Fabric-Reinforced, Raglan Sleeves**


- ▶ SMS Barrier Fabric
- ▶ Raglan Sleeves
- ▶ Adjustable Comfort-Soft Hook-and-Loop Closures
- ▶ Fabric Reinforcement
- ▶ Sleeve-N-Sleeve Construction

BASIC

| Size |  Sterile | Case Qty |  Non-Sterile | Case Qty |  Non-Sterile, Handi-Bin | Case Qty |
|------|-------------------------------------------------------------------------------------------|----------|-----------------------------------------------------------------------------------------------|----------|------------------------------------------------------------------------------------------------------------|----------|
| L    | 91-4320-S1                                                                                | 24       |                                                                                               |          | 91-4320-HB                                                                                                 | 378      |
| XL   | 91-4330-S1                                                                                | 22       |                                                                                               |          | 91-4330-HB                                                                                                 | 324      |
| XXL  | 91-4340-S1                                                                                | 22       |                                                                                               |          | 91-4340-HB                                                                                                 | 324      |

SCRUB NURSE

| Size |  Non-Sterile, Handi-Bin | Case Qty |
|------|------------------------------------------------------------------------------------------------------------|----------|
| L    | 91-4320-W1-HB (Hand Towel & Overwrap)                                                                      | 324      |
| XL   | 91-4330-W1-HB (Hand Towel & Overwrap)                                                                      | 324      |

**2000 SERIES****Surgical Gown, Non-Reinforced, Set-In Sleeves**


- ▶ SMS Barrier Fabric
- ▶ Set-In Sleeves
- ▶ Adjustable Comfort-Soft Hook-and-Loop Closures
- ▶ Non-Reinforced

## BASIC

| Size |  Sterile | Case Qty |  Non-Sterile | Case Qty |  Non-Sterile, Handi-Bin | Case Qty |
|------|-------------------------------------------------------------------------------------------|----------|-----------------------------------------------------------------------------------------------|----------|------------------------------------------------------------------------------------------------------------|----------|
| S    | 90-2000-S                                                                                 | 20       | 90-2000                                                                                       | 48       |                                                                                                            |          |
|      | 90-2000-S1 (1 Hand Towel)                                                                 | 20       |                                                                                               |          |                                                                                                            |          |
|      | 90-2000-S2 (2 Hand Towels)                                                                | 20       |                                                                                               |          |                                                                                                            |          |
| M    | 90-2010-S                                                                                 | 20       | 90-2010                                                                                       | 48       |                                                                                                            |          |
|      | 90-2010-S1 (1 Hand Towel)                                                                 | 20       |                                                                                               |          |                                                                                                            |          |
|      | 90-2010-S2 (2 Hand Towels)                                                                | 20       |                                                                                               |          |                                                                                                            |          |
| L    | 90-2020-S                                                                                 | 24       | 90-2020                                                                                       | 48       | 90-2020-HB                                                                                                 | 486      |
|      | 90-2020-S1 (1 Hand Towel)                                                                 | 24       |                                                                                               |          |                                                                                                            |          |
|      | 90-2020-S2 (2 Hand Towels)                                                                | 24       |                                                                                               |          |                                                                                                            |          |
| XL   | 90-2030-S                                                                                 | 22       | 90-2030                                                                                       | 48       | 90-2030-HB                                                                                                 | 486      |
|      | 90-2030-S1 (1 Hand Towel)                                                                 | 22       |                                                                                               |          |                                                                                                            |          |
|      | 90-2030-S2 (2 Hand Towels)                                                                | 22       |                                                                                               |          |                                                                                                            |          |
| XXL  | 90-2040-S                                                                                 | 24       | 90-2040                                                                                       | 44       |                                                                                                            |          |
|      | 90-2040-S1 (1 Hand Towel)                                                                 | 24       |                                                                                               |          |                                                                                                            |          |
|      | 90-2040-S2 (2 Hand Towels)                                                                | 24       |                                                                                               |          |                                                                                                            |          |

## SPECIALTY

| Size       |  Sterile | Case Qty |  Non-Sterile | Case Qty |  Non-Sterile, Handi-Bin | Case Qty |
|------------|---------------------------------------------------------------------------------------------|----------|-------------------------------------------------------------------------------------------------|----------|--------------------------------------------------------------------------------------------------------------|----------|
| L, XLong   | 90-2021-S                                                                                   | 24       | 90-2021                                                                                         | 48       | 90-2021-HB                                                                                                   | 486      |
|            | 90-2021-S1 (1 Hand Towel)                                                                   | 24       |                                                                                                 |          |                                                                                                              |          |
|            | 90-2021-S2 (2 Hand Towels)                                                                  | 24       |                                                                                                 |          |                                                                                                              |          |
| XL, XLong  | 90-2031-S                                                                                   | 32       | 90-2031                                                                                         | 48       | 90-2031-HB                                                                                                   | 486      |
|            | 90-2031-S1 (1 Hand Towel)                                                                   | 32       |                                                                                                 |          |                                                                                                              |          |
|            | 90-2031-S2 (2 Hand Towels)                                                                  | 32       |                                                                                                 |          |                                                                                                              |          |
| XXL, XLong | 90-2041-S                                                                                   | 24       | 90-2041                                                                                         | 40       |                                                                                                              |          |
|            | 90-2041-S1 (1 Hand Towel)                                                                   | 24       |                                                                                                 |          |                                                                                                              |          |
|            | 90-2041-S2 (2 Hand Towels)                                                                  | 24       |                                                                                                 |          |                                                                                                              |          |



## 2200 SERIES

### Surgical Gown, Film-Reinforced, Set-In Sleeves


- ▶ SMS Barrier Fabric
- ▶ Set-In Sleeves
- ▶ Adjustable Comfort-Soft Hook-and-Loop Closures
- ▶ Film Reinforced
- ▶ Sleeve-N-Sleeve Construction
- ▶ Extra Wide Film Reinforcement

#### BASIC

| Size |  Sterile | Case Qty |  Non-Sterile | Case Qty |  Non-Sterile, Handi-Bin | Case Qty |
|------|-------------------------------------------------------------------------------------------|----------|-----------------------------------------------------------------------------------------------|----------|------------------------------------------------------------------------------------------------------------|----------|
| S    | 90-2200-S                                                                                 | 18       | 90-2200                                                                                       | 44       |                                                                                                            |          |
|      | 90-2200-S1 (1 Hand Towel)                                                                 | 18       |                                                                                               |          |                                                                                                            |          |
|      | 90-2200-S2 (2 Hand Towels)                                                                | 18       |                                                                                               |          |                                                                                                            |          |
| M    | 90-2210-S                                                                                 | 18       | 90-2210                                                                                       | 40       |                                                                                                            |          |
|      | 90-2210-S1 (1 Hand Towel)                                                                 | 18       |                                                                                               |          |                                                                                                            |          |
|      | 90-2210-S2 (2 Hand Towels)                                                                | 18       |                                                                                               |          |                                                                                                            |          |
| L    | 90-2220-S                                                                                 | 22       | 90-2220                                                                                       | 44       | 90-2220-HB                                                                                                 | 432      |
|      | 90-2220-S1 (1 Hand Towel)                                                                 | 22       |                                                                                               |          |                                                                                                            |          |
|      | 90-2220-S2 (2 Hand Towels)                                                                | 22       |                                                                                               |          |                                                                                                            |          |
| XL   | 90-2230-S                                                                                 | 22       | 90-2230                                                                                       | 44       | 90-2230-HB                                                                                                 | 432      |
|      | 90-2230-S1 (1 Hand Towel)                                                                 | 22       |                                                                                               |          |                                                                                                            |          |
|      | 90-2230-S2 (2 Hand Towels)                                                                | 22       |                                                                                               |          |                                                                                                            |          |
| XXL  | 90-2240-S                                                                                 | 22       | 90-2240                                                                                       | 36       |                                                                                                            |          |
|      | 90-2240-S1 (1 Hand Towel)                                                                 | 22       |                                                                                               |          |                                                                                                            |          |
|      | 90-2240-S2 (2 Hand Towels)                                                                | 22       |                                                                                               |          |                                                                                                            |          |

#### SPECIALTY

| Size       |  Sterile | Case Qty |  Non-Sterile | Case Qty |  Non-Sterile, Handi-Bin | Case Qty |
|------------|---------------------------------------------------------------------------------------------|----------|-------------------------------------------------------------------------------------------------|----------|--------------------------------------------------------------------------------------------------------------|----------|
| L, XLong   | 90-2221-S                                                                                   | 22       | 90-2221                                                                                         | 44       | 90-2221-HB                                                                                                   | 432      |
|            | 90-2221-S1 (1 Hand Towel)                                                                   | 22       |                                                                                                 |          |                                                                                                              |          |
|            | 90-2221-S2 (2 Hand Towels)                                                                  | 22       |                                                                                                 |          |                                                                                                              |          |
| XL, XLong  | 90-2231-S                                                                                   | 22       | 90-2231                                                                                         | 44       | 90-2231-HB                                                                                                   | 432      |
|            | 90-2231-S1 (1 Hand Towel)                                                                   | 22       |                                                                                                 |          |                                                                                                              |          |
|            | 90-2231-S2 (2 Hand Towels)                                                                  | 22       |                                                                                                 |          |                                                                                                              |          |
| XXL, XLong | 90-2241-S                                                                                   | 22       | 90-2241                                                                                         | 36       |                                                                                                              |          |
|            | 90-2241-S1 (1 Hand Towel)                                                                   | 22       |                                                                                                 |          |                                                                                                              |          |
|            | 90-2241-S2 (2 Hand Towels)                                                                  | 22       |                                                                                                 |          |                                                                                                              |          |

**2300 SERIES****Surgical Gown, Fabric-Reinforced, Set-In Sleeves**

- ▶ SMS Barrier Fabric
- ▶ Set-In Sleeves
- ▶ Adjustable Comfort-Soft Hook-and-Loop Closures
- ▶ Fabric Reinforced
- ▶ Sleeve-N-Sleeve Construction
- ▶ Extra Wide Reinforcement

## BASIC

| Size |  Sterile | Case Qty |  Non-Sterile | Case Qty |  Non-Sterile, Handi-Bin | Case Qty |
|------|-------------------------------------------------------------------------------------------|----------|-----------------------------------------------------------------------------------------------|----------|------------------------------------------------------------------------------------------------------------|----------|
| S    | 90-2300-S1 (1 Hand Towel)                                                                 | 18       | 90-2300                                                                                       | 44       |                                                                                                            |          |
|      | 90-2300-S2 (2 Hand Towels)                                                                | 18       |                                                                                               |          |                                                                                                            |          |
| M    | 90-2310-S1 (1 Hand Towel)                                                                 | 18       | 90-2310                                                                                       | 44       |                                                                                                            |          |
|      | 90-2310-S2 (2 Hand Towels)                                                                | 18       |                                                                                               |          |                                                                                                            |          |
| L    | 90-2320-S1 (1 Hand Towel)                                                                 | 22       | 90-2320                                                                                       | 44       |                                                                                                            |          |
|      | 90-2320-S2 (2 Hand Towels)                                                                | 22       |                                                                                               |          |                                                                                                            |          |
| XL   | 90-2330-S1 (1 Hand Towel)                                                                 | 22       | 90-2330                                                                                       | 44       |                                                                                                            |          |
|      | 90-2330-S2 (2 Hand Towels)                                                                | 22       |                                                                                               |          |                                                                                                            |          |
| XXL  | 90-2340-S1 (1 Hand Towel)                                                                 | 20       | 90-2340                                                                                       | 36       |                                                                                                            |          |
|      | 90-2340-S2 (2 Hand Towels)                                                                | 20       |                                                                                               |          |                                                                                                            |          |

## SPECIALTY

| Size       |  Sterile | Case Qty |  Non-Sterile | Case Qty |  Non-Sterile, Handi-Bin | Case Qty |
|------------|---------------------------------------------------------------------------------------------|----------|-------------------------------------------------------------------------------------------------|----------|--------------------------------------------------------------------------------------------------------------|----------|
| L, XLong   | 90-2321-S1 (1 Hand Towel)                                                                   | 22       | 90-2321                                                                                         | 44       |                                                                                                              |          |
|            | 90-2321-S2 (2 Hand Towels)                                                                  | 22       |                                                                                                 |          |                                                                                                              |          |
| XL, XLong  | 90-2331-S1 (1 Hand Towel)                                                                   | 22       | 90-2331                                                                                         | 44       |                                                                                                              |          |
|            | 90-2331-S2 (2 Hand Towels)                                                                  | 22       |                                                                                                 |          |                                                                                                              |          |
| XXL, XLong | 90-2341-S1 (1 Hand Towel)                                                                   | 20       | 90-2341                                                                                         | 36       |                                                                                                              |          |
|            | 90-2341-S2 (2 Hand Towels)                                                                  | 20       |                                                                                                 |          |                                                                                                              |          |



# SURGICAL DRAPES

Reinvention of the surgical drape means that every millimeter presents an opportunity to improve performance and in turn, improve success during surgery. GRI-Alleset's SteriGuard® Surgical Drapes are meticulously designed to provide optimal ease of use and unmatched surgical site isolation. Design customization capabilities provide flexibility to fulfill the complex needs of every procedure.



## PROTECTION

Advanced non-woven fabric technologies are leveraged to deliver excellent barrier, strength, durability, and adhesion properties, protecting against infection. Low lint properties help eliminate one source of post-operative patient complications. Flame resistant and anti-static treatment help prevent potential source ignition.



## SURGICAL EDGING

GRI-Alleset's EdgeGuard technology provides adhesive edging to maintain surgical site isolation and ensure drapes stay in place during surgery. 2 in / 5 cm medical tape provides improved drape placement and drape security.



## ABSORBENCY

GRI-Alleset's HydraGuard® advanced absorbency and barrier fabric provides superior management of blood and bodily fluids. Precision placement and integration of HydraGuard® materials ensures that high-absorbency needed areas are fully covered.



## EASE OF USE

Intuitive designs with unfolding instructions clearly stamped on each drape provide assistance in the draping process. Soft and pliable fabrics and adhesive materials increase ease of application and improve drapability.



DELIVERING PRODUCT SOLUTIONS™

# GRI-Alleset's SteriGuard® Surgical Drapes set a new mark for the best in strength, protection, absorbency, and ease of use.

GRI-Alleset's SteriGuard® Surgical Drapes are a breakthrough in design, delivering optimal ease of use, unmatched protection and surgical site isolation for creating and maintaining a sterile field. No matter the procedure type, SteriGuard® combines superior barrier definition, absorption, and strength, to ensure that placement is easy and that drapes remain in place during surgery.

Combining advanced fabric technologies and unique manufacturing processes with strategic, intuitive design is helping SteriGuard® redefine the meaning of quality in the market. Making drapes with the best cost-to-quality value available, GRI-Alleset is taking Surgical Drapes to the next level.



## Advanced SMS Protects Patients From Infection

Superior protection in the critical zones provides excellent barrier protection from fluid strikethrough and minimizes the risk of cross-contamination and infection.



## Low Lint Technology

Each drape utilizes a low linting, three-layer polypropylene non-woven fabric to minimize lint particle generation, which can contribute to bacterial transmission or foreign body contamination.



## 2 in / 5 cm Tape

SteriGuard® utilizes 2 in / 5 cm medical tape for improved drape placement and drape security.



## Precision Adhesive

GRI-Alleset's EdgeGuard technology allows medical-grade adhesion to be applied to the very edge of the fenestration, ensuring precise placement, security of the surgical drape, and fluid management.



## Strong and Durable SMS

GRI-Alleset surgical fabric is constructed with three layers of SMS, making it highly resistant to tears and abrasion.



## Advanced Fluid Management

GRI-Alleset's HydraGuard® reinforcement fabric provides high fluid absorbency around the fenestration to promote fluid management. GRI-Alleset surgical drape pouch and pocket designs allow for advanced fluid management and efficient clean up.



## Low Flammability

State-of-the-art fabric technology provides superior resistance to flame protection that can be caused by the use of high-energy medical devices during surgery.

## Easy To Use

Starting with carefully coordinated packaging, SteriGuard® drapes are easily and aseptically opened to manage and maintain a sterile field.






# ANGIOGRAPHY SURGICAL DRAPE



 **STERIGUARD®**

*grī* | alle  <sup>®</sup>

DELIVERING PRODUCT SOLUTIONS™




### Pediatric Femoral Angiography Drape with Two Clear Windows

Size: 76 x 124 in / 193 x 315 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-01139G-S  | <span style="color: red;">S</span> Sterile                 | 12       |
| 80-01139G    | <span style="color: teal;">N</span> Non-Sterile            | N/A      |
| 80-01139G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Two Circular Adhesive Fenestrations: 1.7 in / 4.3 cm each
- ▶ Two Circular Adhesive Fenestrations: 3.5 in / 9 cm each
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Two Clear Plastic Windows: 27 x 124 in / 69 x 315 cm
- ▶ Patient Comfort Layer




### Femoral Angiography Drape with Two Clear Windows and Absorbent Impervious Reinforcement

Size: 75 x 125 in / 191 x 318 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-01144-S  | <span style="color: red;">S</span> Sterile                 | 10       |
| 80-01144    | <span style="color: teal;">N</span> Non-Sterile            | 12       |
| 80-01144-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Two Circular Adhesive Fenestrations: 3.5 in / 8.9 cm each
- ▶ Absorbent Impervious Reinforcement: 32 x 133 in / 81 x 338 cm
- ▶ Two Clear Plastic Windows: 27 x 125 in / 69 x 318 cm
- ▶ Patient Comfort Layer




### Femoral Angiography Drape (Extra Long)

Size: 87 x 135 in / 221 x 343 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-01105-S  | <span style="color: red;">S</span> Sterile                 | 14       |
| 80-01105    | <span style="color: teal;">N</span> Non-Sterile            | 12       |
| 80-01105-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | 200      |

- ▶ Two Circular Adhesive Fenestrations: 3.5 in / 8.9 cm each
- ▶ Absorbent Impervious Reinforcement: 32 x 133 in / 81 x 338 cm
- ▶ Two Clear Plastic Windows: 27 x 135 in / 69 x 343 cm
- ▶ Patient Comfort Layer




### Femoral Angiography Drape (Extra Long) with Two Clear Windows

Size: 87 x 135 in / 221 x 343 cm

| Code        | Summary                                                                                                  | Case Qty |
|-------------|----------------------------------------------------------------------------------------------------------|----------|
| 80-01108-S  |  Sterile                | 14       |
| 80-01108    |  Non-Sterile            | 12       |
| 80-01108-HB |  Non-Sterile, Handi-Bin | 200      |

- ▶ Two Circular Adhesive Fenestrations: 3.5 in / 8.9 cm each
- ▶ Absorbent Impervious Reinforcement: 32 x 133 in / 81 x 338 cm
- ▶ Two Clear Plastic Windows: 27 x 135 in / 69 x 343 cm
- ▶ Extended Absorbent Fabric Reinforcement: 27 x 31 in / 69 x 79 cm
- ▶ Patient Comfort Layer




### Femoral Angiography Drape (Extra Long) with Two Clear Windows and Pouch

Size: 87 x 135 in / 221 x 343 cm

| Code        | Summary                                                                                                    | Case Qty |
|-------------|------------------------------------------------------------------------------------------------------------|----------|
| 80-01151-S  |  Sterile                | 14       |
| 80-01151    |  Non-Sterile            | 12       |
| 80-01151-HB |  Non-Sterile, Handi-Bin | N/A      |

- ▶ Two Circular Adhesive Fenestrations: 3.5 in / 8.9 cm each
- ▶ Absorbent Impervious Reinforcement
- ▶ Fluid Collection Pouch
- ▶ Patient Comfort Layer




### Femoral Angiography Drape (Extra Long) with Two Clear Windows and Two Pouches

Size: 87 x 135 in / 221 x 343 cm

| Code        | Summary                                                                                                    | Case Qty |
|-------------|------------------------------------------------------------------------------------------------------------|----------|
| 80-01152-S  |  Sterile                | 14       |
| 80-01152    |  Non-Sterile            | 12       |
| 80-01152-HB |  Non-Sterile, Handi-Bin | N/A      |

- ▶ Two Circular Adhesive Fenestrations: 3.5 in / 8.9 cm each
- ▶ Absorbent Impervious Reinforcement
- ▶ Two Fluid Collection Pouches
- ▶ Patient Comfort Layer




### Femoral Angiography Drape with Two Clear Windows and Clear Foot Panel

Size: 94 x 160 in / 239 x 406 cm

| Code        | Summary                                                                                                  | Case Qty |
|-------------|----------------------------------------------------------------------------------------------------------|----------|
| 80-01147-S  |  Sterile                | 12       |
| 80-01147    |  Non-Sterile            | 19       |
| 80-01147-HB |  Non-Sterile, Handi-Bin | N/A      |

- ▶ Two Circular Adhesive Fenestrations: 3.5 in / 8.9 cm each
- ▶ Absorbent Impervious Reinforcement: 32 x 133 in / 81 x 338 cm
- ▶ Two Clear Plastic Windows: 31 x 160 in / 79 x 406 cm
- ▶ Clear Plastic Foot Window: 27 x 90 in / 69 x 229 cm
- ▶ Patient Comfort Layer




### Femoral Angiography Drape with Two Clear Windows

Size: 87 x 124 in / 221 x 315 cm

| Code         | Summary                                                                                                    | Case Qty |
|--------------|------------------------------------------------------------------------------------------------------------|----------|
| 80-01103G-S  |  Sterile                | 16       |
| 80-01103G    |  Non-Sterile            | 44       |
| 80-01103G-HB |  Non-Sterile, Handi-Bin | 252      |

- ▶ Two Circular Adhesive Fenestrations: 3.5 in / 8.9 cm each
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Two Clear Plastic Windows: 27 x 124 in / 69 x 315 cm each
- ▶ Patient Comfort Layer




### Femoral Angiography Drape (Extra Long) with Two Clear Windows

Size: 89 x 135 in / 226 x 343 cm

| Code         | Summary                                                                                                    | Case Qty |
|--------------|------------------------------------------------------------------------------------------------------------|----------|
| 80-01106G-S  |  Sterile                | 16       |
| 80-01106G    |  Non-Sterile            | 18       |
| 80-01106G-HB |  Non-Sterile, Handi-Bin | 200      |

- ▶ Two Circular Adhesive Fenestrations: 3.5 in / 8.9 cm each
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Two Clear Plastic Windows: 27 x 135 in / 69 x 343 cm each
- ▶ Patient Comfort Layer




### Femoral Angiography Drape

Size: 89 x 135 in / 225 x 343 cm

| Code         | Summary                                                                                                  | Case Qty |
|--------------|----------------------------------------------------------------------------------------------------------|----------|
| 80-01181G-S  |  Sterile                | N/A      |
| 80-01181G    |  Non-Sterile            | N/A      |
| 80-01181G-HB |  Non-Sterile, Handi-Bin | N/A      |

- ▶ Two Circular Adhesive Fenestrations: 3.5 in / 8.9 cm each
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Clear PE Film




### Femoral Angiography Drape

Size: 89 x 144 in / 225 x 366 cm

| Code         | Summary                                                                                                    | Case Qty |
|--------------|------------------------------------------------------------------------------------------------------------|----------|
| 80-01135G-S  |  Sterile                | N/A      |
| 80-01135G    |  Non-Sterile            | N/A      |
| 80-01135G-HB |  Non-Sterile, Handi-Bin | N/A      |

- ▶ Two Circular Adhesive Fenestrations: 3.5 in / 8.9 cm each
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Clear PE Film




### Femoral Angiography Drape with Absorbent Pouch

Size: 82 x 124 in / 208 x 315 cm

| Code         | Summary                                                                                                    | Case Qty |
|--------------|------------------------------------------------------------------------------------------------------------|----------|
| 80-01131G-S  |  Sterile                | 16       |
| 80-01131G    |  Non-Sterile            | 16       |
| 80-01131G-HB |  Non-Sterile, Handi-Bin | N/A      |

- ▶ Two Circular Adhesive Fenestrations: 5 in / 12.7 cm each
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ One Clear Plastic Window: 27 x 124 in / 69 x 315 cm
- ▶ Patient Comfort Layer




### Femoral Angiography Drape with One Clear Window

Size: 89 x 120 in / 226 x 305 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-01102G-S  | <span style="color: red;">S</span> Sterile                 | 14       |
| 80-01102G    | <span style="color: teal;">N</span> Non-Sterile            | 36       |
| 80-01102G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | 210      |

- ▶ Two Circular Adhesive Fenestrations: 3.5 in / 8.9 cm each
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ One Clear Plastic Window: 27 x 124 in / 89 x 305 cm
- ▶ Patient Comfort Layer




### Femoral Angiography Drape

Size: 89 x 120 in / 226 x 305 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-01101G-S  | <span style="color: red;">S</span> Non-Sterile             | 22       |
| 80-01101G    | <span style="color: teal;">N</span> Non-Sterile            | 28       |
| 80-01101G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | 168      |

- ▶ Two Circular Adhesive Fenestrations: 3.5 in / 8.9 cm each
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Patient Comfort Layer




### Radial Femoral Angiography Drape (Extra Long) with Two Clear Windows

Size: 95 x 150 in / 241 x 381 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-01132G-S  | <span style="color: red;">S</span> Sterile                 | 12       |
| 80-01132G    | <span style="color: teal;">N</span> Non-Sterile            | 14       |
| 80-01132G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Two Adhesive Radial Fenestrations: 4.7 x 3 in / 11.9 x 7.6 cm each
- ▶ Two Oval Adhesive Fenestrations: 5 in / 12.7 cm each
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Two Clear Plastic Windows: 25 x 150 in / 64 x 381 cm each
- ▶ Patient Comfort Layer




### Radial Femoral Angiography Drape with Clear Windows

Size: 95 x 150 in / 241 x 381 cm

| Code        | Summary                                                     | Case Qty |
|-------------|-------------------------------------------------------------|----------|
| 80-01121-S  | <span style="color: red;">S</span> Sterile                  | 10       |
| 80-01121    | <span style="color: green;">N</span> Non-Sterile            | 12       |
| 80-01121-HB | <span style="color: green;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Two Circular Adhesive Fenestrations: 3.5 in / 8.9 cm each
- ▶ Absorbent Impervious Reinforcement: 32 x 133 in / 81 x 338 cm
- ▶ Two Clear Plastic Windows: 31 x 133 in / 79 x 338 cm
- ▶ Clear Plastic Foot Window: 27 x 90 in / 69 x 229 cm
- ▶ Patient Comfort Layer




### Transradial Femoral Angiography Drape

Size: 87 x 135 in / 221 x 343 cm

| Code        | Summary                                                     | Case Qty |
|-------------|-------------------------------------------------------------|----------|
| 80-01166-S  | <span style="color: red;">S</span> Sterile                  | 10       |
| 80-01166    | <span style="color: green;">N</span> Non-Sterile            | 12       |
| 80-01166-HB | <span style="color: green;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Two Adhesive Radial Fenestrations: 4.7 in / 11.9 cm each
- ▶ Two Oval Adhesive Fenestrations: 2.8 x 4.7 in / 7.1 x 11.9 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Two Clear Plastic Windows: 24 x 194 in / 61 x 493 cm
- ▶ Adjustable Trough
- ▶ Hook-and-Loop Tabs




### Multi Access Angiography Drape

Size: 89 x 137 in / 226 x 348 cm

| Code        | Summary                                                     | Case Qty |
|-------------|-------------------------------------------------------------|----------|
| 80-01165-S  | <span style="color: red;">S</span> Sterile                  | 9        |
| 80-01165    | <span style="color: green;">N</span> Non-Sterile            | 10       |
| 80-01165-HB | <span style="color: green;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Two Circular Femoral Adhesive Fenestrations: 4.7 in / 12 cm each
- ▶ Two Radial Adhesive Fenestrations: 2.7 x 4.7 in / 7 x 12 cm each
- ▶ Two Clear Plastic Windows: 24 x 173 in / 60 x 439 cm
- ▶ Adjustable Absorbent Pocket
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Patient Comfort Layer




### Radial Femoral Angiography Drape

Size: 94 x 148 in / 239 x 376 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-01140G-S  | <span style="color: red;">S</span> Sterile                 | 12       |
| 80-01140G    | <span style="color: teal;">N</span> Non-Sterile            | 12       |
| 80-01140G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Two Circular Adhesive Fenestrations: 3.5 in / 8.9 cm each
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ One Clear Plastic Window: 27 x 124 in / 89 x 305 cm
- ▶ Patient Comfort Layer




### Radial Femoral Angiography Drape (Extra Long) with Two Clear Windows

Size: 94 x 173 in / 239 x 439 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-01141G-S  | <span style="color: red;">S</span> Sterile                 | 10       |
| 80-01141G    | <span style="color: teal;">N</span> Non-Sterile            | 10       |
| 80-01141G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Two Circular Adhesive Fenestrations: 3.5 in / 8.9 cm each
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Patient Comfort Layer



### Radial Femoral Angiography Drape with Two Clear Windows, Absorbent Adjustable Pouch

Size: 96 x 194 in / 244 x 493 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-01143G-S  | <span style="color: red;">S</span> Sterile                 | 10       |
| 80-01143G    | <span style="color: teal;">N</span> Non-Sterile            | 10       |
| 80-01143G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Two Adhesive Radial Fenestrations: 4.7 in / 11.9 cm each
- ▶ Two Oval Adhesive Fenestrations: 2.8 x 4.7 in / 7.1 x 11.9 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Two Clear Plastic Windows: 24 x 194 in / 61 x 493 cm
- ▶ Adjustable Trough
- ▶ Hook-and-Loop Tabs




# SPECIAL PROCEDURES SURGICAL DRAPE



*grī* | alle **set**®

DELIVERING PRODUCT SOLUTIONS™




### Brachial Angiography Drape

Size: 38 x 44 in / 97 x 112 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-01164G-S  | <span style="color: red;">S</span> Sterile                 | 32       |
| 80-01164G    | <span style="color: teal;">N</span> Non-Sterile            | 38       |
| 80-01164G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Circular Adhesive Fenestration: 3 in / 7.6 cm
- ▶ HydraGuard® Absorbent Reinforcement




### Brachial Angiography Drape with Fluid Collection Pouch

Size: 46 x 52.5 in / 117 x 133 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-01163-S  | <span style="color: red;">S</span> Sterile                 | 38       |
| 80-01163    | <span style="color: teal;">N</span> Non-Sterile            | 50       |
| 80-01163-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Oval Adhesive Fenestration: 4.5 x 3.5 in / 11.4 x 8.9
- ▶ Clear Window for checking the patient's fingertips for circulation
- ▶ U-Shaped Fluid Collection Pouch




### Brachial Angiography Drape with Absorbent Impervious Reinforcement

Size: 30 x 42 in / 76 x 107 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-01148-S  | <span style="color: red;">S</span> Sterile                 | 38       |
| 80-01148    | <span style="color: teal;">N</span> Non-Sterile            | 50       |
| 80-01148-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Oval Fenestration: 4.5 x 3.5 in / 11.4 x 8.9 cm
- ▶ Absorbent Impervious Reinforcement
- ▶ Patient Comfort Layer



### Pacemaker Drape


Size: 87 x 135 in / 221 x 343 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-02109G-S  | <span style="color: red;">S</span> Sterile                 | 16       |
| 80-02109G    | <span style="color: teal;">N</span> Non-Sterile            | 20       |
| 80-02109G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Two Bilateral Aperture Fenestrations: 6 x 8 in / 15.2 x 20.3 cm
- ▶ Bilateral Pockets
- ▶ Adhesive Tape with Release Liner
- ▶ Adjustable Breather Bar
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Two Clear Plastic Windows

## SPECIAL PROCEDURES

### Surgical Drapes




#### Pacemaker Drape

Size: 91 x 135 in / 231 x 343 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-02113G-S  | <span style="color: red;">S</span> Sterile                 | 14       |
| 80-02113G    | <span style="color: teal;">N</span> Non-Sterile            | 18       |
| 80-02113G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Two Bilateral Aperture Fenestrations: 3 x 5 in / 7.6 x 12.7 cm
- ▶ Absorbent Bilateral Pockets
- ▶ Adhesive Incise with Release Liner
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Two Clear Plastic Windows




#### Pacemaker Drape

Size: 91 x 150 in / 231 x 381 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-02114G-S  | <span style="color: red;">S</span> Sterile                 | 14       |
| 80-02114G    | <span style="color: teal;">N</span> Non-Sterile            | 18       |
| 80-02114G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Two Bilateral Aperture Fenestrations: 3 x 5 in / 7.6 x 12.7 cm
- ▶ Absorbent Bilateral Pockets
- ▶ Adhesive Incise with Release Liner
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Two Clear Plastic Windows




#### Pacemaker Drape

Size: 101 x 79 x 145 in / 256 x 200 x 368 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-02120G-S  | <span style="color: red;">S</span> Sterile                 | 8        |
| 80-02120G    | <span style="color: teal;">N</span> Non-Sterile            | 16       |
| 80-02120G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Incise Fenestration: 15 x 23 in / 38.1 x 58.4 cm
- ▶ Absorbent Troughs
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Hook and loop Line Control
- ▶ Clear Instrumentation Window




#### Peripherally Inserted Central Catheter (PICC) Drape

Size: 78 x 111 in / 198 x 282 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-01149G-S  | <span style="color: red;">S</span> Sterile                 | 20       |
| 80-01149G    | <span style="color: teal;">N</span> Non-Sterile            | 24       |
| 80-01149G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Two Fenestrations: 6 x 5 in / 15.2 x 12.7 cm
- ▶ Adhesive Tape with Release Liner
- ▶ HydraGuard® Absorbent Reinforcement




### TAVI Drape (Extra Long) with Two Clear Windows and Absorbent Impervious Reinforcement

Size: 92 x 180 in / 234 x 457 cm

| Code         | Summary                                                                                                  | Case Qty |
|--------------|----------------------------------------------------------------------------------------------------------|----------|
| 80-02107G-S  |  Sterile                | 8        |
| 80-02107G    |  Non-Sterile            | 8        |
| 80-02107G-HB |  Non-Sterile, Handi-Bin | N/A      |

- ▶ Transcatheter Aortic Valve Replacement Drape
- ▶ One Rectangular Adhesive Fenestration: 15 x 12 in / 38.1 x 33 cm each
- ▶ Two Oval Adhesive Fenestrations: 5 in / 12.7 cm each
- ▶ Two Clear Plastic Windows
- ▶ Instrument Pouch



### Botto Transpedal Access Draping System

| Code        | Summary                                                                                                  | Case Qty |
|-------------|----------------------------------------------------------------------------------------------------------|----------|
| 86-01607-S  |  Sterile                | 5        |
| 86-01607    |  Non-Sterile            | N/A      |
| 86-01607-HB |  Non-Sterile, Handi-Bin | N/A      |

- 1 CV Split Drape with 2 in Positioning Tape (80-02125G)
  - ▶ 115 x 85 x 114 in / 292 x 216 x 290 cm
  - ▶ Longitudinal Line Control Flaps
  - ▶ Clear Film Anesthesia
  - ▶ HydraGuard® Absorbent Reinforcement

- 1 Groin Drape (80-02101G)
  - ▶ 84 x 80 in / 213 x 203 cm
  - ▶ HydraGuard® Absorbent Reinforcement
  - ▶ Impervious Perineal Cover
  - ▶ Umbilicus Cover



# BASIC DRAPES SURGICAL DRAPES



 **STERIGUARD**

*grī* | alle **set**

DELIVERING PRODUCT SOLUTIONS™




### Half Drape

Size: 44 x 60 in / 112 x 152 cm

| Code        | Summary                                                                                                  | Case Qty |
|-------------|----------------------------------------------------------------------------------------------------------|----------|
| 80-12401-S  |  Sterile                | 22       |
| 80-12401    |  Non-Sterile            | 74       |
| 80-12401-HB |  Non-Sterile, Handi-Bin | 960      |

► SMS Fabric




### Medium Drape

Size: 44 x 76 in / 112 x 193 cm

| Code        | Summary                                                                                                  | Case Qty |
|-------------|----------------------------------------------------------------------------------------------------------|----------|
| 80-12402-S  |  Sterile                | 40       |
| 80-12402    |  Non-Sterile            | 58       |
| 80-12402-HB |  Non-Sterile, Handi-Bin | N/A      |

► SMS Fabric



### Medium Drape with Tape

Size: 44 x 76 in / 112 x 193 cm


| Code        | Summary                                                                                                    | Case Qty |
|-------------|------------------------------------------------------------------------------------------------------------|----------|
| 80-12406-S  |  Sterile                | N/A      |
| 80-12406    |  Non-Sterile            | 50       |
| 80-12406-HB |  Non-Sterile, Handi-Bin | N/A      |

► SMS Fabric

► 2 in Positioning Tape

## BASIC DRAPES

### Surgical Drapes




#### Large Drape

Size: 60 x 76 in / 152 x 193 cm

| Code        | Summary                                                                                                  | Case Qty |
|-------------|----------------------------------------------------------------------------------------------------------|----------|
| 80-12403-S  |  Sterile                | 20       |
| 80-12403    |  Non-Sterile            | 24       |
| 80-12403-HB |  Non-Sterile, Handi-Bin | 525      |

► SMS Fabric



#### Large Drape, Reverse Fold

Size: 60 x 76 in / 152 x 193 cm

| Code        | Summary                                                                                                  | Case Qty |
|-------------|----------------------------------------------------------------------------------------------------------|----------|
| 80-12408-S  |  Sterile                | 24       |
| 80-12408    |  Non-Sterile            | 24       |
| 80-12408-HB |  Non-Sterile, Handi-Bin | N/A      |

► SMS Fabric




#### X-Large Drape

Size: 76 x 100 in / 193 x 254 cm

| Code        | Summary                                                                                                    | Case Qty |
|-------------|------------------------------------------------------------------------------------------------------------|----------|
| 80-12404-S  |  Sterile                | 28       |
| 80-12404    |  Non-Sterile            | 30       |
| 80-12404-HB |  Non-Sterile, Handi-Bin | N/A      |

► SMS Fabric




### Utility Drape with Tape

Size: 15 x 26 in / 38 x 66 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-05126-S  | <span style="color: red;">S</span> Sterile                 | 200      |
| 80-05126    | <span style="color: teal;">N</span> Non-Sterile            | 180      |
| 80-05126-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | 4800     |

► SMS Fabric




### Utility Drape

Size: 30 x 35 in / 75 x 90 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-05803-S  | <span style="color: red;">S</span> Sterile                 | N/A      |
| 80-05803    | <span style="color: teal;">N</span> Non-Sterile            | N/A      |
| 80-05803-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

► SMS Fabric



### Three-Quarter Drape

Size: 60 x 76 in / 152 x 193 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-12110G-S  | <span style="color: red;">S</span> Sterile                 | 30       |
| 80-12110G    | <span style="color: teal;">N</span> Non-Sterile            | 38       |
| 80-12110G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | 192      |


► HydraGuard® Absorbent Reinforcement



# CARDIOVASCULAR SURGICAL DRAPE



DELIVERING PRODUCT SOLUTIONS™




### TECH70 Cardiovascular Split Drape

Size: 108 x 80 x 126 in / 274 x 203 x 320 cm

| Code         | Summary                | Case Qty |
|--------------|------------------------|----------|
| 80-02102G-S  | Sterile                | 10       |
| 80-02102G    | Non-Sterile            | 8        |
| 80-02102G-HB | Non-Sterile, Handi-Bin | N/A      |



- ▶ TECH70 Performance SMS Fabric
- ▶ Split: 6 x 84 in / 15.2 x 213.4 cm
- ▶ Longitudinal Pockets with Hook-and-Loop Closure
- ▶ 2 in / 5 cm Adhesive Tape
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Armboard Covers
- ▶ Tube Holders




### TECH70 Cardiovascular Split Drape with Clear Anesthesia Screen

Size: 115 x 85 x 151 in / 292 x 216 x 384 cm

| Code         | Summary                | Case Qty |
|--------------|------------------------|----------|
| 80-02110G-S  | Sterile                | 6        |
| 80-02110G    | Non-Sterile            | 7        |
| 80-02110G-HB | Non-Sterile, Handi-Bin | N/A      |



- ▶ TECH70 Performance SMS Fabric
- ▶ Split: 4 x 7 x 97 in / 10.2 x 17.8 x 246.4 cm
- ▶ 2 in / 5 cm Adhesive Tape
- ▶ Clear Anesthesia Screen
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Armboard Covers




### Cardiovascular Drape

Size: 108 x 126 in / 275 x 320 cm

| Code         | Summary                | Case Qty |
|--------------|------------------------|----------|
| 80-02103G-S  | Sterile                | 10       |
| 80-02103G    | Non-Sterile            | 10       |
| 80-02103G-HB | Non-Sterile, Handi-Bin | N/A      |

- ▶ Chest Adhesive Incise: 14 x 17 in / 35.6 x 43.2 cm
- ▶ Leg Adhesive Incise: 20 x 41 in / 50 x 103 cm
- ▶ Longitudinal Pockets with Hook-and-Loop Closure
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Armboard Covers
- ▶ Tube Holders




### Cardiovascular Incise Drape

Size: 110 x 80 x 128 in / 279 x 203 x 325 cm

| Code         | Summary                | Case Qty |
|--------------|------------------------|----------|
| 80-02136G-S  | Sterile                | N/A      |
| 80-02136G    | Non-Sterile            | 6        |
| 80-02136G-HB | Non-Sterile, Handi-Bin | N/A      |

- ▶ Chest Adhesive Incise: 14 x 17 in / 35.6 x 43.2 cm
- ▶ Leg Adhesive Incise: 20 x 41 in / 50.8 x 104.1 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders
- ▶ Armboard Covers
- ▶ Longitudinal Flaps






### Bilateral Underleg Perineal Drape

Size: 84 x 80 in / 213 x 203 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-02101G-S  | <span style="color: red;">S</span> Sterile                 | 24       |
| 80-02101G    | <span style="color: teal;">N</span> Non-Sterile            | 58       |
| 80-02101G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Impervious Perineal Cover
- ▶ Umbilicus Cover




### Thoracic Drape

Size: 94 x 132 in / 239 x 335 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-02108G-S  | <span style="color: red;">S</span> Sterile                 | N/A      |
| 80-02108G    | <span style="color: teal;">N</span> Non-Sterile            | 14       |
| 80-02108G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Fenestration: 20 x 16 in / 50 x 40 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders




### Chest Drape

Size: 100 x 72 x 126 in / 254 x 183 x 320 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-05122G-S  | <span style="color: red;">S</span> Sterile                 | 12       |
| 80-05122G    | <span style="color: teal;">N</span> Non-Sterile            | 14       |
| 80-05122G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Fenestration: 13 x 14 in / 33 x 35.6 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Instrument Pad
- ▶ Tube Holders
- ▶ Armboard Covers




### TAVI Drape (Extra Long) with Two Clear Windows and Absorbent Impervious Reinforcement

Size: 92 x 180 in / 234 x 457 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-02107G-S  | <span style="color: red;">S</span> Sterile                 | N/A      |
| 80-02107G    | <span style="color: teal;">N</span> Non-Sterile            | N/A      |
| 80-02107G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Transcatheter Aortic Valve Replacement Drape
- ▶ Rectangular Incise Fenestration: 15 x 12 in / 38.1 x 30.5 cm
- ▶ Two Oval Adhesive Fenestrations: 5 in / 12.7 cm each
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Two Clear Plastic Windows
- ▶ Instrument Pouch




### Cardiovascular Split Drape with Clear Anesthesia Screen

Size: 109 x 84.5 x 109 in / 277 x 215 x 277 cm

| Code         | Summary                                                                                                  | Case Qty |
|--------------|----------------------------------------------------------------------------------------------------------|----------|
| 80-01205G-S  |  Sterile                | 7        |
| 80-01205G    |  Non-Sterile            | 7        |
| 80-01205G-HB |  Non-Sterile, Handi-Bin | N/A      |

- ▶ Split: 6 x 84 in / 15.2 x 213.4 cm
- ▶ Longitudinal Pockets with Hook-and-Loop Closure
- ▶ 2 in / 5 cm Adhesive Tape
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Armboard Covers
- ▶ Tube Holders



### Heavy-Duty Universal Draping System

| Code        | Summary                                                                                                  | Case Qty |
|-------------|----------------------------------------------------------------------------------------------------------|----------|
| 86-01603-S  |  Sterile                | 5        |
| 86-01603    |  Non-Sterile            | 5        |
| 86-01603-HB |  Non-Sterile, Handi-Bin | N/A      |

- ▶ TECH70 Performance SMS Fabric
- ▶ 1 Top Drape, 100 x 60 in / 254 x 152 cm
- ▶ 1 Bottom Drape, 76 x 60 in / 193 x 152 cm
- ▶ 2 Side Drapes, 38 x 84 in / 97 x 213 cm
- ▶ Longitudinal Flaps
- ▶ HydraGuard® Absorbent Reinforcement

**Tech70**  
PERFORMANCE




# EENT / NEUROLOGY SURGICAL DRAPE



 **STERIGUARD**

*grī* | alle 

DELIVERING PRODUCT SOLUTIONS™




### EENT Split Drape

Size: 76 x 124 in / 193 x 315 cm

| Code         | Summary                                                     | Case Qty |
|--------------|-------------------------------------------------------------|----------|
| 80-09143G-S  | <span style="color: red;">S</span> Sterile                  | 22       |
| 80-09143G    | <span style="color: green;">N</span> Non-Sterile            | 20       |
| 80-09143G-HB | <span style="color: green;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Split: 3 x 25 in / 7.6 x 63.5 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders




### EENT Split Drape

Size: 76 x 124 in / 193 x 315 cm

| Code         | Summary                                                     | Case Qty |
|--------------|-------------------------------------------------------------|----------|
| 80-09144G-S  | <span style="color: red;">S</span> Sterile                  | 26       |
| 80-09144G    | <span style="color: green;">N</span> Non-Sterile            | 30       |
| 80-09144G-HB | <span style="color: green;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Split: 3 x 28 in / 7.6 x 71.1 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Instrument Pad
- ▶ Tube Holders




### EENT Split Drape

Size: 76 x 124 in / 193 x 315 cm

| Code         | Summary                                                     | Case Qty |
|--------------|-------------------------------------------------------------|----------|
| 80-09142G-S  | <span style="color: red;">S</span> Sterile                  | 16       |
| 80-09142G    | <span style="color: green;">N</span> Non-Sterile            | 36       |
| 80-09142G-HB | <span style="color: green;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Split: 3 x 31 in / 7.6 x 78.7 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders




### Split Sheet Drape

Size: 76 x 120 in / 193 x 305 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-09118G-S  | <span style="color: red;">S</span> Sterile                 | 14       |
| 80-09118G    | <span style="color: teal;">N</span> Non-Sterile            | 24       |
| 80-09118G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Split: 4 x 40 in / 10 x 102 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders




### Head Drape without Turban

Size: 50 x 36 in / 127 x 91 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-04104-S  | <span style="color: red;">S</span> Sterile                 | 42       |
| 80-04104    | <span style="color: teal;">N</span> Non-Sterile            | 60       |
| 80-04104-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Turban: 33 x 20 in / 83.8 x 50.8 cm




### Head Drape with Turban

Size: 41 x 32 in / 104 x 81 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-04105-S  | <span style="color: red;">S</span> Sterile                 | 40       |
| 80-04105    | <span style="color: teal;">N</span> Non-Sterile            | 56       |
| 80-04105-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Detachable Turban: 36 x 24 in / 91 x 61 cm



### Craniotomy Drape with Pouch

Size: 122 x 74 x 134 in / 310 x 188 x 340 cm

| Code         | Summary                                                                                                  | Case Qty |
|--------------|----------------------------------------------------------------------------------------------------------|----------|
| 80-06001G-S  |  Sterile                | 26       |
| 80-06001G    |  Non-Sterile            | 20       |
| 80-06001G-HB |  Non-Sterile, Handi-Bin | N/A      |

- ▶ Oval Fenestration: 8 x 12 in / 20.3 x 30.5 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Fluid Collection Pouch with Filtering Screen and Suction Port
- ▶ Clear Gusseted Anesthesia Panels




# GENERAL DRAPES SURGICAL DRAPES



 **STERIGUARD**

*grī* | alle **set**<sup>®</sup>

DELIVERING PRODUCT SOLUTIONS™




### Bar Drape with Armboard Covers

Size: 106 x 62 in / 269 x 157 cm

| Code         | Summary                                                                                                  | Case Qty |
|--------------|----------------------------------------------------------------------------------------------------------|----------|
| 80-05113G-S  |  Sterile                | 32       |
| 80-05113G    |  Non-Sterile            | 20       |
| 80-05113G-HB |  Non-Sterile, Handi-Bin | N/A      |

- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders
- ▶ Armboard Covers




### Orthopedic Bar Drape

Size: 100 x 60 in / 254 x 152 cm

| Code         | Summary                                                                                                  | Case Qty |
|--------------|----------------------------------------------------------------------------------------------------------|----------|
| 80-05115G-S  |  Sterile                | 32       |
| 80-05115G    |  Non-Sterile            | 28       |
| 80-05115G-HB |  Non-Sterile, Handi-Bin | N/A      |

- ▶ 2 in / 5 cm Adhesive Tape
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders




### Top Drape

Size: 99 x 53 in / 251 x 135 cm

| Code         | Summary                                                                                                    | Case Qty |
|--------------|------------------------------------------------------------------------------------------------------------|----------|
| 80-05152G-S  |  Sterile                | 28       |
| 80-05152G    |  Non-Sterile            | 32       |
| 80-05152G-HB |  Non-Sterile, Handi-Bin | N/A      |

- ▶ 2 in / 5 cm Adhesive Tape
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders



### Side Drape


Size: 38 x 76 in / 97 x 193 cm

| Code         | Summary                                                                                                    | Case Qty |
|--------------|------------------------------------------------------------------------------------------------------------|----------|
| 80-05153G-S  |  Sterile                | 30       |
| 80-05153G    |  Non-Sterile            | 24       |
| 80-05153G-HB |  Non-Sterile, Handi-Bin | N/A      |

- ▶ 2 in / 5 cm Adhesive Tape
- ▶ HydraGuard® Absorbent Reinforcement

## GENERAL DRAPES

### Surgical Drapes




#### Bottom Drape

Size: 75 x 72 in / 191 x 183 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-05154G-S  | <span style="color: red;">S</span> Sterile                 | 18       |
| 80-05154G    | <span style="color: teal;">N</span> Non-Sterile            | 22       |
| 80-05154G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ 2 in / 5 cm Adhesive Tape
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders




#### Split Drape with Absorbent Reinforcement

Size: 53 x 76 in / 135 x 193 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-09135G-S  | <span style="color: red;">S</span> Sterile                 | 26       |
| 80-09135G    | <span style="color: teal;">N</span> Non-Sterile            | 36       |
| 80-09135G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Split: 2.5 x 30 in / 6.4 x 76.2 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Hook and Loop Tube and Line Controls




#### U-Drape with Split

Size: 76 x 120 in / 193 x 305 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-09107G-S  | <span style="color: red;">S</span> Sterile                 | 18       |
| 80-09107G    | <span style="color: teal;">N</span> Non-Sterile            | 20       |
| 80-09107G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | 112      |

- ▶ Split: 4 x 40 in / 10.2 x 101.6 cm
- ▶ 2 in / 5 cm Adhesive Tape
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders




#### U-Drape without Tape

Size: 76 x 120 in / 193 x 305 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-09177G-S  | <span style="color: red;">S</span> Sterile                 | 18       |
| 80-09177G    | <span style="color: teal;">N</span> Non-Sterile            | 20       |
| 80-09177G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Split: 4 x 40 in / 10.2 x 101.6 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders




### U-Drape

Size: 80 x 124 in / 203 x 315 cm

| Code         | Summary                | Case Qty |
|--------------|------------------------|----------|
| 80-09149G-S  | Sterile                | 20       |
| 80-09149G    | Non-Sterile            | 24       |
| 80-09149G-HB | Non-Sterile, Handi-Bin | N/A      |

- ▶ Split: 4 x 44 in / 10.2 x 111.8 cm
- ▶ 2 in / 5 cm Adhesive Tape
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders




### Split Sheet

Size: 77 x 117 in / 196 x 297 cm

| Code        | Summary                | Case Qty |
|-------------|------------------------|----------|
| 80-09173-S  | Sterile                | 24       |
| 80-09173    | Non-Sterile            | 20       |
| 80-09173-HB | Non-Sterile, Handi-Bin | N/A      |

- ▶ Split: 3.5 x 43 in / 9 x 109 cm
- ▶ 1 in / 2.5 cm Adhesive Tape
- ▶ Absorbent Reinforcement
- ▶ Hook and Loop Tube Holder




### Fenestrated Drape

Size: 24 x 30 in / 61 x 76 cm

| Code         | Summary                | Case Qty |
|--------------|------------------------|----------|
| 80-05209G-S  | Sterile                | 80       |
| 80-05209G    | Non-Sterile            | 270      |
| 80-05209G-HB | Non-Sterile, Handi-Bin | N/A      |

- ▶ Circular Fenestration: 4 in / 10 cm
- ▶ Adhesive Tape and Release Liner
- ▶ HydraGuard® Absorbent Reinforcement



### Small Minor Procedure Fenestrated Drape


Size: 30 x 30 in / 76 x 76 cm

| Code         | Summary                | Case Qty |
|--------------|------------------------|----------|
| 80-05173G-S  | Sterile                | 60       |
| 80-05173G    | Non-Sterile            | 70       |
| 80-05173G-HB | Non-Sterile, Handi-Bin | N/A      |

- ▶ Rectangular Fenestration: 1.5 x 2 in / 3.8 x 5.1 cm
- ▶ HydraGuard® Absorbent Reinforcement

## GENERAL DRAPES

### Surgical Drapes




#### Medium Minor Procedure Fenestrated Drape

Size: 48 x 60 in / 122 x 152 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-05178G-S  | <span style="color: red;">S</span> Sterile                 | 40       |
| 80-05178G    | <span style="color: teal;">N</span> Non-Sterile            | 60       |
| 80-05178G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Circular Fenestration: 3 in / 7.6 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders




#### Large Minor Procedure Fenestrated Drape

Size: 76 x 124 in / 193 x 315 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-05103G-S  | <span style="color: red;">S</span> Sterile                 | 20       |
| 80-05103G    | <span style="color: teal;">N</span> Non-Sterile            | 24       |
| 80-05103G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Square Fenestration: 6 x 6 in / 15 x 15 cm
- ▶ HydraGuard® Absorbent Reinforcement




#### Thyroid Drape

Size: 76 x 124 in / 193 x 315 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-05180G-S  | <span style="color: red;">S</span> Sterile                 | 13       |
| 80-05180G    | <span style="color: teal;">N</span> Non-Sterile            | 16       |
| 80-05180G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Diamond Fenestration: 5 x 5 in / 12.7 x 12.7 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Instrument Pad
- ▶ Tube Holders




#### Thyroid Drape with Arm Board Covers

Size: 100 x 72 x 126 in / 254 x 183 x 320 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-05211G-S  | <span style="color: red;">S</span> Sterile                 | 10       |
| 80-05211G    | <span style="color: teal;">N</span> Non-Sterile            | 10       |
| 80-05211G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Diamond Fenestration: 5 x 5 in / 12.7 x 12.7 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Instrument Pad
- ▶ Tube Holders
- ▶ Arm Board Covers




### Breast Drape

Size: 100 x 76 x 126 in / 254 x 193 x 320 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-05214G-S  | <span style="color: red;">S</span> Sterile                 | 10       |
| 80-05214G    | <span style="color: teal;">N</span> Non-Sterile            | 12       |
| 80-05214G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Square Fenestration: 6 x 6 in / 15 x 15 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Hook and Loop Tube and Line Control
- ▶ Armboard Covers




### Chest Drape

Size: 100 x 72 x 126 in / 254 x 183 x 320 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-05122G-S  | <span style="color: red;">S</span> Sterile                 | 12       |
| 80-05122G    | <span style="color: teal;">N</span> Non-Sterile            | 14       |
| 80-05122G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Instrument Pad
- ▶ Tube Holders
- ▶ Armboard Holders




### Heavy-Duty Universal Draping System

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 86-01603-S  | <span style="color: red;">S</span> Sterile                 | 5        |
| 86-01603    | <span style="color: teal;">N</span> Non-Sterile            | 6        |
| 86-01603-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

**Tech70**  
PERFORMANCE

- ▶ TECH70 Performance SMS Fabric
- ▶ 1 Top Drape, 100 x 60 in / 254 x 152 cm
- ▶ 1 Bottom Drape, 76 x 60 in / 193 x 152 cm
- ▶ 2 Side Drapes, 38 x 84 in / 97 x 213 cm
- ▶ Longitudinal Flaps
- ▶ HydraGuard® Absorbent Reinforcement



### Utility Drape with Tape

Size: 15 x 26 in / 38 x 66 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-05126-S  | <span style="color: red;">S</span> Sterile                 | 200      |
| 80-05126    | <span style="color: teal;">N</span> Non-Sterile            | 180      |
| 80-05126-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | 4800     |

- ▶ SMS Fabric




# LAPAROTOMY SURGICAL DRAPE



*grí* | alle<sup>set</sup><sup>®</sup>

DELIVERING PRODUCT SOLUTIONS™




### Pediatric Laparotomy Drape

Size: 74 x 124 in / 188 x 315 cm

| Code         | Summary                | Case Qty |
|--------------|------------------------|----------|
| 80-05119G-S  | Sterile                | 13       |
| 80-05119G    | Non-Sterile            | 19       |
| 80-05119G-HB | Non-Sterile, Handi-Bin | N/A      |

- ▶ Fenestration: 4 x 6 in / 10.1 x 15.2 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders




### Pediatric Laparotomy Drape

Size: 60 x 108 in / 152 x 274 cm

| Code         | Summary                | Case Qty |
|--------------|------------------------|----------|
| 80-05157G-S  | Sterile                | 15       |
| 80-05157G    | Non-Sterile            | 25       |
| 80-05157G-HB | Non-Sterile, Handi-Bin | N/A      |

- ▶ Fenestration: 3 x 8 in / 7.6 x 20.3 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders




### Pediatric Laparotomy Drape with Armboard Covers

Size: 100 x 72 x 124 in / 254 x 183 x 315 cm

| Code         | Summary                | Case Qty |
|--------------|------------------------|----------|
| 80-05212G-S  | Sterile                | 10       |
| 80-05212G    | Non-Sterile            | 10       |
| 80-05212G-HB | Non-Sterile, Handi-Bin | N/A      |

- ▶ Fenestration: 4 x 6 in / 10 x 15 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Armboard Covers



### Laparotomy Drape with Armboard Covers


Size: 104 x 76 x 120 in / 264 x 193 x 305 cm

| Code         | Summary                | Case Qty |
|--------------|------------------------|----------|
| 80-05170G-S  | Sterile                | 13       |
| 80-05170G    | Non-Sterile            | 18       |
| 80-05170G-HB | Non-Sterile, Handi-Bin | N/A      |

- ▶ Fenestration: 4 x 12 in / 10.2 x 30.5 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders
- ▶ Armboard Covers

## LAPAROTOMY

### Surgical Drapes




#### Laparotomy Drape

Size: 76 x 120 in / 193 x 305 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-05171G-S  | <span style="color: red;">S</span> Sterile                 | 16       |
| 80-05171G    | <span style="color: teal;">N</span> Non-Sterile            | 20       |
| 80-05171G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Fenestration: 4 x 12 in / 10.2 x 30.5 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders




#### Laparotomy Drape with Absorbent Troughs

Size: 100 x 72 x 124 in / 254 x 183 x 315 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-05162G-S  | <span style="color: red;">S</span> Sterile                 | 10       |
| 80-05162G    | <span style="color: teal;">N</span> Non-Sterile            | 12       |
| 80-05162G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Fenestration: 4 x 12 in / 10.2 x 30.5 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Absorbent Troughs
- ▶ Tube Holders
- ▶ Armboard Covers




#### Laparotomy Drape with Pouches

Size: 100 x 72 x 124 in / 254 x 183 x 315 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-05164G-S  | <span style="color: red;">S</span> Sterile                 | 10       |
| 80-05164G    | <span style="color: teal;">N</span> Non-Sterile            | 13       |
| 80-05164G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Fenestration: 4 x 12 in / 10.2 x 30.5 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Two Longitudinal Fluid Collection Pouches
- ▶ Tube Holders
- ▶ Armboard Covers




#### Laparotomy Drape with Adjustable Fenestration

Size: 104 x 76 x 120 in / 264 x 193 x 305 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-05123G-S  | <span style="color: red;">S</span> Sterile                 | 10       |
| 80-05123G    | <span style="color: teal;">N</span> Non-Sterile            | 13       |
| 80-05123G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Adjustable Fenestration: 12 x 13.5 in / 30.5 x 34.3 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Hook and Loop Tube and Line Control
- ▶ Armboard Covers




### Laparoscopic Abdominal Drape with Pouches

Size: 102 x 76 x 120 in / 259 x 193 x 305 cm

| Code         | Summary                                                                                                  | Case Qty |
|--------------|----------------------------------------------------------------------------------------------------------|----------|
| 80-05155G-S  |  Sterile                | 10       |
| 80-05155G    |  Non-Sterile            | 9        |
| 80-05155G-HB |  Non-Sterile, Handi-Bin | N/A      |

- ▶ Abdominal Fenestration: 12 x 13 in / 30.5 x 33 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Clear Instrument Pouches
- ▶ Hook and Loop Tube and Line Control
- ▶ Armboard Covers




### Laparoscopic Cholecystectomy Drape

Size: 102 x 76 x 120 in / 259 x 193 x 305 cm

| Code         | Summary                                                                                                  | Case Qty |
|--------------|----------------------------------------------------------------------------------------------------------|----------|
| 80-05125G-S  |  Sterile                | 11       |
| 80-05125G    |  Non-Sterile            | 9        |
| 80-05125G-HB |  Non-Sterile, Handi-Bin | N/A      |

- ▶ Abdominal Fenestration: 12 x 13 in / 30.5 x 33 cm
- ▶ Longitudinal Line Control Flaps
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders
- ▶ Armboard Covers




### Laparoscopic/Cholecystectomy Drape with Instrument Pouch

Size: 102 x 76 x 120 in / 259 x 193 x 305 cm

| Code         | Summary                                                                                                    | Case Qty |
|--------------|------------------------------------------------------------------------------------------------------------|----------|
| 80-05163G-S  |  Sterile                | 10       |
| 80-05163G    |  Non-Sterile            | 11       |
| 80-05163G-HB |  Non-Sterile, Handi-Bin | N/A      |

- ▶ Abdominal Fenestration: 12 x 13 in / 30.5 x 33 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Instrument Pouch
- ▶ Tube Holders
- ▶ Armboard Covers
- ▶ Line Control Flaps



### Transverse Laparotomy Drape


Size: 100 x 72 x 120 in / 254 x 183 x 305 cm

| Code         | Summary                                                                                                    | Case Qty |
|--------------|------------------------------------------------------------------------------------------------------------|----------|
| 80-05172G-S  |  Sterile                | 14       |
| 80-05172G    |  Non-Sterile            | 13       |
| 80-05172G-HB |  Non-Sterile, Handi-Bin | N/A      |

- ▶ Transverse Fenestration: 14 x 4 in / 35.6 x 10.2 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders
- ▶ Armboard Covers

## LAPAROTOMY

### Surgical Drapes




#### Laparotomy T-Drape with Fenestrations

Size: 98 x 41 x 99 in / 249 x 104 x 251 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-05121G-S  | <span style="color: red;">S</span> Sterile                 | 22       |
| 80-05121G    | <span style="color: teal;">N</span> Non-Sterile            | 26       |
| 80-05121G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Abdominal Fenestration: 6 x 8 in / 15 x 20 cm
- ▶ Perineal Fenestration: 3 x 6 in / 7.6 x 15 cm
- ▶ Fenestration Flap
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders




#### Laparoscopy / Pelviscopy Drape

Size: 104 x 104 in / 264 x 264 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-05124G-S  | <span style="color: red;">S</span> Sterile                 | 20       |
| 80-05124G    | <span style="color: teal;">N</span> Non-Sterile            | 24       |
| 80-05124G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Trapezoidal Abdominal Fenestration: 7 x 9 in / 18 x 23 cm
- ▶ Perineal Fenestration: 3 x 6 in / 7.6 x 15 cm
- ▶ Fenestration Flap
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Hook and Loop Tube and Line Controls
- ▶ Arm Board Covers




#### Laparotomy / Pelviscopy Drape with Attached Leggings

Size: 100 x 30 x 123 in / 254 x 76 x 312 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-05182G-S  | <span style="color: red;">S</span> Sterile                 | 9        |
| 80-05182G    | <span style="color: teal;">N</span> Non-Sterile            | 10       |
| 80-05182G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Trapezoidal Abdominal Fenestration: 7 x 9 in / 18 x 23 cm
- ▶ Perineal Fenestration: 3 x 8 in / 7.6 x 20 cm
- ▶ Fenestration Flap
- ▶ Attached Leggings
- ▶ Hook and Loop Tube and Line Control
- ▶ HydraGuard® Absorbent Reinforcement




#### Laparoscopic Pelviscopy Drape with Fluid Collection Pouch

Size: 112 x 30 x 123 in / 284 x 76 x 312 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-05174G-S  | <span style="color: red;">S</span> Sterile                 | 10       |
| 80-05174G    | <span style="color: teal;">N</span> Non-Sterile            | 10       |
| 80-05174G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Trapezoidal Abdominal Fenestration: 7 x 9 in / 18 x 23 cm
- ▶ Perineal Fenestration: 3 x 8 in / 7.6 x 20 cm
- ▶ Attached Leggings
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Hook and Loop Tube and Line Control
- ▶ Armboard Covers



Robotic Drape

Size: 104 x 160 in / 264 x 406 cm

| Code         | Summary                                                                                                  | Case Qty |
|--------------|----------------------------------------------------------------------------------------------------------|----------|
| 80-05181G-S  |  Sterile                | 10       |
| 80-05181G    |  Non-Sterile            | 12       |
| 80-05181G-HB |  Non-Sterile, Handi-Bin | N/A      |

- ▶ Abdominal Fenestration: 12 x 14 in / 30.5 x 35.6 cm
- ▶ Perineal Fenestration: 9 x 4 x 8 in / 22.9 x 10.2 x 20.3 cm
- ▶ Fenestration Flap
- ▶ Clear Instrument Pouches
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Hook and Loop Tube and Line Control




# LITHOTOMY / UROLOGY SURGICAL DRAPE



 **STERIGUARD®**

*grī* | alle **set** 

DELIVERING PRODUCT SOLUTIONS™




### Lithotomy T Sheet Drape

Size: 60 x 31 x 75 in / 152 x 79 x 191 cm

| Code         | Summary                                                                                                  | Case Qty |
|--------------|----------------------------------------------------------------------------------------------------------|----------|
| 80-11102G-S  |  Sterile                | 52       |
| 80-11102G    |  Non-Sterile            | 64       |
| 80-11102G-HB |  Non-Sterile, Handi-Bin | N/A      |

- ▶ Perineal Fenestration: 3 x 6 in / 7.6 x 15.2 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders




### Lithotomy Drape with Attached Leggings

Size: 100 x 30 x 87 in / 254 x 76 x 221 cm

| Code         | Summary                                                                                                    | Case Qty |
|--------------|------------------------------------------------------------------------------------------------------------|----------|
| 80-11101G-S  |  Sterile                | 20       |
| 80-11101G    |  Non-Sterile            | 26       |
| 80-11101G-HB |  Non-Sterile, Handi-Bin | N/A      |

- ▶ Perineal Fenestration: 3 x 6 in / 7.6 x 15 cm
- ▶ Attached Leggings
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders




### Lithotomy Drape with Fluid Collection Pouch

Size: 100 x 30 x 87 in / 254 x 76 x 221 cm

| Code         | Summary                                                                                                    | Case Qty |
|--------------|------------------------------------------------------------------------------------------------------------|----------|
| 80-11111G-S  |  Sterile                | 13       |
| 80-11111G    |  Non-Sterile            | 26       |
| 80-11111G-HB |  Non-Sterile, Handi-Bin | N/A      |

- ▶ Perineal Fenestration: 3 x 6 in / 7.6 x 15.2 cm
- ▶ Fluid Collection Pouch with Graduated Markings
- ▶ Attached Leggings
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders




### GYN Urology Drape with Pouch

Size: 112 x 63 x 131 in / 284 x 160 x 333 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-11114G-S  | <span style="color: red;">S</span> Sterile                 | 14       |
| 80-11114G    | <span style="color: teal;">N</span> Non-Sterile            | 16       |
| 80-11114G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Perineal Fenestration: 3 x 6 in / 7.6 x 15.2 cm
- ▶ Fluid Collection Pouch with Filtering Screen and Suction Port
- ▶ Attached Leggings
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders
- ▶ Armboard Covers




### TUR Drape

Size: 63 x 85 x 102 in / 160 x 217 x 260 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-11106-S  | <span style="color: red;">S</span> Sterile                 | N/A      |
| 80-11106    | <span style="color: teal;">N</span> Non-Sterile            | N/A      |
| 80-11106-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ 6cm Circular Fenestration
- ▶ Fluid Collection Pouch with Filtering Screen and Suction Port
- ▶ Fingercot




### Cystoscopy T-Drape with Fenestrations

Size: 62 x 39 x 102 in / 158 x 100 x 260 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-03103G-S  | <span style="color: red;">S</span> Sterile                 | 32       |
| 80-03103G    | <span style="color: teal;">N</span> Non-Sterile            | 42       |
| 80-03103G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Perineal Fenestration: 3 x 6 in / 7.6 x 15.2 cm
- ▶ Filtering Screen: 9 x 9 in / 22.9 x 22.9 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders




### Cystoscopy T-Drape with Fenestrations

Size: 62 x 39 x 102 in / 158 x 100 x 260 cm

| Code         | Summary                                                                                                  | Case Qty |
|--------------|----------------------------------------------------------------------------------------------------------|----------|
| 80-03102G-S  |  Sterile                | 26       |
| 80-03102G    |  Non-Sterile            | 36       |
| 80-03102G-HB |  Non-Sterile, Handi-Bin | N/A      |

- ▶ Perineal Fenestration: 3 x 6 in / 7.6 x 15.2 cm
- ▶ Filtering Screen: 9 x 9 in / 22.9 x 22.9 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders



### Cystoscopy Drape with Attached Leggings

Size: 100 x 30 x 87 in / 254 x 76 x 220 cm

| Code         | Summary                                                                                                    | Case Qty |
|--------------|------------------------------------------------------------------------------------------------------------|----------|
| 80-03101G-S  |  Sterile                | 18       |
| 80-03101G    |  Non-Sterile            | 20       |
| 80-03101G-HB |  Non-Sterile, Handi-Bin | N/A      |

- ▶ Perineal Fenestration: 3 x 6 in / 7.6 x 15.2 cm
- ▶ Filtering Screen: 9 x 9 in / 22.9 x 22.9 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders




# OB / GYN SURGICAL DRAPE



 **STERIGUARD®**

*grī* | alle  **set**®

DELIVERING PRODUCT SOLUTIONS™




### C-Section Drape with Fenestrated Incise and Parent Viewing Window

Size: 110 x 72 x 127 in / 279 x 183 x 323 cm

| Code        | Summary                                                     | Case Qty |
|-------------|-------------------------------------------------------------|----------|
| 80-07109-S  | <span style="color: red;">S</span> Sterile                  | 10       |
| 80-07109    | <span style="color: green;">N</span> Non-Sterile            | 13       |
| 80-07109-HB | <span style="color: green;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Adhesive Incise Area with Triangular Fenestration: 15 x 13 in / 38.1 x 33 cm
- ▶ Fluid Collection Pouch with Suction Port
- ▶ SMS and Clear Anesthesia Screen for Patient Viewing
- ▶ Malleable Banded Pouch
- ▶ Tube Holders
- ▶ Armboard Covers




### C-Section Drape with Full Incise and Parent Viewing Window

Size: 110 x 72 x 127 in / 279 x 183 x 323 cm

| Code        | Summary                                                     | Case Qty |
|-------------|-------------------------------------------------------------|----------|
| 80-07110-S  | <span style="color: red;">S</span> Sterile                  | 10       |
| 80-07110    | <span style="color: green;">N</span> Non-Sterile            | 13       |
| 80-07110-HB | <span style="color: green;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Full Adhesive Incise Area: 13 x 13 in / 33 x 33 cm
- ▶ Fluid Collection Pouch with Suction Port
- ▶ SMS and Clear Anesthesia Screen for Patient Viewing
- ▶ Malleable Banded Pouch
- ▶ Tube Holders
- ▶ Armboard Covers




### C-Section Drape with Fenestrated Incise and Fluid Collection Pouch

Size: 100 x 72 x 120 in / 254 x 183 x 305 cm

| Code        | Summary                                                     | Case Qty |
|-------------|-------------------------------------------------------------|----------|
| 80-07101-S  | <span style="color: red;">S</span> Sterile                  | 10       |
| 80-07101    | <span style="color: green;">N</span> Non-Sterile            | 13       |
| 80-07101-HB | <span style="color: green;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Adhesive Incise Area with Triangular Fenestration: 15 x 13 in / 38.1 x 33 cm
- ▶ Fluid Collection Pouch with Suction Port
- ▶ Malleable Banded Pouch
- ▶ Tube Holders
- ▶ Armboard Covers




### C-Section Drape with Fenestrated Incise and Fluid Collection Pouch

Size: 100 x 72 x 120 in / 254 x 183 x 305 cm

| Code        | Summary                                                     | Case Qty |
|-------------|-------------------------------------------------------------|----------|
| 80-07102-S  | <span style="color: red;">S</span> Sterile                  | N/A      |
| 80-07102    | <span style="color: green;">N</span> Non-Sterile            | 12       |
| 80-07102-HB | <span style="color: green;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Adhesive Incise Area with Triangular Fenestration: 15 x 13 in / 38.1 x 33 cm
- ▶ Fluid Collection Pouch with Suction Port
- ▶ Foam Banded Pouch
- ▶ Tube Holders
- ▶ Armboard Covers




### C-Section Drape with Full Incise and Fluid Collection Pouch

Size: 100 x 72 x 120 in / 254 x 183 x 305 cm

| Code        | Summary                                                     | Case Qty |
|-------------|-------------------------------------------------------------|----------|
| 80-07103-S  | <span style="color: red;">S</span> Sterile                  | 10       |
| 80-07103    | <span style="color: green;">N</span> Non-Sterile            | 12       |
| 80-07103-HB | <span style="color: green;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Full Adhesive Incise Area: 13 x 13 in / 33 x 33 cm
- ▶ Fluid Collection Pouch with Suction Port
- ▶ Foam Banded Pouch
- ▶ Tube Holders
- ▶ Armboard Covers




### C-Section Drape with Full Incise and Fluid Collection Pouch

Size: 100 x 72 x 120 in / 254 x 183 x 305 cm

| Code        | Summary                                                     | Case Qty |
|-------------|-------------------------------------------------------------|----------|
| 80-07104-S  | <span style="color: red;">S</span> Sterile                  | 11       |
| 80-07104    | <span style="color: green;">N</span> Non-Sterile            | 13       |
| 80-07104-HB | <span style="color: green;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Full Adhesive Incise Area: 13 x 13 in / 33 x 33 cm
- ▶ Fluid Collection Pouch with Suction Port
- ▶ Malleable Banded Pouch
- ▶ Tube Holders
- ▶ Armboard Covers




### C-Section Drape with Teardrop Fenestration

Size: 100 x 72 x 120 in / 254 x 183 x 305 cm

| Code         | Summary                                                     | Case Qty |
|--------------|-------------------------------------------------------------|----------|
| 80-07107G-S  | <span style="color: red;">S</span> Sterile                  | N/A      |
| 80-07107G    | <span style="color: green;">N</span> Non-Sterile            | 16       |
| 80-07107G-HB | <span style="color: green;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Teardrop Fenestration: 8 x 8 in / 20 x 20 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders




### Abdominal Drape

Size: 60 x 38 in / 152 x 97 cm

| Code         | Summary                                                     | Case Qty |
|--------------|-------------------------------------------------------------|----------|
| 80-11110G-S  | <span style="color: red;">S</span> Sterile                  | 24       |
| 80-11110G    | <span style="color: green;">N</span> Non-Sterile            | 72       |
| 80-11110G-HB | <span style="color: green;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ 2 in / 5 cm Adhesive Tape
- ▶ HydraGuard® Absorbent Reinforcement




### Leggings

Size: 31 x 54 in / 79 x 137 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-11105-S  | <span style="color: red;">S</span> Sterile                 | 50       |
| 80-11105    | <span style="color: teal;">N</span> Non-Sterile            | 72       |
| 80-11105-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ SMS Fabric
- ▶ 6 in / 15 cm Cuff




### Clear Leggings

Size: 30 x 35 in / 75 x 90 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-11120-S  | <span style="color: red;">S</span> Sterile                 | N/A      |
| 80-11120    | <span style="color: teal;">N</span> Non-Sterile            | N/A      |
| 80-11120-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Clear PE Film




### Clear Leggings

Size: 33 x 49 in / 84 x 124 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-11115-S  | <span style="color: red;">S</span> Sterile                 | 100      |
| 80-11115    | <span style="color: teal;">N</span> Non-Sterile            | 120      |
| 80-11115-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Clear PE Film




#### Under Buttocks Drape with Fluid Collection Pouch

Size: 40 x 44 in / 102 x 112 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-14106-S  | <span style="color: red;">S</span> Sterile                 | 40       |
| 80-14106    | <span style="color: teal;">N</span> Non-Sterile            | 42       |
| 80-14106-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ SMS Fabric Reinforcement
- ▶ Fluid Collection Pouch with Malleable Band and Suction Port
- ▶ Filtering Screen
- ▶ Graduated Markings




#### Under Buttocks Drape with Fluid Collection Pouch

Size: 40 x 44 in / 102 x 112 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-14112-S  | <span style="color: red;">S</span> Sterile                 | 70       |
| 80-14112    | <span style="color: teal;">N</span> Non-Sterile            | 70       |
| 80-14112-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ SMS Fabric Reinforcement
- ▶ Fluid Collection Pouch




#### Under Buttocks Drape with Fluid Collection Pouch

Size: 40 x 44 in / 102 x 112 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-14113-S  | <span style="color: red;">S</span> Sterile                 | 70       |
| 80-14113    | <span style="color: teal;">N</span> Non-Sterile            | 76       |
| 80-14113-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ SMS Fabric Reinforcement
- ▶ Fluid Collection Pouch
- ▶ Suction Port
- ▶ Graduated Markings
- ▶ Filtering Screen



#### Under Buttocks Drape with Fluid Collection Pouch

Size: 40 x 44 in / 102 x 112 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-14111-S  | <span style="color: red;">S</span> Sterile                 | 80       |
| 80-14111    | <span style="color: teal;">N</span> Non-Sterile            | 80       |
| 80-14111-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ SMS Fabric Reinforcement
- ▶ Fluid Collection Pouch




### Under Buttocks Drape

Size: 40 x 44 in / 102 x 112 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-14108-S  | <span style="color: red;">S</span> Sterile                 | 100      |
| 80-14108    | <span style="color: teal;">N</span> Non-Sterile            | 100      |
| 80-14108-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

► PE Film



### Under Buttocks Drape with Reinforcement

Size: 38 x 54 in / 97 x 137 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-14104-S  | <span style="color: red;">S</span> Sterile                 | 30       |
| 80-14104    | <span style="color: teal;">N</span> Non-Sterile            | 100      |
| 80-14104-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |


► Film Fabric Reinforcement



# OPHTHALMIC SURGICAL DRAPES



DELIVERING PRODUCT SOLUTIONS™




### Ophthalmic Drape with Aperture, Incise and Fluid Collection Pouches

Size: 40 x 48 in / 102 x 122 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-08126-S  | <span style="color: red;">S</span> Sterile                 | 80       |
| 80-08126    | <span style="color: teal;">N</span> Non-Sterile            | 200      |
| 80-08126-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Rectangular Adhesive Incise Area: 5 x 4 in / 12.7 x 10.2 cm
- ▶ Oval Fenestration: 3 x 2 in / 7.6 x 5.1 cm
- ▶ Fluid Collection Pouch
- ▶ Clear Micro Embossed Film




### Ophthalmic Drape with Embossed Film, Full Incise and Fluid Pouch

Size: 40 x 48 in / 102 x 122 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-08109-S  | <span style="color: red;">S</span> Sterile                 | 80       |
| 80-08109    | <span style="color: teal;">N</span> Non-Sterile            | 200      |
| 80-08109-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Rectangular Adhesive Incise Area: 5 x 4 in / 12.7 x 10.2 cm
- ▶ Fluid Collection Pouch
- ▶ Clear Micro Embossed Film




### Ophthalmic Drape with Full Incise and Fluid Pouch

Size: 40 x 40 in / 102 x 102 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-04110-S  | <span style="color: red;">S</span> Sterile                 | 50       |
| 80-04110    | <span style="color: teal;">N</span> Non-Sterile            | 100      |
| 80-04110-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |


- ▶ Adhesive Incise Area: 5 x 3.5 in / 12.7 x 8.9 cm
- ▶ Fluid Collection Pouch
- ▶ Eye Wick



### Ophthalmic Drape

Size: 59 x 59 in / 150 x 150 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-04407-S  | <span style="color: red;">S</span> Sterile                 | N/A      |
| 80-04407    | <span style="color: teal;">N</span> Non-Sterile            | N/A      |
| 80-04407-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |




### Ophthalmic Drape with Aperture and Fluid Collection Pouches

Size: 40 x 48 in / 102 x 122 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-08127-S  | <span style="color: red;">S</span> Sterile                 | 54       |
| 80-08127    | <span style="color: teal;">N</span> Non-Sterile            | 66       |
| 80-08127-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Rectangular Adhesive Incise Area: 5 x 4 in / 12.7 x 10.2 cm
- ▶ Oval Fenestration: 3 x 2 in / 7.6 x 5.1 cm
- ▶ Fluid Collection Pouch




### Ophthalmic Drape with Aperture, Eye Wick and Fluid Collection Pouches

Size: 47 x 47 in / 119 x 119 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-08124-S  | <span style="color: red;">S</span> Sterile                 | 54       |
| 80-08124    | <span style="color: teal;">N</span> Non-Sterile            | 66       |
| 80-08124-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Rectangular Adhesive Incise Area: 5 x 4 in / 12.7 x 10.2 cm
- ▶ Fluid Collection Pouch
- ▶ Malleable Breather Bar




### Ophthalmic Drape with Aperture, Eye Wick and Fluid Collection Pouches

Size: 35 x 35 in / 89 x 89 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-08123-S  | <span style="color: red;">S</span> Sterile                 | 60       |
| 80-08123    | <span style="color: teal;">N</span> Non-Sterile            | 180      |
| 80-08123-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Rectangular Adhesive Incise Area: 5 x 4 in / 12.7 x 10.2 cm
- ▶ Oval Fenestration: 3 x 2 in / 7.6 x 5.1 cm
- ▶ Fluid Collection Pouch
- ▶ Malleable Breather Bar




### Ophthalmic Drape with Oval Fenestration

Size: 56 x 52 in / 142 x 132 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-04103-S  | <span style="color: red;">S</span> Sterile                 | 44       |
| 80-04103    | <span style="color: teal;">N</span> Non-Sterile            | 60       |
| 80-04103-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Adhesive Incise Area: 5 x 3.5 in / 12.7 x 8.9 cm
- ▶ Oval Fenestration: 3 x 2 in / 7.6 x 5.1 cm
- ▶ Fluid Collection Pouch
- ▶ Eye Wick




### Ophthalmic Drape with Full Incise

Size: 56 x 52 in / 142 x 132 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-04106-S  | <span style="color: red;">S</span> Sterile                 | 44       |
| 80-04106    | <span style="color: teal;">N</span> Non-Sterile            | 60       |
| 80-04106-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Adhesive Incise Area: 5 x 3.5 in / 12.7 x 8.9 cm
- ▶ Fluid Collection Pouch
- ▶ Eye Wick
- ▶ Malleable Breather Bar




### Ophthalmic Drape with Aperture and Fluid Collection Pouches

Size: 48 x 55 in / 122 x 140 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-08117-S  | <span style="color: red;">S</span> Sterile                 | 54       |
| 80-08117    | <span style="color: teal;">N</span> Non-Sterile            | 66       |
| 80-08117-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Adhesive Fenestrated Incise Area: 5 x 3.5 in / 12.7 x 8.9 cm
- ▶ Oval Fenestration: 3 x 2 in / 7.6 x 5.1 cm
- ▶ Fluid Collection Pouch
- ▶ Malleable Breather Bar




### Ophthalmic Drape with Full Incise and Fluid Collection Pouches

Size: 48 x 55 in / 122 x 140 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-08118-S  | <span style="color: red;">S</span> Sterile                 | 54       |
| 80-08118    | <span style="color: teal;">N</span> Non-Sterile            | 66       |
| 80-08118-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Adhesive Incise Area: 5 x 3.5 in / 12.7 x 8.9 cm
- ▶ Fluid Collection Pouch
- ▶ Malleable Breather Bar




### Ophthalmic Drape with Aperture

Size: 10.6 x 16 in / 27 x 41 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-08125-S  | <span style="color: red;">S</span> Sterile                 | N/A      |
| 80-08125    | <span style="color: teal;">N</span> Non-Sterile            | 250      |
| 80-08125-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Oval Fenestration: 3 x 2 in / 7.6 x 5.1 cm
- ▶ Clear Micro Embossed Film




### Ophthalmic Drape with Aperture and Fluid Collection Pouches

Size: 55 x 71 in / 140 x 180 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-08119-S  | <span style="color: red;">S</span> Sterile                 | 17       |
| 80-08119    | <span style="color: teal;">N</span> Non-Sterile            | 50       |
| 80-08119-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Adhesive Fenestrated Incise Area: 5 x 3.5 in / 12.7 x 8.9 cm
- ▶ Oval Fenestration: 3 x 2 in / 7.6 x 5.1 cm
- ▶ Saddlebag Fluid Collection Pouch
- ▶ Malleable Breather Bar




### Ophthalmic Drape with Full Incise and Fluid Collection Pouch

Size: 55 x 71 in / 140 x 180 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-08120-S  | <span style="color: red;">S</span> Sterile                 | 17       |
| 80-08120    | <span style="color: teal;">N</span> Non-Sterile            | 50       |
| 80-08120-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |


- ▶ Adhesive Incise Area: 5 x 3.5 in / 12.7 x 8.9 cm
- ▶ Saddlebag Fluid Collection Pouch
- ▶ Malleable Breather Bar



### Ophthalmic Drape with Full Incise and Fluid Collection Pouches

Size: 63 x 71 in / 160 x 180 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-04115-S  | <span style="color: red;">S</span> Sterile                 | N/A      |
| 80-04115    | <span style="color: teal;">N</span> Non-Sterile            | N/A      |
| 80-04115-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |




### Ophthalmic Drape with Oval Incise and Fluid Pouch

Size: 77 x 107 in / 196 x 272 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-08110-S  | <span style="color: red;">S</span> Sterile                 | 12       |
| 80-08110    | <span style="color: teal;">N</span> Non-Sterile            | 40       |
| 80-08110-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Adhesive Incise Area: 5 x 3.5 in / 12.7 x 8.9 cm
- ▶ Oval Fenestration: 3 x 2 in / 7.6 x 5.1 cm
- ▶ Saddlebag Fluid Collection Pouch
- ▶ Malleable Breather Bar
- ▶ Full Body Coverage




### Ophthalmic Drape with Aperture and Fluid Collection Pouch

Size: 71 x 100 in / 180 x 254 cm

| Code        | Summary                                                                                                  | Case Qty |
|-------------|----------------------------------------------------------------------------------------------------------|----------|
| 80-08121-S  |  Sterile                | 14       |
| 80-08121    |  Non-Sterile            | 32       |
| 80-08121-HB |  Non-Sterile, Handi-Bin | N/A      |

- ▶ Adhesive Fenestrated Incise Area: 5 x 3.5 in / 12.7 x 8.9 cm
- ▶ Oval Fenestration: 3 x 2 in / 7.6 x 5.1 cm
- ▶ Saddlebag Fluid Collection Pouch
- ▶ Malleable Breather Bar
- ▶ Full Body Coverage
- ▶ Tube and Line Control Holders




### Ophthalmic with Full Incise and Fluid Collection Pouch

Size: 71 x 100 in / 180 x 254 cm

| Code        | Summary                                                                                                  | Case Qty |
|-------------|----------------------------------------------------------------------------------------------------------|----------|
| 80-08122-S  |  Sterile                | 14       |
| 80-08122    |  Non-Sterile            | 32       |
| 80-08122-HB |  Non-Sterile, Handi-Bin | N/A      |

- ▶ Adhesive Incise Area: 5 x 3.5 in / 12.7 x 8.9 cm
- ▶ Saddlebag Fluid Collection Pouch
- ▶ Malleable Breather Bar
- ▶ Malleable Breather Bar
- ▶ Full Body Coverage
- ▶ Tube and Line Control Holders



### Ophthalmic Drape with Fluid Collection Pouch

Size: 72 x 115 in / 183 x 292 cm

| Code         | Summary                                                                                                    | Case Qty |
|--------------|------------------------------------------------------------------------------------------------------------|----------|
| 80-08102G-S  |  Sterile                | 14       |
| 80-08102G    |  Non-Sterile            | 12       |
| 80-08102G-HB |  Non-Sterile, Handi-Bin | N/A      |

- ▶ Adhesive Incise Area: 4.5 x 4.5 in / 11.4 x 11.4 cm
- ▶ Optional Incise Drape: 5 x 6 in / 12.7 x 15.2 cm
- ▶ 270 Degree Fluid Collection Pouch
- ▶ Malleable Breather Bar
- ▶ Full Body Coverage
- ▶ Line Control and Tube Holders




# ORTHOPEDIC SURGICAL DRAPE



*grī* | alle **set**<sup>®</sup>

DELIVERING PRODUCT SOLUTIONS™




### Bar Drape with Armboard Covers

Size: 106 x 62 in / 269 x 157 cm

| Code         | Summary                | Case Qty |
|--------------|------------------------|----------|
| 80-05113G-S  | Sterile                | 32       |
| 80-05113G    | Non-Sterile            | 20       |
| 80-05113G-HB | Non-Sterile, Handi-Bin | N/A      |

- ▶ 2 in / 5 cm Adhesive Tape
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders
- ▶ Armboard Covers




### Orthopedic Bar Drape

Size: 100 x 60 in / 254 x 152 cm

| Code         | Summary                | Case Qty |
|--------------|------------------------|----------|
| 80-05115G-S  | Sterile                | 32       |
| 80-05115G    | Non-Sterile            | 28       |
| 80-05115G-HB | Non-Sterile, Handi-Bin | N/A      |

- ▶ 2 in / 5 cm Adhesive Tape
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders




### TECH70 Heavy-Duty Bar Drape

Size: 112 x 72 in / 284 x 183 cm

| Code         | Summary                | Case Qty |
|--------------|------------------------|----------|
| 80-05215G-S  | Sterile                | 22       |
| 80-05215G    | Non-Sterile            | 22       |
| 80-05215G-HB | Non-Sterile, Handi-Bin | N/A      |

- ▶ TECH70 Performance SMS Fabric




### Impervious Split Drape

Size: 60 x 70 in / 152 x 178 cm

| Code        | Summary                | Case Qty |
|-------------|------------------------|----------|
| 80-09106-S  | Sterile                | 92       |
| 80-09106    | Non-Sterile            | 120      |
| 80-09106-HB | Non-Sterile, Handi-Bin | N/A      |

- ▶ Split: 4 x 21 in / 10.2 x 53.3 cm
- ▶ 2 in / 5 cm Adhesive Tape




### Large Impervious Split Drape

Size: 60 x 81 in / 152 x 206 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-09147-S  | <span style="color: red;">S</span> Sterile                 | 84       |
| 80-09147    | <span style="color: teal;">N</span> Non-Sterile            | 100      |
| 80-09147-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Split: 4 x 21 in / 10 x 53 cm
- ▶ Blue PE Film




### U-Drape with Pouch

Size: 35 x 30 in / 89 x 76 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-09169-S  | <span style="color: red;">S</span> Sterile                 | 18       |
| 80-09169    | <span style="color: teal;">N</span> Non-Sterile            | 50       |
| 80-09169-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Adhesive Split: 6 x 15 in / 16 x 40 cm
- ▶ U-Shaped Fluid Collection Pouch
- ▶ 2 Suction Ports




### Orthopedic Large U Drape

Size: 88 x 125 in / 224 x 318 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-09166G-S  | <span style="color: red;">S</span> Sterile                 | 17       |
| 80-09166G    | <span style="color: teal;">N</span> Non-Sterile            | 12       |
| 80-09166G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Split: 6 x 54 in / 15 x 137 cm
- ▶ HydraGuard® Absorbent Reinforcement




### TECH70 Body Split Drape

Size: 88 x 116 in / 224 x 295 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-09105G-S  | <span style="color: red;">S</span> Sterile                 | 12       |
| 80-09105G    | <span style="color: teal;">N</span> Non-Sterile            | 20       |
| 80-09105G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ TECH70 Performance SMS Fabric
- ▶ Split: 4 X 27 in / 10 x 69 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders

**Tech70**  
PERFORMANCE




### TECH70 Orthopedic U-Drape

Size: 88 x 108 in / 224 x 274 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-09155G-S  | <span style="color: red;">S</span> Sterile                 | 10       |
| 80-09155G    | <span style="color: teal;">N</span> Non-Sterile            | 13       |
| 80-09155G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

**Tech70**  
PERFORMANCE

- ▶ TECH70 Performance SMS Fabric
- ▶ Split: 4 x 40 in / 10.2 x 101.6 cm
- ▶ HydraGuard® Absorbent Reinforcement




### Universal Extremity Drape

Size: 90 x 124 in / 229 x 315 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-09101G-S  | <span style="color: red;">S</span> Sterile                 | 10       |
| 80-09101G    | <span style="color: teal;">N</span> Non-Sterile            | 13       |
| 80-09101G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Conformable Circular Fenestration: 2.5 in / 6.4 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders




### Lower Extremity Drape

Size: 88 x 131 in / 223 x 332 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-09157G-S  | <span style="color: red;">S</span> Sterile                 | 10       |
| 80-09157G    | <span style="color: teal;">N</span> Non-Sterile            | 13       |
| 80-09157G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Conformable Fenestration: 2.5 in / 6.4 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders




### Lower Extremity Drape

Size: 108 x 80 x 124 in / 274 x 203 x 315 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-09158G-S  | <span style="color: red;">S</span> Sterile                 | 10       |
| 80-09158G    | <span style="color: teal;">N</span> Non-Sterile            | 13       |
| 80-09158G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Conformable Fenestration: 2.5 in / 6.4 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders
- ▶ Armboard Covers




### Lower Extremity Drape (Extra Long)

Size: 114 x 80 x 144 in / 290 x 203 x 366 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-09151G-S  | <span style="color: red;">S</span> Sterile                 | 10       |
| 80-09151G    | <span style="color: teal;">N</span> Non-Sterile            | 13       |
| 80-09151G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Conformable Fenestration: 2.5 in / 6.4 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders
- ▶ Armboard Covers




### TECH70 Lower Extremity Drape

Size: 114 x 88 x 126 in / 290 x 224 x 320 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-09152G-S  | <span style="color: red;">S</span> Sterile                 | 10       |
| 80-09152G    | <span style="color: teal;">N</span> Non-Sterile            | 12       |
| 80-09152G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ TECH70 Performance SMS Fabric
- ▶ Conformable Fenestration: 2.5 in / 6.4 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders
- ▶ Armboard Covers




### Bilateral Limb Drape

Size: 116 x 88 x 144 in / 295 x 224 x 366 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-09138G-S  | <span style="color: red;">S</span> Sterile                 | 12       |
| 80-09138G    | <span style="color: teal;">N</span> Non-Sterile            | 12       |
| 80-09138G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Two Conformable Circular Fenestration: 2.5 in / 6.4 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders
- ▶ Armboard Covers




### Knee Arthroscopy Drape

Size: 90 x 124 in / 229 x 315 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-09102-S  | <span style="color: red;">S</span> Sterile                 | 11       |
| 80-09102    | <span style="color: teal;">N</span> Non-Sterile            | 12       |
| 80-09102-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Conformable Circular Fenestration: 2.5 in / 6.4 cm
- ▶ Fluid Collection Pouch with Suction Port
- ▶ Tube Holders




### Knee Arthroscopy Drape with Armboard Covers

Size: 114 x 88 x 124 in / 290 x 224 x 315 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-09156-S  | <span style="color: red;">S</span> Sterile                 | 8        |
| 80-09156    | <span style="color: teal;">N</span> Non-Sterile            | 10       |
| 80-09156-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Fluid Collection Pouch
- ▶ Tube Holders
- ▶ Armboard Covers




### TECH70 Knee Arthroscopy Drape with Absorbent Reinforcement

Size: 114 x 88 x 116 in / 290 x 224 x 295 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-09153G-S  | <span style="color: red;">S</span> Sterile                 | 8        |
| 80-09153G    | <span style="color: teal;">N</span> Non-Sterile            | 10       |
| 80-09153G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ TECH70 Performance SMS Fabric
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Fluid Collection Pouch with Suction Port
- ▶ Tube Holders
- ▶ Armboard Covers




### Hip Drape with Clear Leg Pouches

Size: 116 x 88 x 136 in / 295 x 224 x 345 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-09154G-S  | <span style="color: red;">S</span> Sterile                 | 7        |
| 80-09154G    | <span style="color: teal;">N</span> Non-Sterile            | 9        |
| 80-09154G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Conformable Teardrop Fenestration: 8.25 x 6.75 in / 21 x 17.1 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Clear Leg Pouches
- ▶ Hook and Loop Tube Holders




### TECH70 Orthopedic Hip Drape with Leg Pockets

Size: 114 x 88 x 126 in / 290 x 224 x 320 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-09104G-S  | <span style="color: red;">S</span> Sterile                 | 6        |
| 80-09104G    | <span style="color: teal;">N</span> Non-Sterile            | 8        |
| 80-09104G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ TECH70 Performance SMS Fabric
- ▶ Conformable Teardrop Fenestration: 8.25 x 6.75 in / 21 x 17.1 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders
- ▶ Armboard Covers
- ▶ Side Leg Pockets






### TECH70 Orthopedic Hip Drape with Teardrop Fenestration

Size: 114 x 88 x 126 in / 290 x 224 x 320 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-09119G-S  | <span style="color: red;">S</span> Sterile                 | 7        |
| 80-09119G    | <span style="color: teal;">N</span> Non-Sterile            | 7        |
| 80-09119G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |



- ▶ TECH70 Performance SMS Fabric
- ▶ Conformable Teardrop Fenestration: 8.25 x 6.75 in / 21 x 17.1 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders
- ▶ Armboard Covers




### Shoulder Arthroscopy Drape

Size: 100 x 60 in / 254 x 152 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-09109-S  | <span style="color: red;">S</span> Sterile                 | 12       |
| 80-09109    | <span style="color: teal;">N</span> Non-Sterile            | 10       |
| 80-09109-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Fluid Collection Pouch
- ▶ Heavy-Duty Fabric




### TECH70 Heavy-Duty Shoulder Arthroscopy Drape

Size: 100 x 60 in / 254 x 152 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-09172-S  | <span style="color: red;">S</span> Sterile                 | 12       |
| 80-09172    | <span style="color: teal;">N</span> Non-Sterile            | 10       |
| 80-09172-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |



- ▶ TECH70 Performance SMS Fabric
- ▶ Fluid Collection Pouch




### Shoulder Arthroscopy with Fluid Collection Pouch

Size: 102 x 160 in / 259 x 406 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-09123G-S  | <span style="color: red;">S</span> Sterile                 | 10       |
| 80-09123G    | <span style="color: teal;">N</span> Non-Sterile            | 10       |
| 80-09123G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Conformable Teardrop Fenestration: 2.5 in / 6.4 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Attached Fluid Collection Pouch




### Beach Chair Shoulder Drape

Size: 100 x 172 in / 254 x 437 cm

| Code         | Summary                                                    | Case Qty |
|--------------|------------------------------------------------------------|----------|
| 80-09111G-S  | <span style="color: red;">S</span> Sterile                 | 7        |
| 80-09111G    | <span style="color: teal;">N</span> Non-Sterile            | 16       |
| 80-09111G-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ HydraGuard® Absorbent Reinforcement
- ▶ 2 Pouch Attachment Areas
- ▶ Hook and Loop Tube and Line Controls




### Attachable Fluid Collection Pouch

Size: 35 x 43 in / 89 x 109 cm

| Code        | Summary                                                    | Case Qty |
|-------------|------------------------------------------------------------|----------|
| 80-14122-S  | <span style="color: red;">S</span> Sterile                 | N/A      |
| 80-14122    | <span style="color: teal;">N</span> Non-Sterile            | 90       |
| 80-14122-HB | <span style="color: teal;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ X-Large Pleated Fluid Collection Pouch
- ▶ Suction Port




### Shoulder Arthroscopy Drape

Size: 100 x 62 in / 254 x 157 cm

| Code        | Summary                                                     | Case Qty |
|-------------|-------------------------------------------------------------|----------|
| 80-09159-S  | <span style="color: red;">S</span> Sterile                  | N/A      |
| 80-09159    | <span style="color: green;">N</span> Non-Sterile            | 30       |
| 80-09159-HB | <span style="color: green;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Fluid Collection Pouch with Suction Ports
- ▶ Tube Holders




### Shoulder Split Drape with Fluid Collection Pouch

Size: 98 x 54 in / 249 x 137 cm

| Code         | Summary                                                     | Case Qty |
|--------------|-------------------------------------------------------------|----------|
| 80-09167G-S  | <span style="color: red;">S</span> Sterile                  | 20       |
| 80-09167G    | <span style="color: green;">N</span> Non-Sterile            | 22       |
| 80-09167G-HB | <span style="color: green;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Pear-shaped Incise Aperture Fenestration: 11 x 14 in / 28 x 36 cm
- ▶ Fluid Collection Pouch with Suction Ports
- ▶ Tube Holders




### 2-Piece Shoulder Arthroscopy Draping System

| Code        | Summary                                                     | Case Qty |
|-------------|-------------------------------------------------------------|----------|
| 86-01527-S  | <span style="color: red;">S</span> Sterile                  | 12       |
| 86-01527    | <span style="color: green;">N</span> Non-Sterile            | N/A      |
| 86-01527-HB | <span style="color: green;">N</span> Non-Sterile, Handi-Bin | N/A      |

**Tech70**  
PERFORMANCE

- 1 Body Split Drape (80-09105G)
  - ▶ 88 x 116 in / 224 x 295 cm
  - ▶ TECH70 Performance SMS Fabric
  - ▶ Split: 4 X 27 in / 10 x 69 cm
  - ▶ HydraGuard® Absorbent Reinforcement
  - ▶ Tube Holders

- 1 Shoulder Arthroscopy Drape
  - ▶ 62 x 100 in / 157 x 254 cm
  - ▶ Fluid Collection Pouch with Suction Ports
  - ▶ Tube Holders




### Spine Drape with Armboard Covers

Size: 112 x 99 x 138 in / 284 x 251 x 351 cm

| Code         | Summary                                                     | Case Qty |
|--------------|-------------------------------------------------------------|----------|
| 80-09168G-S  | <span style="color: red;">S</span> Sterile                  | 10       |
| 80-09168G    | <span style="color: green;">N</span> Non-Sterile            | 10       |
| 80-09168G-HB | <span style="color: green;">N</span> Non-Sterile, Handi-Bin | N/A      |

- ▶ Fenestration: 4 x 12 in / 10 x 30 cm
- ▶ Clear Instrument Pouches
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Hook and Loop Tube and Line Controls
- ▶ Armboard Covers




### TECH70 Hand Drape

Size: 76 x 150 in / 193 x 381 cm

| Code         | Summary                                                                                                  | Case Qty |
|--------------|----------------------------------------------------------------------------------------------------------|----------|
| 80-09108G-S  |  Sterile                | 14       |
| 80-09108G    |  Non-Sterile            | 13       |
| 80-09108G-HB |  Non-Sterile, Handi-Bin | N/A      |



- ▶ TECH70 Performance SMS Fabric
- ▶ Conformable Fenestration: 1.25 in / 3.2 cm
- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Tube Holders
- ▶ Heavy-Duty Fabric
- ▶ Gusseted Hand Table Cover: 62 x 58 in / 157.5 x 147.3 cm



### Vertical Isolation Drape

Size: 130 x 91 in / 330 x 230 cm

| Code        | Summary                                                                                                  | Case Qty |
|-------------|----------------------------------------------------------------------------------------------------------|----------|
| 80-09133-S  |  Sterile                | 15       |
| 80-09133    |  Non-Sterile            | N/A      |
| 80-09133-HB |  Non-Sterile, Handi-Bin | N/A      |


- ▶ Clear PE Film
- ▶ Wide Adhesive Incise Area
- ▶ Hook and loop tube holders
- ▶ Central Release Adhesive Strips for Horizontal Rod Attachment
- ▶ Fluid Collection Pouch with Suction Port



# EQUIPMENT COVERS SURGICAL DRAPE



DELIVERING PRODUCT SOLUTIONS™




### Absorbent Heavy Duty Back Table Cover

Size: 60 x 90 in / 152 x 229 cm

| Code        | Summary                                                                                                               | Case Qty |
|-------------|-----------------------------------------------------------------------------------------------------------------------|----------|
| 80-13322-S  |  Sterile, Fan-Folded                 | 26       |
| 80-13322    |  Non-Sterile, Fan-Folded             | 14       |
| 80-13322-HB |  Non-Sterile, Soft-Folded, Handi-Bin | 140      |

- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Heavy-Duty Film




### Absorbent Heavy Duty Back Table Cover

Size: 80 x 110 in / 203 x 279 cm

| Code        | Summary                                                                                                               | Case Qty |
|-------------|-----------------------------------------------------------------------------------------------------------------------|----------|
| 80-13323-S  |  Sterile, Fan-Folded                 | 20       |
| 80-13323    |  Non-Sterile, Fan-Folded             | 18       |
| 80-13323-HB |  Non-Sterile, Soft-Folded, Handi-Bin | 240      |

- ▶ HydraGuard® Absorbent Reinforcement
- ▶ Heavy-Duty Film




### Fully Reinforced Heavy Duty Back Table Cover

Size: 70 x 110 in / 178 x 279 cm

| Code        | Summary                                                                                                                 | Case Qty |
|-------------|-------------------------------------------------------------------------------------------------------------------------|----------|
| 80-13317-S  |  Sterile, Fan-Folded                 | 16       |
| 80-13317    |  Non-Sterile, Fan-Folded             | 20       |
| 80-13317-HB |  Non-Sterile, Soft-Folded, Handi-Bin | 250      |

- ▶ SMS Reinforced Fabric
- ▶ Advanced eCoat Technology



### Fully Reinforced Heavy Duty Back Table Cover


Size: 60 x 90 in / 152 x 229 cm

| Code        | Summary                                                                                                                 | Case Qty |
|-------------|-------------------------------------------------------------------------------------------------------------------------|----------|
| 80-13336-S  |  Sterile, Fan-Folded                 | 32       |
| 80-13336    |  Non-Sterile, Fan-Folded             | 40       |
| 80-13336-HB |  Non-Sterile, Soft-Folded, Handi-Bin | 400      |

- ▶ SMS Reinforced Fabric
- ▶ Advanced eCoat Technology

## EQUIPMENT COVERS

### Surgical Drapes




#### Fully Reinforced Back Table Cover

Size: 60 x 90 in / 152 x 229 cm

| Code        | Summary                                                                                                               | Case Qty |
|-------------|-----------------------------------------------------------------------------------------------------------------------|----------|
| 80-13303-S  |  Sterile, Fan-Folded                 | N/A      |
| 80-13303    |  Non-Sterile, Fan-Folded             | 20       |
| 80-13303-HB |  Non-Sterile, Soft-Folded, Handi-Bin | 252      |

- ▶ SMS Reinforced Fabric
- ▶ Advanced eCoat Technology




#### Fully Reinforced Back Table Cover

Size: 75 x 79 in / 191 x 201 cm

| Code        | Summary                                                                                                               | Case Qty |
|-------------|-----------------------------------------------------------------------------------------------------------------------|----------|
| 80-13304-S  |  Sterile, Fan-Folded                 | 20       |
| 80-13304    |  Non-Sterile, Fan-Folded             | 32       |
| 80-13304-HB |  Non-Sterile, Soft-Folded, Handi-Bin | N/A      |

- ▶ SMS Reinforced Fabric
- ▶ Advanced eCoat Technology




#### Zone Reinforced Heavy Duty Back Table Cover

Size: 44 x 75 in / 112 x 191 cm

| Code        | Summary                                                                                                                 | Case Qty |
|-------------|-------------------------------------------------------------------------------------------------------------------------|----------|
| 80-13101-S  |  Sterile, Fan-Folded                 | 40       |
| 80-13101    |  Non-Sterile, Fan-Folded             | 52       |
| 80-13101-HB |  Non-Sterile, Soft-Folded, Handi-Bin | 320      |

- ▶ SMS Reinforced Fabric
- ▶ Advanced eCoat Technology




#### Zone Reinforced Heavy Duty Back Table Cover

Size: 50 x 90 in / 127 x 229 cm

| Code        | Summary                                                                                                                 | Case Qty |
|-------------|-------------------------------------------------------------------------------------------------------------------------|----------|
| 80-13103-S  |  Sterile, Fan-Folded                 | 26       |
| 80-13103    |  Non-Sterile, Fan-Folded             | 40       |
| 80-13103-HB |  Non-Sterile, Soft-Folded, Handi-Bin | 320      |

- ▶ SMS Reinforced Fabric
- ▶ Advanced eCoat Technology




#### Zone Reinforced Heavy Duty Back Table Cover

Size: 60 x 90 in / 152 x 229 cm

| Code        | Summary                                                                                                               | Case Qty |
|-------------|-----------------------------------------------------------------------------------------------------------------------|----------|
| 80-13109-S  |  Sterile, Fan-Folded                 | 26       |
| 80-13109    |  Non-Sterile, Soft-Folded            | 36       |
| 80-13109-HB |  Non-Sterile, Soft-Folded, Handi-Bin | N/A      |

- ▶ SMS Reinforced Fabric
- ▶ Advanced eCoat Technology




#### Zone Reinforced Heavy Duty Back Table Cover

Size: 60 x 75 in / 152 x 190 cm

| Code        | Summary                                                                                                               | Case Qty |
|-------------|-----------------------------------------------------------------------------------------------------------------------|----------|
| 80-13108-S  |  Sterile, Fan-Folded                 | 28       |
| 80-13108    |  Non-Sterile, Fan-Folded             | 42       |
| 80-13108-HB |  Non-Sterile, Soft-Folded, Handi-Bin | 300      |

- ▶ SMS Reinforced Fabric
- ▶ Advanced eCoat Technology




#### Zone Reinforced Back Table Cover

Size: 44 x 90 in / 112 x 229 cm

| Code        | Summary                                                                                                                 | Case Qty |
|-------------|-------------------------------------------------------------------------------------------------------------------------|----------|
| 80-13310-S  |  Sterile, Fan-Folded                 | 20       |
| 80-13310    |  Non-Sterile, Fan-Folded             | 60       |
| 80-13310-HB |  Non-Sterile, Soft-Folded, Handi-Bin | 480      |

- ▶ SMS Reinforced Fabric
- ▶ Advanced eCoat Technology



#### Zone Reinforced Back Table Cover


Size: 44 x 44 in / 112 x 112 cm

| Code        | Summary                                                                                                                 | Case Qty |
|-------------|-------------------------------------------------------------------------------------------------------------------------|----------|
| 80-13311-S  |  Sterile, Fan-Folded                 | 80       |
| 80-13311    |  Non-Sterile, Fan-Folded             | 100      |
| 80-13311-HB |  Non-Sterile, Soft-Folded, Handi-Bin | 1040     |

- ▶ SMS Reinforced Fabric
- ▶ Advanced eCoat Technology

## EQUIPMENT COVERS

### Surgical Drapes




#### Zone Reinforced Back Table Cover

Size: 70 x 110 in / 178 x 279 cm

| Code        | Summary                                                                                                               | Case Qty |
|-------------|-----------------------------------------------------------------------------------------------------------------------|----------|
| 80-13312-S  |  Sterile, Fan-Folded                 | 16       |
| 80-13312    |  Non-Sterile, Fan-Folded             | 20       |
| 80-13312-HB |  Non-Sterile, Soft-Folded, Handi-Bin | 240      |

- ▶ SMS Reinforced Fabric
- ▶ Advanced eCoat Technology




#### Zone Reinforced Back Table Cover

Size: 60 x 90 in / 152 x 229 cm

| Code        | Summary                                                                                                               | Case Qty |
|-------------|-----------------------------------------------------------------------------------------------------------------------|----------|
| 80-13313-S  |  Sterile, Fan-Folded                 | 32       |
| 80-13313    |  Non-Sterile, Fan-Folded             | 48       |
| 80-13313-HB |  Non-Sterile, Soft-Folded, Handi-Bin | 280      |

- ▶ SMS Reinforced Fabric
- ▶ Advanced eCoat Technology




#### Zone Reinforced Back Table Cover

Size: 50 x 90 in / 127 x 229 cm

| Code        | Summary                                                                                                                 | Case Qty |
|-------------|-------------------------------------------------------------------------------------------------------------------------|----------|
| 80-13316-S  |  Sterile, Fan-Folded                 | 32       |
| 80-13316    |  Non-Sterile, Fan-Folded             | 50       |
| 80-13316-HB |  Non-Sterile, Soft-Folded, Handi-Bin | 320      |

- ▶ SMS Reinforced Fabric
- ▶ Advanced eCoat Technology



#### Zone Reinforced Back Table Cover

Size: 44 x 78 in / 112 x 198 cm

| Code        | Summary                                                                                                                 | Case Qty |
|-------------|-------------------------------------------------------------------------------------------------------------------------|----------|
| 80-13327-S  |  Sterile, Fan-Folded                 | 46       |
| 80-13327    |  Non-Sterile, Fan-Folded             | 52       |
| 80-13327-HB |  Non-Sterile, Soft-Folded, Handi-Bin | 480      |

- ▶ SMS Reinforced Fabric
- ▶ Advanced eCoat Technology

23 in / 58 cm

54 in / 137 cm



### Mayo Stand Cover

Size: 23 x 54 in / 58 x 137 cm

| Code        | Summary                                                                                                  | Case Qty |
|-------------|----------------------------------------------------------------------------------------------------------|----------|
| 80-13105-S  |  Sterile                | 80       |
| 80-13105    |  Non-Sterile            | 180      |
| 80-13105-HB |  Non-Sterile, Handi-Bin | 1280     |

► SMS Reinforced Fabric

30 in / 76 cm

57 in / 145 cm



### Mayo Stand Cover, X-Large

Size: 30 x 57 in / 76 x 145 cm

| Code        | Summary                                                                                                  | Case Qty |
|-------------|----------------------------------------------------------------------------------------------------------|----------|
| 80-13106-S  |  Sterile                | 60       |
| 80-13106    |  Non-Sterile            | 140      |
| 80-13106-HB |  Non-Sterile, Handi-Bin | 1120     |

► SMS Reinforced Fabric

31.5 in / 80 cm

57 in / 145 cm



### Mayo Stand Cover

Size: 31.5 x 57 in / 80 x 145 cm

| Code        | Summary                                                                                                    | Case Qty |
|-------------|------------------------------------------------------------------------------------------------------------|----------|
| 80-13126-S  |  Sterile                | 80       |
| 80-13126    |  Non-Sterile            | 180      |
| 80-13126-HB |  Non-Sterile, Handi-Bin | N/A      |

► SMS Reinforced Fabric

### Basin Liner Stand Cover, Single

| Code        | Summary                                                                                                    | Case Qty |
|-------------|------------------------------------------------------------------------------------------------------------|----------|
| 80-13506-S  |  Sterile                | N/A      |
| 80-13506    |  Non-Sterile            | 50       |
| 80-13506-HB |  Non-Sterile, Handi-Bin | N/A      |

► Round Basin - Diameter 13.5 in / 34 cm  
► Film Skirt Size - 52 x 52 in / 132 x 132 cm

### Basin Liner Stand Cover, Single, Flat Fold

| Code        | Summary                                                                                                    | Case Qty |
|-------------|------------------------------------------------------------------------------------------------------------|----------|
| 80-13507-S  |  Sterile                | N/A      |
| 80-13507    |  Non-Sterile            | 50       |
| 80-13507-HB |  Non-Sterile, Handi-Bin | N/A      |


### Basin Liner Stand Cover, Double

| Code        | Summary                                                                                                    | Case Qty |
|-------------|------------------------------------------------------------------------------------------------------------|----------|
| 80-13509-S  |  Sterile                | N/A      |
| 80-13509    |  Non-Sterile            | 50       |
| 80-13509-HB |  Non-Sterile, Handi-Bin | N/A      |

► Round Basin - Diameter 13.5 in / 34 cm  
► Film Skirt Size - 52 x 69 in / 132 x 175 cm

## EQUIPMENT COVERS

### Surgical Drapes




#### Rectangular Band Bags

- ▶ Natural Rubber - Latex-Free
- ▶ CPSC Class 1 Flammability Rating
- ▶ ECO2 line made with organic additive to enhance ability of plastic to break down in landfills\*

**Sterile:** 25/Case

**Non-Sterile:** 100/Case

| Width          | Line     | Summary     | 21 in / 71 cm | 30 in / 76 cm | 36 in / 91 cm | 36 in / 91 cm<br>w/Rubber Bands | 40 in / 102 cm | 44 in / 112 cm | 50 in / 127 cm |
|----------------|----------|-------------|---------------|---------------|---------------|---------------------------------|----------------|----------------|----------------|
| 12 in / 30 cm  | Standard | Sterile     | 40-2812-S     | -             | -             | -                               | -              | -              | -              |
|                |          | Non-Sterile | 40-2812       | -             | -             | -                               | -              | -              | -              |
|                | ECO2     | Sterile     | e40-2812-S    | -             | -             | -                               | -              | -              | -              |
|                |          | Non-Sterile | e40-2812      | -             | -             | -                               | -              | -              | -              |
| 20 in / 51 cm  | Standard | Sterile     | -             | 40-3020-S     | 40-3620-S     | -                               | 40-4020-S      | -              | -              |
|                |          | Non-Sterile | -             | 40-3020       | 40-3620       | -                               | 40-4020        | -              | -              |
|                | ECO2     | Sterile     | -             | e40-3020-S    | e40-3620-S    | -                               | e40-4020-S     | -              | -              |
|                |          | Non-Sterile | -             | e40-3020      | e40-3620      | -                               | e40-4020       | -              | -              |
| 30 in / 76 cm  | Standard | Sterile     | -             | 40-3030-S     | 40-3630-S     | -                               | 40-4030-S      | -              | 40-5030-S      |
|                |          | Non-Sterile | -             | 40-3030       | 40-3630       | -                               | 40-4030        | -              | 40-5030        |
|                | ECO2     | Sterile     | -             | -             | -             | -                               | -              | -              | -              |
|                |          | Non-Sterile | -             | -             | -             | -                               | -              | -              | -              |
| 36 in / 91 cm  | Standard | Sterile     | -             | 40-3036-S     | 40-3636-S     | -                               | -              | 40-4036-S      | -              |
|                |          | Non-Sterile | -             | 40-3036       | 40-3636       | -                               | -              | 40-4036        | -              |
|                | ECO2     | Sterile     | -             | e40-3036-S    | e40-3636-S    | -                               | -              | e40-4036-S     | -              |
|                |          | Non-Sterile | -             | e40-3036      | e40-3636      | -                               | -              | e40-4036       | -              |
| 28 in / 71 cm  | Standard | Sterile     | -             | -             | 40-3628-S     | -                               | -              | -              | -              |
|                |          | Non-Sterile | -             | -             | 40-3628       | -                               | -              | -              | -              |
|                | ECO2     | Sterile     | -             | -             | e40-3628-S    | -                               | -              | -              | -              |
|                |          | Non-Sterile | -             | -             | e40-3628      | -                               | -              | -              | -              |
| 40 in / 102 cm | Standard | Sterile     | -             | -             | 40-3640-S     | -                               | 40-4040-S      | -              | 40-5040-S      |
|                |          | Non-Sterile | -             | -             | 40-3640       | -                               | 40-4040        | -              | 40-5040        |
|                | ECO2     | Sterile     | -             | -             | e40-3640-S    | -                               | e40-4040-S     | -              | e40-5040-S     |
|                |          | Non-Sterile | -             | -             | e40-3640      | -                               | e40-4040       | -              | e40-5040       |
| 50 in / 127 cm | Standard | Sterile     | -             | -             | 40-3650-S     | 40-3650R-S                      | -              | -              | -              |
|                |          | Non-Sterile | -             | -             | 40-3650       | 40-3650R                        | -              | -              | -              |
|                | ECO2     | Sterile     | -             | -             | e40-3650-S    | e40-3650R-S                     | -              | -              | -              |
|                |          | Non-Sterile | -             | -             | e40-3650      | e40-3650R                       | -              | -              | -              |
| 54 in / 137 cm | Standard | Sterile     | -             | -             | 40-3654-S     | 40-3654R-S                      | -              | -              | -              |
|                |          | Non-Sterile | -             | -             | 40-3654       | 40-3654R                        | -              | -              | -              |
|                | ECO2     | Sterile     | -             | -             | e40-3654-S    | e40-3654R-S                     | -              | -              | -              |
|                |          | Non-Sterile | -             | -             | e40-3654      | e40-3654R                       | -              | -              | -              |



### Dome Band Bags


- ▶ Natural Rubber - Latex-Free
- ▶ CPSC Class 1 Flammability Rating
- ▶ ECO2 line made with organic additive to enhance ability of plastic to break down in landfills\*

| Size           | SUMMARY     | Standard  | ECO2       | Case Qty |
|----------------|-------------|-----------|------------|----------|
| 14 in / 36 cm  | Sterile     | 40-0014-S | e40-0014-S | 25       |
|                | Non-sterile | 40-0014   | e40-0014   | 100      |
| 18 in / 46 cm  | Sterile     | 40-0018-S | e40-0018-S | 100      |
|                | Non-sterile | 40-0018   | e40-0018   | 25       |
| 22 in / 56 cm  | Sterile     | 40-0022-S | e40-0022-S | 100      |
|                | Non-sterile | 40-0022   | e40-0022   | 25       |
| 26 in / 66 cm  | Sterile     | 40-0026-S | e40-0026-S | 100      |
|                | Non-sterile | 40-0026   | e40-0026   | 25       |
| 27 in / 69 cm  | Sterile     | 40-0027-S | e40-0027-S | 100      |
|                | Non-sterile | 40-0027   | e40-0027   | 25       |
| 30 in / 76 cm  | Sterile     | 40-0030-S | e40-0030-S | 100      |
|                | Non-sterile | 40-0030   | e40-0030   | 25       |
| 40 in / 102 cm | Sterile     | 40-0040-S | e40-0040-S | 25       |
|                | Non-sterile | 40-0040   | e40-0040   | 25       |
| 50 in / 127 cm | Sterile     | 40-0050-S | e40-0050-S | 25       |
|                | Non-sterile | 40-0050   | e40-0050   | 25       |
| 60 in / 152 cm | Sterile     | 40-0060-S | e40-0060-S | 25       |
|                | Non-sterile | 40-0060   | e40-0060   | 25       |

\* Representative products were tested according to ASTM D5511: Standard Test Method for Determining Anaerobic Biodegradation of Plastic Materials Under High-Solids Anaerobic-Digestion Conditions. The results showed a 31% biodegradation of the polyethylene film after 566 days. Actual break-down times will differ based on many variables such as exposure to moisture, air and sunlight.

## EQUIPMENT COVERS

### Surgical Drapes




#### X-Ray Cassette Cover, Medium - Cuffed with Adhesive Closure

Size: 20 x 24 in / 51 x 61 cm

| Code         | Summary                                         | Case Qty |
|--------------|-------------------------------------------------|----------|
| 46-CAS2024-S | <span style="color: red;">S</span> Sterile      | 50       |
| 46-CAS2024   | <span style="color: teal;">N</span> Non-Sterile | 100      |

- ▶ Cuffed
- ▶ Adhesive Closure




#### X-Ray Cassette Cover, Large - Cuffed with Adhesive Closure

Size: 24 x 25 in / 61 x 64 cm

| Code         | Summary                                         | Case Qty |
|--------------|-------------------------------------------------|----------|
| 46-CAS2425-S | <span style="color: red;">S</span> Sterile      | 50       |
| 46-CAS2425   | <span style="color: teal;">N</span> Non-Sterile | 100      |

- ▶ Cuffed
- ▶ Adhesive Closure




#### X-Ray Cassette Cover, Medium - Cuffed Release Liner

Size: 20 x 24 in / 51 x 61 cm

| Code          | Summary                                         | Case Qty |
|---------------|-------------------------------------------------|----------|
| 46-CAS2024L-S | <span style="color: red;">S</span> Sterile      | 50       |
| 46-CAS2024L   | <span style="color: teal;">N</span> Non-Sterile | 100      |

- ▶ Cuffed Release Liner




#### X-Ray Cassette Cover, Large - Cuffed Release Liner

Size: 24 x 25 in / 61 x 64 cm

| Code          | Summary                                         | Case Qty |
|---------------|-------------------------------------------------|----------|
| 46-CAS2425L-S | <span style="color: red;">S</span> Sterile      | 50       |
| 46-CAS2425L   | <span style="color: teal;">N</span> Non-Sterile | 100      |

- ▶ Cuffed Release Liner



#### X-Ray Cassette Cover, Large - Double Cuffed

Size: 24 x 40 in / 61 x 102 cm

| Code          | Summary                                         | Case Qty |
|---------------|-------------------------------------------------|----------|
| 46-CAS2040D-S | <span style="color: red;">S</span> Sterile      | 50       |
| 46-CAS2040D   | <span style="color: teal;">N</span> Non-Sterile | 100      |

- ▶ Double Cuffed

# LEGAL INFORMATION

The following are trademarks of Global Resources International, Inc. in the U.S. and other countries.

The ® symbol refers to a U.S. registration:

- Alleset®
- ComfortGuard®
- DELIVERING PRODUCT SOLUTIONS™
- GRI™
- HydraGuard®
- SteriGuard®
- TECH70™

NOTE: These products may not be available for legal sale in every market. Please check with your sales or account representative for availability of your market of interest.

Copyright © 2017 GRI-Alleset. All rights reserved.

# FACILITIES GRI-ALLESET



LEGEND: ● Office ● Manufacturing ● Distribution ● Corporate Headquarters

## Corporate Office

GRI-Alleset  
4142 Industry Way  
Flowery Branch,  
Georgia 30542  
USA  
Tel: 1.678.866.0550  
Email: [info@gri-alleset.com](mailto:info@gri-alleset.com)

## Europe Office

GRI-Alleset Limited BV  
Kingsfordweg 151  
1043 GR  
Amsterdam  
The Netherlands  
Tel: 31.20.491.9163  
Email: [info@gri-alleset.com](mailto:info@gri-alleset.com)

## India Office

GRI-Alleset Pvt. Ltd.  
331 Arun Chambers  
80 Tardeo Rd.  
Mumbai - 400034  
Maharashtra, India  
Tel: 91.22.23516311  
Email: [info@gri-alleset.com](mailto:info@gri-alleset.com)

## Asia Pacific Office

GRI-Alleset Limited  
27/F, World Wide House,  
19 Des Voeux Road  
Central  
Hong Kong  
Tel: 852.3975.2898  
Email: [info@gri-alleset.com](mailto:info@gri-alleset.com)

## China Office

GRI-Alleset China  
1805 North Honggao Road  
Xiuzhou Industry Zone  
Jiaxing, Zhejiang,  
China 314031  
Tel: 86.573.83916501  
Email: [info@gri-alleset.com](mailto:info@gri-alleset.com)



DELIVERING PRODUCT SOLUTIONS™

[www.gri-alleset.com](http://www.gri-alleset.com)